

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\GRASE\_VASO\_BOLD-TI1\_1460\_TI2\_2600

TA: 0:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L0.7 P39.8 F16.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	3000 ms
TE	59 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
TI 1	1400 ms
TI 2	0 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Pause after meas. 1	0 s
Pause after meas. 2	0 s
Pause after meas. 3	0 s
Pause after meas. 4	0 s
Pause after meas. 5	0 s
Pause after meas. 6	0 s
Pause after meas. 7	0 s
Pause after meas. 8	0 s
Pause after meas. 9	0 s
Pause after meas. 10	0 s
Pause after meas. 11	0 s
Pause after meas. 12	0 s
Pause after meas. 13	0 s
Pause after meas. 14	0 s
Pause after meas. 15	0 s
Pause after meas. 16	0 s
Pause after meas. 17	0 s

Pause after meas. 18	0 s
Pause after meas. 19	0 s
Pause after meas. 20	0 s
Pause after meas. 21	0 s
Pause after meas. 22	0 s
Pause after meas. 23	0 s
Pause after meas. 24	0 s
Multiple series	Off

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm
Position	L0.7 P39.8 F16.9
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.7 P39.8 F16.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	41 mm
F >> H	7 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1144 Hz/Px
<hr/>	
Turbo factor	5
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\BP\_grase\_clean\_VASO\_V07

TA: 0:10 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	3350 ms
TE	58.6 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
TI 1	1400 ms
TI 2	2600 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Off

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
--------	-------------

## Sat. region 1

Thickness	40 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	41 mm
F >> H	7 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1144 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	Mosaic
prepscans	0

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\GRASE\_VASO\_BOLD-TI1\_1400\_TI2\_2500\_20170808

TA: 16:40 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R3.4 P47.4 F27.3
Orientation	T > C-38.3
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	3000 ms
TE	60 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
TI 1	1400 ms
TI 2	2500 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	100
Pause after meas.	0 s
Multiple series	Off

## Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	

Thickness	40 mm
Position	L0.0 P57.8 F14.0
Orientation	C > T38.3
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R3.4 P47.4 F27.3
Orientation	T > C-38.3
Rotation	0.00 deg
R >> L	130 mm
A >> P	41 mm
F >> H	8 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1102 Hz/Px
Echo spacing	1 ms

Turbo factor	6
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast

flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\BP\_grase\_clean\_VASO\_V07\_170912

TA: 3.0 s PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07\_170912

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.8 %
Slice thickness	1.0 mm
TR	3000 ms
TE	52.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
TI 1	1400 ms
TI 2	2600 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

## Resolution

Base resolution	130
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm

Position	Isocenter
Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	40 mm
F >> H	8 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1040 Hz/Px
Echo spacing	1.1 ms
Turbo factor	5
EPI factor	40
RF pulse type	Normal
Gradient mode	Fast
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\BP\_grase\_clean\_VASO\_V07\_10132017

TA: 3.0 s PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07\_101320

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	35.4 %
Slice thickness	1.0 mm
TR	3000 ms
TE	46.5 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1200 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	130
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm
Position	Isocenter

Orientation  
Special sat. Coronal  
None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 Off  
B4 Off  
M3 Off  
V32 Off  
Positioning mode REF  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
? Ref. amplitude 1H 0.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position Isocenter  
Orientation Transversal  
Rotation 0.00 deg  
R >> L 130 mm  
A >> P 46 mm  
F >> H 8 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1424 Hz/Px  
Echo spacing 0.8 ms  
Turbo factor 5  
EPI factor 46  
RF pulse type Normal  
Gradient mode Fast  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\BP\_grase\_clean\_VASO\_V07\_func\_10132017

TA: 2:44 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V07\_101320

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	35.4 %
Slice thickness	1.0 mm
TR	4000 ms
TE	46.5 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1450 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	41
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s

Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Pause after meas. 40	0.0 s
Multiple series	Off

## Resolution

Base resolution	130
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	46 mm
F >> H	8 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1424 Hz/Px
Echo spacing	0.8 ms
<hr/>	
Turbo factor	5
EPI factor	46
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0



# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE\_VASO\_BOLD\BP\_grase\_clean\_VASO\_V08\_func\_10132017

TA: 2:40 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V08\_101320

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.8 %
Slice thickness	1.0 mm
TR	4000 ms
TE	49.6 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1200.0 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s

Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

## Resolution

Base resolution	130
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	66 mm
Position	Isocenter
Orientation	Transversal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	40 mm
F >> H	8 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1098 Hz/Px
Echo spacing	1 ms
<hr/>	
Turbo factor	5
EPI factor	40
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
Phase skip	30
Ampl	100
BWDTH	300 3.1kHz
thickness	100
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	Mosaic
Phase skip	30

## Table of contents

\\USER

Feinberglab

Tanja

GRASE\_VASO\_BOLD

GRASE\_VASO\_BOLD-TI1\_1460\_TI2\_2600

BP\_grase\_clean\_VASO\_V07

GRASE\_VASO\_BOLD-TI1\_1400\_TI2\_2500\_20170808

BP\_grase\_clean\_VASO\_V07\_170912

BP\_grase\_clean\_VASO\_V07\_10132017

BP\_grase\_clean\_VASO\_V07\_func\_10132017

BP\_grase\_clean\_VASO\_V08\_func\_10132017