

## Table of contents

\\USER

VAARC

Joseph

## UC\_Berkeley\_Comparison

localizer  
b1map\_230V  
ep2d\_bold\_sos\_0pt8\_TE29\_ax  
ep2d\_bold\_sos\_1pt0\_TE22\_ax  
ep2d\_bold\_sos\_0pt8\_TE29\_sag  
ep2d\_bold\_sos\_1pt0\_TE22\_sag  
t2\_swi\_tra\_p4\_pt25\_ipat15\_minTE  
T2\_3D\_FLAIR  
T2\_3D\_FLAIR\_p15  
t2\_spc\_TE146\_p15  
t2\_spc\_TE146  
mp2rage\_0.7mm\_TR4500  
ep2d\_bold\_dpg\_1105c\_trig  
ep2d\_bold\_pat3\_sms3\_1.3x1.3x1.3\_dpg  
ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV3\_dpg  
ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV2\_dpg  
ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_dpg  
ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_68PF\_dpg  
ep2d\_bold\_pat5\_sms3\_0.6x0.6x0.6\_dpg  
ep2d\_bold\_pat4\_sms3\_0.6x0.6x0.6\_68pf\_dpg  
ep2d\_bold\_pat4\_sms2\_0.6x0.6x0.6\_68pf\_dpg  
ep2d\_bold\_pat4\_sms3\_0.5x0.5x0.6\_fov3\_58PF\_dpg

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\localizer

TA: 0:28 PM: REF Voxel size: 0.5×0.5×5.0 mmPAT: 3 Rel. SNR: 1.00 : qfl

**Properties**

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	8.6 ms
TE	3.69 ms
Averages	2
Concatenations	15
Filter	Elliptical filter
Coil elements	A32

**Contrast - Common**

TR	8.6 ms
TE	3.69 ms
TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	2
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

**Contrast - Dynamic**

Multiple series	Each measurement
-----------------	------------------

**Resolution - Common**

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	On

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	On

**Geometry - Common**

Slice group	1
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	5
Dist. factor	100 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	8.6 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	15

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P

**Geometry - AutoAlign**

Slice group	2
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**Geometry - Tim CT**

Tim CT mode	Off
Slices	5
Slice thickness	5.0 mm
Dist. factor	100 %
FoV read	250 mm
FoV phase	100.0 %
Segments	1

**System - Miscellaneous**

Positioning mode	REF
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm

**System - Adjust Volume**

R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	8.6 ms
Concatenations	15
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	250 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	15

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	20 deg
Measurements	1
Contrasts	1

**Inline - MapIt**

TR	8.6 ms
TE	3.69 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Contrasts	1
Flow comp.	No
Multi-slice mode	Sequential
Bandwidth	320 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	Active
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

**Sequence - Nuclei**

TX/RX Nucleus	<sup>1</sup> H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\b1map\_230V

TA: 9.2 s PM: FIX Voxel size: 4.0×4.0×4.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	25
Dist. factor	100 %
Position	L0.0 A18.9 F10.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4000.0 ms
TE	1.72 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4000.0 ms
TE	1.72 ms
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	16
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	25
Dist. factor	100 %
Position	L0.0 A18.9 F10.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	4000.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 A18.9 F10.8 mm
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L0.0 A18.9 F10.8
L	0.0 mm
A	18.9 mm
F	10.8 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
--------------	---------

**System - Adjustments**

B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	Off
Dimension	2D
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Interleaved
Echo spacing	3.9 ms
Bandwidth	490 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On
Turbo factor	64

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\lep2d\_bold\_sos\_opt8\_TE29\_ax

TA: 2:56 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 9 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	141
Dist. factor	0 %
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
TR	3990 ms
TE	29.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	3990 ms
TE	29.0 ms
MTC	Off
Flip angle exc	84 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
Base resolution	250
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	141
Dist. factor	0 %
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
TR	3990 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.2 A7.8 H7.8
L	1.2 mm
A	7.8 mm
H	7.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-14.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

**System - Miscellaneous**

Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Rotation	0.00 deg
A >> P	197 mm
R >> L	200 mm
F >> H	113 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3990 ms
Concatenations	1

**BOLD**

Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1112 Hz/Px

**Sequence - Part 2**

EPI factor	246
RF pulse type	Fast
Gradient mode	Fast

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active



\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\lep2d\_bold\_sos\_1pt0\_TE22\_ax

TA: 1:55 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 9 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	114
Dist. factor	0 %
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2640 ms
TE	22.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2640 ms
TE	22.0 ms
MTC	Off
Flip angle exc	70 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
Base resolution	200
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	114
Dist. factor	0 %
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2640 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L1.2 A7.8 H7.8
L	1.2 mm
A	7.8 mm
H	7.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-14.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P

**System - Miscellaneous**

Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.2 A7.8 H7.8 mm
Orientation	T > C-14.4
Rotation	0.00 deg
A >> P	201 mm
R >> L	201 mm
F >> H	114 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2640 ms
Concatenations	1

**BOLD**

Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.83 ms
Bandwidth	1388 Hz/Px

**Sequence - Part 2**

EPI factor	200
RF pulse type	Fast
Gradient mode	Normal

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_sos\_0pt8\_TE29\_sag

TA: 2:56 PM: REF Voxel size: 0.8×0.8×0.8 mmPAT: 9 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	141
Dist. factor	0 %
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
TR	3990 ms
TE	29.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	3990 ms
TE	29.0 ms
MTC	Off
Flip angle exc	84 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
Base resolution	250
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	141
Dist. factor	0 %
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	98.4 %
Slice thickness	0.8 mm
TR	3990 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L2.4 A15.7 F23.5
L	2.4 mm
A	15.7 mm
F	23.5 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares

**System - Miscellaneous**

Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	197 mm
F >> H	200 mm
R >> L	113 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3990 ms
Concatenations	1

**BOLD**

Spatial filter	Off
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1112 Hz/Px

**Sequence - Part 2**

EPI factor	246
RF pulse type	Fast
Gradient mode	Fast

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_sos\_1pt0\_TE22\_sag

TA: 1:51 PM: FIX Voxel size: 1.0×1.0×1.0 mmPAT: 9 Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	114
Dist. factor	0 %
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2530 ms
TE	22.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2530 ms
TE	22.0 ms
MTC	Off
Flip angle exc	70 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
Base resolution	200
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	114
Dist. factor	0 %
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
FoV read	201 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2530 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	L2.4 A15.7 F23.5
L	2.4 mm
A	15.7 mm
F	23.5 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares

**System - Miscellaneous**

Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L2.4 A15.7 F23.5 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	201 mm
F >> H	201 mm
R >> L	114 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2530 ms
Concatenations	1

**BOLD**

Spatial filter	Off
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.83 ms
Bandwidth	1388 Hz/Px

**Sequence - Part 2**

EPI factor	200
RF pulse type	Fast
Gradient mode	Normal

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\t2\_swi\_tra\_p4\_pt25\_ipat15\_minTE

TA: 5:24 PM: FIX Voxel size: 0.4×0.4×1.5 mmPAT: 15 Rel. SNR: 1.00 : fl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Contrast - Common**

Flip angle	15 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

**Contrast - Dynamic**

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

**Routine**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R3.0 A22.3 F3.6 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	12.5 %
Slices per slab	128
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	42.0 ms
TE 1	3.10 ms
TE 2	5.98 ms
TE 3	8.86 ms
TE 4	11.74 ms
TE 5	14.62 ms
TE 6	17.50 ms
TE 7	20.38 ms
TE 8	23.26 ms
TE 9	26.14 ms
TE 10	29.02 ms
TE 11	31.90 ms
TE 12	34.78 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Resolution - Common**

FoV read	220 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
Base resolution	572
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	5
Ref. lines PE	48
Accel. factor 3D	3
Ref. lines 3D	48
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R3.0 A22.3 F3.6 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
Slice oversampling	12.5 %
Slices per slab	128
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	42.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Contrast - Common**

TR	42.0 ms
TE 1	3.10 ms
TE 2	5.98 ms
TE 3	8.86 ms
TE 4	11.74 ms
TE 5	14.62 ms
TE 6	17.50 ms
TE 7	20.38 ms
TE 8	23.26 ms
TE 9	26.14 ms
TE 10	29.02 ms
TE 11	31.90 ms
TE 12	34.78 ms
MTC	Off
Magn. preparation	None

**Geometry - AutoAlign**

Slab group	1
Position	R3.0 A22.3 F3.6 mm
Orientation	T > C-14.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R3.0 A22.3 F3.6
R	3.0 mm
A	22.3 mm
F	3.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-14.4
> S	0.0

**Geometry - Saturation**

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**Geometry - Tim CT**

Tim CT mode	Off
Slabs	1
Slices per slab	128
Slice thickness	1.50 mm
Dist. factor	20 %
FoV read	220 mm
FoV phase	100.0 %
Segments	1

**System - Miscellaneous**

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R3.0 A22.3 F3.6 mm
Orientation	T > C-14.4
Rotation	0.00 deg
A >> P	220 mm
R >> L	220 mm

**System - Adjust Volume**

F >> H	192 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	42.0 ms
Concatenations	1
Segments	1

**Physio - Cardiac**

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	220 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Soft Tissue**

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	None
Flip angle	15 deg
Measurements	1
Contrasts	12
TR	42.0 ms



**Inline - MapIt**

TE 1	3.10 ms
TE 2	5.98 ms
TE 3	8.86 ms
TE 4	11.74 ms
TE 5	14.62 ms
TE 6	17.50 ms
TE 7	20.38 ms
TE 8	23.26 ms
TE 9	26.14 ms
TE 10	29.02 ms
TE 11	31.90 ms
TE 12	34.78 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	12
Flow comp. 1	No
Readout mode	Bipolar
Multi-slice mode	Interleaved
Bandwidth 1	400 Hz/Px
Bandwidth 2	400 Hz/Px
Bandwidth 3	400 Hz/Px
Bandwidth 4	400 Hz/Px
Bandwidth 5	400 Hz/Px
Bandwidth 6	400 Hz/Px
Bandwidth 7	400 Hz/Px
Bandwidth 8	400 Hz/Px
Bandwidth 9	400 Hz/Px
Bandwidth 10	400 Hz/Px
Bandwidth 11	400 Hz/Px
Bandwidth 12	400 Hz/Px

**Sequence - Part 2**

Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On

**Sequence - Nuclei**

TX/RX Nucleus	<sup>1</sup> H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\T2\_3D\_FLAIR

TA: 6:50 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 4 Rel. SNR: 1.00 : spcir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	8000 ms
TE	399 ms
Averages	1.0
Concatenations	1
Filter	Raw filter, Image Filter, B1 filter
Coil elements	A32

**Contrast - Common**

TR	8000 ms
TE	399 ms
MTC	Off
Magn. preparation	Non-sel. T2-IR
TI 1	2150 ms
Fat suppr.	None
Blood suppr.	Off
Restore magn.	Off

**Contrast - Dynamic**

Averages	1.0
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	2
Ref. lines 3D	36
Reordering Shift 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	On
Intensity	Medium
Edge Enhancement	3
Smoothing	3
Unfiltered images	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	On

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	8000 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	R3.0 A22.3 F25.6
R	3.0 mm
A	22.3 mm
F	25.6 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Restore magn.	Off
Special sat.	None

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	157 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	2.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	8000 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. T2-IR
TI 1	2150 ms
Fat suppr.	None
Dark blood	Off
FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	On
Reordering	Linear
Flow comp.	No
Echo spacing	3.3 ms
Adiabatic-mode	Off
Bandwidth	781 Hz/Px

**Sequence - Part 2**

Echo train duration	891 ms
RF pulse type	Normal
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
Turbo factor	300

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

## \\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\T2\_3D\_FLAIR\_p15

TA: 2:58 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 15 Rel. SNR: 1.00 : spcir

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	7.1 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	8000 ms
TE	399 ms
Averages	1.0
Concatenations	1
Filter	Raw filter, Image Filter, B1 filter
Coil elements	A32

**Contrast - Common**

TR	8000 ms
TE	399 ms
MTC	Off
Magn. preparation	Non-sel. T2-IR
TI 1	2150 ms
Fat suppr.	None
Blood suppr.	Off
Restore magn.	Off

**Contrast - Dynamic**

Averages	1.0
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Accel. factor PE	5
Ref. lines PE	48
Accel. factor 3D	3
Ref. lines 3D	48
Reordering Shift 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	On
Intensity	Medium
Edge Enhancement	3
Smoothing	3
Unfiltered images	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	On

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	7.1 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	8000 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	R3.0 A22.3 F25.6
R	3.0 mm
A	22.3 mm
F	25.6 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Restore magn.	Off
Special sat.	None

**Geometry - Navigator**

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R3.0 A22.3 F25.6 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	157 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	2.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	8000 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. T2-IR
TI 1	2150 ms
Fat suppr.	None
Dark blood	Off
FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	On
Reordering	Linear
Flow comp.	No
Echo spacing	3.3 ms
Adiabatic-mode	Off
Bandwidth	781 Hz/Px

**Sequence - Part 2**

Echo train duration	891 ms
RF pulse type	Normal
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
Turbo factor	300

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\t2\_spc\_TE146\_p15

TA: 2:33 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 15 Rel. SNR: 1.00 : spc

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	1 %
Slice oversampling	7.1 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	2290 ms
TE	148 ms
Averages	1.0
Concatenations	1
Filter	Raw filter, B1 filter
Coil elements	A32

**Contrast - Common**

TR	2290 ms
TE	148 ms
MTC	Off
Magn. preparation	None
Fat suppr.	None
Blood suppr.	Off
Restore magn.	Off

**Contrast - Dynamic**

Averages	1.0
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Accel. factor PE	5
Ref. lines PE	48
Accel. factor 3D	3
Ref. lines 3D	48
Reordering Shift 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	7.1 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	2290 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	R3.0 A22.3 F33.7
R	3.0 mm
A	22.3 mm
F	33.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Restore magn.	Off
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	157 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	2.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	2290 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	On
Reordering	Linear
Flow comp.	No
Echo spacing	4.34 ms
Adiabatic-mode	Off
Bandwidth	521 Hz/Px

**Sequence - Part 2**

Echo train duration	352 ms
RF pulse type	Low SAR
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
Turbo factor	96

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\t2\_spc\_TE146

TA: 6:24 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 4 Rel. SNR: 1.00 : spc

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	1 %
Slice oversampling	0.0 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	2290 ms
TE	148 ms
Averages	1.0
Concatenations	1
Filter	Raw filter, B1 filter
Coil elements	A32

**Contrast - Common**

TR	2290 ms
TE	148 ms
MTC	Off
Magn. preparation	None
Fat suppr.	None
Blood suppr.	Off
Restore magn.	Off

**Contrast - Dynamic**

Averages	1.0
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Allowed
Slice partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	CAIPIRINHA
Accel. factor PE	2
Ref. lines PE	48
Accel. factor 3D	2
Ref. lines 3D	48
Reordering Shift 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	On
Unfiltered images	Off

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	224
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	2290 ms
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	R3.0 A22.3 F33.7
R	3.0 mm
A	22.3 mm
F	33.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Saturation**

Fat suppr.	None
Restore magn.	Off
Special sat.	None

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off



**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Off - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R3.0 A22.3 F33.7 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	157 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	2.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
Trigger delay	0 ms
TR	2290 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	On
Reordering	Linear
Flow comp.	No
Echo spacing	4.34 ms
Adiabatic-mode	Off
Bandwidth	521 Hz/Px

**Sequence - Part 2**

Echo train duration	352 ms
RF pulse type	Low SAR
Gradient mode	Fast
Excitation	Non-sel.
Flip angle mode	T2 var
Turbo factor	96

**Sequence - Assistant**

Allowed delay	30 s
---------------	------

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\mp2rage\_0.7mm\_TR4500

TA: 7:50 PM: FIX Voxel size: 0.7×0.7×0.7 mmPAT: 3 Rel. SNR: 1.00 : tfl\_rs

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L1.8 A7.2 F33.1 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	240
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	4500.0 ms
TE	3.37 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4500.0 ms
TE	3.37 ms
Magn. preparation	Non-sel. IR
T1 1	900 ms
T1 2	2750 ms
Flip angle 1	5.0 deg
Flip angle 2	3.0 deg
Fat suppr.	Water excit. fast
Water suppr.	None

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
Base resolution	320
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8

**Resolution - Common**

Slice partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slab group	1
Slabs	1
Dist. factor	50 %
Position	L1.8 A7.2 F33.1 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	240
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	0.70 mm
TR	4500.0 ms
Multi-slice mode	Single shot
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slab group	1
Position	L1.8 A7.2 F33.1 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
AutoAlign	Head > Basis
Initial Position	L1.8 A7.2 F33.1
L	1.8 mm
A	7.2 mm
F	33.1 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Basis
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	L1.8 A7.2 F33.1 mm
Orientation	Sagittal
Rotation	0.00 deg
A >> P	224 mm
F >> H	224 mm
R >> L	168 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	230.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4500.0 ms
Concatenations	1

**Physio - Cardiac**

Magn. preparation	Non-sel. IR
TI 1	900 ms
TI 2	2750 ms
Fat suppr.	Water excit. fast
Dark blood	Off
FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	1

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Inline Composing	Off
Distortion Corr.	Off

**Inline - MapIt**

Save original images	On
MapIt	T1 map
Flip angle 1	5.0 deg
Flip angle 2	3.0 deg
Measurements	1
TR	4500.0 ms
TE	3.37 ms

**Sequence - Part 1**

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Flow comp.	Slice
Multi-slice mode	Single shot
Echo spacing	8.3 ms
Bandwidth	200 Hz/Px

**Sequence - Part 2**

RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	180

**Sequence - Nuclei**

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	A32

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_dpg\_1105c\_trig

TA: 6.1 s PM: REF Voxel size: 3.9×3.9×5.0 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	3
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	196 ms
TE	96.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	196 ms
TE	96.0 ms
MTC	Off
Flip angle exc	90 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	500 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	36
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	3
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	196 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

**System - Miscellaneous**

Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	196 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active

**BOLD**

Motion correction	Off
Spatial filter	Off
Measurements	20
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.52 ms
Bandwidth	752 Hz/Px

**Sequence - Part 2**

EPI factor	128
RF pulse type	Fast
Gradient mode	Fast

**Sequence - Special**

DPG	Off
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_1.3x1.3x1.3\_dpg

TA: 1:08 PM: FIX Voxel size: 1.3x1.3x1.3 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	111
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
TR	2160 ms
TE	23.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2160 ms
TE	23.0 ms
MTC	Off
Flip angle exc	78 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
Base resolution	168
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	111
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
TR	2160 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	210 mm
R >> L	210 mm
F >> H	145 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	170.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2160 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.73 ms
Bandwidth	1654 Hz/Px

**Sequence - Part 2**

EPI factor	168
RF pulse type	Fast
Gradient mode	Normal

**Sequence - Special**

DPG	On
FLEET	On

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV3\_dpg

TA: 1:02 PM: FIX Voxel size: 0.8x0.8x0.8 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	84
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	2470 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2470 ms
TE	30.0 ms
MTC	Off
Flip angle exc	87 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
Base resolution	260
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	84
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	2470 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L



**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	205 mm
R >> L	208 mm
F >> H	68 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2470 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1202 Hz/Px

**Sequence - Part 2**

EPI factor	256
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV2\_dpg

TA: 1:02 PM: FIX Voxel size: 0.8x0.8x0.8 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	84
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	2470 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2470 ms
TE	30.0 ms
MTC	Off
Flip angle exc	87 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
Base resolution	260
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	84
Dist. factor	0 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	2470 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	205 mm
R >> L	208 mm
F >> H	68 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2470 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1202 Hz/Px

**Sequence - Part 2**

EPI factor	256
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\lep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_dpg

TA: 1:46 PM: FIX Voxel size: 0.5x0.5x0.6 mmPAT: 16 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	100
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	2720 ms
TE	26.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2720 ms
TE	26.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	88
Accel. factor slice	4
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	100
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	2720 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	120 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2720 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_68PF\_d  
pg

TA: 2:02 PM: FIX Voxel size: 0.5x0.5x0.6 mmPAT: 16 Rel. SNR: 1.00 : WIPep

### Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

### Routine

Slice group	1
Slices	100
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	3140 ms
TE	42.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

### Contrast - Common

TR	3140 ms
TE	42.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

### Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

### Resolution - Common

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

### Resolution - iPAT

Accel. mode	Slice accel.
-------------	--------------

### Resolution - iPAT

Accel. factor PE	4
Ref. lines PE	88
Accel. factor slice	4
FOV shift factor	3
Reference scan mode	EPI/separate

### Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

### Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

### Geometry - Common

Slice group	1
Slices	100
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	3140 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

### Geometry - AutoAlign

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

### Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

### Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

### System - Miscellaneous

Positioning mode	FIX
Table position	H

**System - Miscellaneous**

Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	120 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3140 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active

**BOLD**

Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat5\_sms3\_0.6x0.6x0.6\_dpg

TA: 1:41 PM: FIX Voxel size: 0.6×0.6×0.6 mmPAT: 15 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	2730 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2730 ms
TE	30.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	5

**Resolution - iPAT**

Ref. lines PE	70
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	2730 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L



**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	108 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2730 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\lep2d\_bold\_pat4\_sms3\_0.6x0.6x0.6\_68pf\_dpg

TA: 1:39 PM: FIX Voxel size: 0.6×0.6×0.6 mmPAT: 12 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	3210 ms
TE	35.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	3210 ms
TE	35.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	56
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	3210 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	108 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3210 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\lep2d\_bold\_pat4\_sms2\_0.6x0.6x0.6\_68pf\_dpg

TA: 1:46 PM: FIX Voxel size: 0.6x0.6x0.6 mmPAT: 8 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	4810 ms
TE	35.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4810 ms
TE	35.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	56
Accel. factor slice	2
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	90
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	4810 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	108 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4810 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	2
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\VAARC\Joseph\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms3\_0.5x0.5x0.6\_fov3\_58PF\_d  
pg

TA: 1:15 PM: FIX Voxel size: 0.5x0.5x1.0 mmPAT: 12 Rel. SNR: 1.00 : WIPep

### Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

### Routine

Slice group	1
Slices	69
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2500 ms
TE	26.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

### Contrast - Common

TR	2500 ms
TE	26.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

### Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

### Resolution - Common

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off

### Resolution - iPAT

Accel. mode	Slice accel.
-------------	--------------

### Resolution - iPAT

Accel. factor PE	4
Ref. lines PE	88
Accel. factor slice	3
FOV shift factor	3
Reference scan mode	EPI/separate

### Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

### Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

### Geometry - Common

Slice group	1
Slices	69
Dist. factor	100 %
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	2500 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

### Geometry - AutoAlign

Slice group	1
Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R2.0 A28.1 F8.5
R	2.0 mm
A	28.1 mm
F	8.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

### Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

### Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

### System - Miscellaneous

Positioning mode	FIX
Table position	H

**System - Miscellaneous**

Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R2.0 A28.1 F8.5 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	137 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210994 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2500 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active

**BOLD**

Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	1
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off