

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\Renzo\_Playground\VASO\_121

TA: 0:10 PAT: Off Voxel size: 9.4x4.7x1.0 mm Rel. SNR: 1.00 USER: VASO\_121

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	32
FoV read	300.0 mm
FoV phase	200.0 %
Slice thickness	1.00 mm
TR	5000.00 ms
TE	52 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Perfusion mode	SS-SI VASO
T12	1100 ms
T11	50 ms
T11s	50 ms
Flip angle	90 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1100.0 ms
Flow limit	100 cm/s

## Resolution

Base resolution	64
Phase resolution	50 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	600 mm
R >> L	300 mm
F >> H	32 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	752 Hz/Px
Free echo spacing	Off
Echo spacing	1.4 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Ampl	100
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	5.2
Renzo: Delta TI	102 ms
EFFECTIVE TR	160000 ms
PatPartitions	32
EPI phase correction	local

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\Renzo\_Playground\VASO\_116\_orig

TA: 0:10 PAT: Off Voxel size: 20.7x6.8x1.0 mm Rel. SNR: 1.00 USER: VASO\_116

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	300.0 mm
FoV phase	300.0 %
Slice thickness	1.00 mm
TR	5000.00 ms
TE	52 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Perfusion mode	Picore Q2TIPS
T12	1100 ms
T11	50 ms
T11s	50 ms
Flip angle	90 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1100.0 ms
Flow limit	100.0 cm/s

## Resolution

Base resolution	44
Phase resolution	33 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	900 mm
R >> L	300 mm
F >> H	8 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	752 Hz/Px
Free echo spacing	Off
Echo spacing	1.43 ms
EPI factor	44
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Ampl	100
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	5.2
Renzo: Delta TI	88 ms
EFFECTIVE TR	40000 ms
PatPartitions	8
EPI phase correction	local

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\Renzo\_Playground\VASO\_116\_phantom

TA: 7.5 s PAT: 2 Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: VASO\_116

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R36.4 P0.0 H40.5
Orientation	T > S31.0
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	32.8 mm
FoV phase	300.0 %
Slice thickness	1.50 mm
TR	1500.00 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Perfusion mode	Picore Q2TIPS
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	90 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

## Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R31.0 A24.3 H24.3
! Orientation	T > S-0.4
! Rotation	180.00 deg
! R >> L	99 mm
! A >> P	82 mm
! F >> H	78 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1062 Hz/Px
Free echo spacing	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Echo spacing	1.07 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	12000 ms
PatPartitions	8
EPI phase correction	local
PAT refscan mode	segm LIN->PAR
use CAIPI	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_PhSk338\_2rep\_V230

TA: 6.0 s PAT: Off Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R29.1 A11.5 H23.0
Orientation	T > S28.5
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.0 %
Slice thickness	1.5 mm
TR	3000 ms
TE	45.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	2
Pause after meas. 1	0.0 s
Multiple series	Off

## Resolution

Base resolution	126
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	32 mm

Position	L1.1 A3.0 H0.5
Orientation	C > S20.0 > T8.1
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R29.1 A11.5 H23.0
Orientation	T > S28.5
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1102 Hz/Px
Echo spacing	1.1 ms

Turbo factor	5
EPI factor	34
RF pulse type	Normal
Gradient mode	Fast*

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR7680

TA: 3.0 s PAT: Off Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR7680

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	28.6 %
Slice thickness	1.5 mm
TR	3000 ms
TE	44.5 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	126
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	66 mm
Position	Isocenter

Orientation  
Special sat. Transversal  
None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	94 mm
A >> P	27 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1168 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	36
RF pulse type	Normal
Gradient mode	Fast
IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0



## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_test

TA: 3.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R29.1 A11.5 H23.0
Orientation	T > S28.5
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	3000 ms
TE	42.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1200 ms
Flip angle	160 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	32 mm
Position	L1.1 A3.0 H0.5

Orientation C > S20.0 > T8.1  
Special sat. NoneTable position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R29.1 A11.5 H23.0
Orientation	T > S28.5
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1086 Hz/Px
Echo spacing	1.1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_func

TA: 12:00 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R34.7 A19.5 H32.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	3000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	240
Pause after meas.	0.0 s
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm

Position	L2.1 A5.8 H0.0
Orientation	C > S20.0
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R34.7 A19.5 H32.0
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\VASO\_116\_phantom

TA: 12:05 PAT: 2 Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: VASO\_116

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R6.6 A8.4 H32.0
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	70.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	35.0 mm
FoV phase	300.0 %
Slice thickness	1.50 mm
TR	1500.00 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Perfusion mode	Picore Q2TIPS
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	483
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

## Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R25.0 A24.3 H16.2
! Orientation	T > S-0.4
! Rotation	180.00 deg
! R >> L	84 mm
! A >> P	82 mm
! F >> H	66 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1062 Hz/Px
Free echo spacing	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Echo spacing	1.07 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	12000 ms
PatPartitions	8
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	44
FlashRef BW	136 Hz/px
FlashRef TE	5000 us
FlashRef FA	5 deg
use CAIPI	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1

TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R34.7 A19.5 H32.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	23 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L2.1 A5.8 H0.0

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R34.7 A19.5 H32.0  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_func

TA: 12:00 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R34.7 A19.5 H32.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	3000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	240
Pause after meas.	0.0 s
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm

Position	L2.1 A5.8 H0.0
Orientation	C > S20.0
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R34.7 A19.5 H32.0
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\VASO\_116\_phantom

TA: 12:05 PAT: 2 Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: VASO\_116

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R6.6 A8.4 H32.0
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	70.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	35.0 mm
FoV phase	300.0 %
Slice thickness	1.50 mm
TR	1500.00 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Perfusion mode	Picore Q2TIPS
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	483
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

## Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

## Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	R25.0 A24.3 H16.2
! Orientation	T > S-0.4
! Rotation	180.00 deg
! R >> L	84 mm
! A >> P	82 mm
! F >> H	66 mm

## Physio

1st Signal/Mode	None
-----------------	------

## BOLD

Motion correction	Off
Spatial filter	Off

## Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1062 Hz/Px
Free echo spacing	Off

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

Echo spacing	1.07 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	12000 ms
PatPartitions	8
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	44
FlashRef BW	136 Hz/px
FlashRef TE	5000 us
FlashRef FA	5 deg
use CAIPI	Off



## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_23

TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	23 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. NoneTable position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast
IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_100					
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400	

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	100 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R32.8 A8.0 H18.5  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_300					
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400	

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	300 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off

Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R32.8 A8.0 H18.5  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_500				
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	500 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R32.8 A8.0 H18.5  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_600				
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	600 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R32.8 A8.0 H18.5  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_700

TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	700 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. NoneTable position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_800  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	800 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast
IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_900				
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	900 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast
IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0



# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1000  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1000 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1 On  
 M2 On  
 B4 On  
 M3 On  
 V32 Off  
 Positioning mode FIX  
 MSMA S - C - T  
 Sagittal R >> L  
 Coronal A >> P  
 Transversal F >> H  
 Save uncombined Off  
 Coil Combine Mode Adaptive Combine  
 AutoAlign ---  
 Auto Coil Select Default

Shim mode Standard  
 Adjust with body coil Off  
 Confirm freq. adjustment Off  
 Assume Silicone Off  
 ! Ref. amplitude 1H 230.000 V  
 Adjustment Tolerance Auto  
 Adjust volume  
 Position R32.8 A8.0 H18.5  
 Orientation Transversal  
 Rotation -20.00 deg  
 R >> L 94 mm  
 A >> P 26 mm  
 F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
 Dimension 3D  
 Reordering Centric  
 Contrasts 2  
 Bandwidth 1116 Hz/Px  
 Echo spacing 1 ms

Turbo factor 5  
 EPI factor 32  
 RF pulse type Normal  
 Gradient mode Fast

IR-RF: 2nd segm phase 338  
 flip angle excit 90  
 phase encoding ON  
 Maxwell compensation Off  
 ICE program single  
 prep scans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1200  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1200 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1 On  
 M2 On  
 B4 On  
 M3 On  
 V32 Off  
 Positioning mode FIX  
 MSMA S - C - T  
 Sagittal R >> L  
 Coronal A >> P  
 Transversal F >> H  
 Save uncombined Off  
 Coil Combine Mode Adaptive Combine  
 AutoAlign ---  
 Auto Coil Select Default

Shim mode Standard  
 Adjust with body coil Off  
 Confirm freq. adjustment Off  
 Assume Silicone Off  
 ! Ref. amplitude 1H 230.000 V  
 Adjustment Tolerance Auto  
 Adjust volume  
 Position R32.8 A8.0 H18.5  
 Orientation Transversal  
 Rotation -20.00 deg  
 R >> L 94 mm  
 A >> P 26 mm  
 F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
 Dimension 3D  
 Reordering Centric  
 Contrasts 2  
 Bandwidth 1116 Hz/Px  
 Echo spacing 1 ms

Turbo factor 5  
 EPI factor 32  
 RF pulse type Normal  
 Gradient mode Fast

IR-RF: 2nd segm phase 338  
 flip angle excit 90  
 phase encoding ON  
 Maxwell compensation Off  
 ICE program single  
 prep scans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1400  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1400 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1 On  
 M2 On  
 B4 On  
 M3 On  
 V32 Off  
 Positioning mode FIX  
 MSMA S - C - T  
 Sagittal R >> L  
 Coronal A >> P  
 Transversal F >> H  
 Save uncombined Off  
 Coil Combine Mode Adaptive Combine  
 AutoAlign ---  
 Auto Coil Select Default

Shim mode Standard  
 Adjust with body coil Off  
 Confirm freq. adjustment Off  
 Assume Silicone Off  
 ! Ref. amplitude 1H 230.000 V  
 Adjustment Tolerance Auto  
 Adjust volume  
 Position R32.8 A8.0 H18.5  
 Orientation Transversal  
 Rotation -20.00 deg  
 R >> L 94 mm  
 A >> P 26 mm  
 F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
 Dimension 3D  
 Reordering Centric  
 Contrasts 2  
 Bandwidth 1116 Hz/Px  
 Echo spacing 1 ms

Turbo factor 5  
 EPI factor 32  
 RF pulse type Normal  
 Gradient mode Fast

IR-RF: 2nd segm phase 338  
 flip angle excit 90  
 phase encoding ON  
 Maxwell compensation Off  
 ICE program single  
 prep scans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1600  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1600 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1 On  
 M2 On  
 B4 On  
 M3 On  
 V32 Off  
 Positioning mode FIX  
 MSMA S - C - T  
 Sagittal R >> L  
 Coronal A >> P  
 Transversal F >> H  
 Save uncombined Off  
 Coil Combine Mode Adaptive Combine  
 AutoAlign ---  
 Auto Coil Select Default

Shim mode Standard  
 Adjust with body coil Off  
 Confirm freq. adjustment Off  
 Assume Silicone Off  
 ! Ref. amplitude 1H 230.000 V  
 Adjustment Tolerance Auto  
 Adjust volume  
 Position R32.8 A8.0 H18.5  
 Orientation Transversal  
 Rotation -20.00 deg  
 R >> L 94 mm  
 A >> P 26 mm  
 F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
 Dimension 3D  
 Reordering Centric  
 Contrasts 2  
 Bandwidth 1116 Hz/Px  
 Echo spacing 1 ms

Turbo factor 5  
 EPI factor 32  
 RF pulse type Normal  
 Gradient mode Fast

IR-RF: 2nd segm phase 338  
 flip angle excit 90  
 phase encoding ON  
 Maxwell compensation Off  
 ICE program single  
 prep scans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1800  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1800 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Rotation	-20.00 deg
R >> L	94 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1116 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

IR-RF: 2nd segm phase	338
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_2000  
 TA: 5.0 s PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_IR6400

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	2000 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
 Special sat. None

Table position H  
 Table position 0 mm  
 Inline Composing Off

## System

T1 On  
 M2 On  
 B4 On  
 M3 On  
 V32 Off  
 Positioning mode FIX  
 MSMA S - C - T  
 Sagittal R >> L  
 Coronal A >> P  
 Transversal F >> H  
 Save uncombined Off  
 Coil Combine Mode Adaptive Combine  
 AutoAlign ---  
 Auto Coil Select Default

Shim mode Standard  
 Adjust with body coil Off  
 Confirm freq. adjustment Off  
 Assume Silicone Off  
 ! Ref. amplitude 1H 230.000 V  
 Adjustment Tolerance Auto  
 Adjust volume  
 Position R32.8 A8.0 H18.5  
 Orientation Transversal  
 Rotation -20.00 deg  
 R >> L 94 mm  
 A >> P 26 mm  
 F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
 Dimension 3D  
 Reordering Centric  
 Contrasts 2  
 Bandwidth 1116 Hz/Px  
 Echo spacing 1 ms

Turbo factor 5  
 EPI factor 32  
 RF pulse type Normal  
 Gradient mode Fast

IR-RF: 2nd segm phase 338  
 flip angle excit 90  
 phase encoding ON  
 Maxwell compensation Off  
 ICE program single  
 prep scans 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASEandRenzosEPI\BP_grase_clean_VASO_V10_IR6400_T1_2200					
TA: 5.0 s	PAT: Off	Voxel size: 0.8x0.8x1.5 mm	Rel. SNR: 1.00	USER: BP_grase_clean_VASO_V10_IR6400	

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R32.8 A8.0 H18.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	27.1 %
Slice thickness	1.5 mm
TR	5000 ms
TE	41.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	2200 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

## Resolution

Base resolution	118
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	L4.0 P5.6 F13.5

Orientation C > S20.0  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode FIX  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
! Ref. amplitude 1H 230.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position R32.8 A8.0 H18.5  
Orientation Transversal  
Rotation -20.00 deg  
R >> L 94 mm  
A >> P 26 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 2  
Bandwidth 1116 Hz/Px  
Echo spacing 1 ms

Turbo factor 5  
EPI factor 32  
RF pulse type Normal  
Gradient mode Fast

IR-RF: 2nd segm phase 338  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0

## Table of contents

\\USER

Feinberglab

Tanja

Renzo\_Playground

VASO\_121

VASO\_116\_orig

VASO\_116\_phantom

GRASEandRenzosEPI

BP\_grase\_clean\_VASO\_V10\_PhSk338\_2rep\_V230

BP\_grase\_clean\_VASO\_V10\_IR7680

BP\_grase\_clean\_VASO\_V10\_IR6400\_test

-----

BP\_grase\_clean\_VASO\_V10\_IR6400\_func

VASO\_116\_phantom

--- 03082018 ---

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1

BP\_grase\_clean\_VASO\_V10\_IR6400\_func

VASO\_116\_phantom

--T1 mapping

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_23

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_100

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_300

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_500

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_600

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_700

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_800

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_900

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1000

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1200

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1400

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1600

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_1800

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_2000

BP\_grase\_clean\_VASO\_V10\_IR6400\_T1\_2200