

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\BP_grase_clean_VASO_V04_pF58_V170

TA: 2.5 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V04

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R67.7 A29.0 F60.0
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	2500 ms
TE	59.78 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1400 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm
Position	R67.7 A24.0 F60.0

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R67.7 A29.0 F60.0
Orientation	Transversal
Rotation	90.00 deg
A >> P	130 mm
R >> L	41 mm
F >> H	7 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1 ms
Turbo factor	5
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast*
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\BP_grase_clean_VASO_V04_pFoff_iPAT2x1_V170

TA: 0:13 PAT: 2 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V04

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L47.0 A29.0 F6.7
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	2500 ms
TE	36.88 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1400 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off

Geometry

Series

Interleaved

Sat. region 1	
Thickness	40 mm
Position	L47.0 A24.0 F6.7
Orientation	Sagittal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L47.0 A29.0 F6.7
Orientation	Transversal
Rotation	90.00 deg
A >> P	130 mm
R >> L	41 mm
F >> H	7 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.1 ms
Turbo factor	8
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast*
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	3

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\BP_grase_clean_VASO

TA: 6.0 s PAT: Off Voxel size: 1.7x1.7x1.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	184 mm
FoV phase	33.3 %
Slice thickness	1.8 mm
TR	2000 ms
TE	35.2 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Off

Resolution

Base resolution	108
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	184 mm
A >> P	62 mm
F >> H	33 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1402 Hz/Px
Echo spacing	0.8 ms
Turbo factor	11
EPI factor	36
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\BP_grase_clean_VASO_V05_pF58_V170

TA: 2.5 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V05

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L45.0 A25.0 H0.0
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	120 mm
FoV phase	30.7 %
Slice thickness	0.8 mm
TR	2500 ms
TE	62.44 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1400 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	150
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm
Position	L45.0 A20.0 H0.0

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L45.0 A25.0 H0.0
Orientation	Transversal
Rotation	90.00 deg
A >> P	120 mm
R >> L	37 mm
F >> H	7 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	952 Hz/Px
Echo spacing	1.2 ms
RF pulse type	Normal
Gradient mode	Fast
Crusher Momentum	40000
Crusher Time	2000
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\pgrs3d_ey

TA: 6.0 s PAT: Off Voxel size: 1.7x1.7x1.8 mm Rel. SNR: 1.00 USER: pgrs3d_ey

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	184 mm
FoV phase	33.3 %
Slice thickness	1.8 mm
TR	2000 ms
TE	76.6 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Off

Resolution

Base resolution	108
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	58 mm
Position	Isocenter
Orientation	Sagittal

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	184 mm
R >> L	62 mm
F >> H	33 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	Off
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1402 Hz/Px
Echo spacing	0.8 ms
Turbo factor	11
EPI factor	37
RF pulse type	Normal
Gradient mode	Fast
Adjust flipangles	Off
Crusher Momentum	40000
Crusher Time	2000
FLIP ANGLES[1]	180 degrees
FLIP ANGLES[2]	180 degrees
FLIP ANGLES[3]	180 degrees
FLIP ANGLES[4]	180 degrees
FLIP ANGLES[5]	180 degrees
FLIP ANGLES[6]	180 degrees
FLIP ANGLES[7]	180 degrees
FLIP ANGLES[8]	180 degrees
FLIP ANGLES[9]	180 degrees
FLIP ANGLES[10]	180 degrees

FLIP ANGLES[11]	180 degrees
FLIP ANGLES[12]	180 degrees
FLIP ANGLES[13]	180 degrees
FLIP ANGLES[14]	180 degrees
FLIP ANGLES[15]	180 degrees
FLIP ANGLES[16]	180 degrees
FLIP ANGLES[17]	180 degrees
FLIP ANGLES[18]	180 degrees
FLIP ANGLES[19]	180 degrees
FLIP ANGLES[20]	180 degrees
FLIP ANGLES[21]	180 degrees
FLIP ANGLES[22]	180 degrees
FLIP ANGLES[23]	180 degrees
FLIP ANGLES[24]	180 degrees
FLIP ANGLES[25]	180 degrees
FLIP ANGLES[26]	180 degrees
FLIP ANGLES[27]	180 degrees
FLIP ANGLES[28]	180 degrees
FLIP ANGLES[29]	180 degrees
FLIP ANGLES[30]	180 degrees
FLIP ANGLES[31]	180 degrees
FLIP ANGLES[32]	180 degrees
Phase Encoding	On
Measurement Number	0
Inversion Time	1500000 us
Inversion Flag	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\fl2d10_retro_1slice_30s

TA: 6:00 PAT: 3 Voxel size: 0.5x0.5x3.0 mm Rel. SNR: 1.00 SIEMENS: CV

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	On
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	On
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	6
Dist. factor	0 %
Position	R1.3 P48.6 F13.5
Orientation	C > T-5.1
Phase enc. dir.	H >> F
Rotation	-89.999728 deg
Auto	On
Phase oversampling	0 %
FoV read	128 mm
FoV phase	100.0 %
Slice thickness	3 mm
TR	1837.82 ms
TE	16.00 ms
Averages	1
Concatenations	6
Filter	Distortion Corr.(2D), Image Filter
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
Magn. preparation	None
Flip angle	60 deg
Fat suppr.	None
Restore magn.	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each slice

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Trajectory	Cartesian
View sharing	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	On
! Intensity	Medium
Edge Enhancement	1
Smoothing	3

Unfiltered images	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	Pulse/Trigger
Average cycle	No Signal ms
Captured cycle	-not set-
Acquisition window	30000 ms
Trigger pulse	1
Trigger delay	0 ms
Segments	86
Phases	16
Tagging	None
Dark blood	Off
Cine	On
Dummy heartbeats	1
Inline ventricular function	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Resp. control	Off
---------------	-----

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	149 Hz/Px
Flow comp.	Yes
Optimization	Min. TR
Allowed delay	0 s
Echo spacing	21.4 ms
Sequence type	Gre
Define	Segments
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slice-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\fl3d_retro_pt5mm_FC

TA: 10:30 PAT: 3 Voxel size: 0.5x0.5x3.0 mm Rel. SNR: 1.00 SIEMENS: CV

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	On
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	On
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R1.3 P38.5 F10.8
Orientation	Coronal
Phase enc. dir.	H >> F
Rotation	-90.00 deg
Auto	On
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	128 mm
FoV phase	100.0 %
Slice thickness	3.00 mm
TR	1842.12 ms
TE	16.00 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Image Filter
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	30 deg
Fat suppr.	None
Restore magn.	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Image Filter	On
! Intensity	Medium
Edge Enhancement	1
Smoothing	3
Unfiltered images	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	ECG/Trigger
Average cycle	No Signal ms
Captured cycle	-not set-
Acquisition window	30000 ms
Trigger pulse	1
Trigger delay	0 ms
Segments	84
Phases	16
Tagging	None

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Dark blood	Off
Cine	On
Dummy heartbeats	1
Inline ventricular function	Off
<hr/>	
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	

Sequence

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	149 Hz/Px
Flow comp.	Yes
Optimization	Min. TR
Allowed delay	0 s
Echo spacing	21.9 ms
Sequence type	Gre
<hr/>	
Define	Segments
RF pulse type	Low SAR
Gradient mode	Normal
Excitation	Slab-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\fl3d_retro_pt5mm_smallFlip

TA: 10:30 PAT: 3 Voxel size: 0.5x0.5x3.0 mm Rel. SNR: 1.00 SIEMENS: CV

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	On
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	On
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R1.3 P38.5 F10.8
Orientation	Coronal
Phase enc. dir.	H >> F
Rotation	-90.00 deg
Auto	On
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	128 mm
FoV phase	100.0 %
Slice thickness	3.00 mm
TR	1816.92 ms
TE	16.00 ms
Averages	1
Concatenations	1
Filter	Distortion Corr.(2D), Image Filter
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	5 deg
Fat suppr.	None
Restore magn.	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Image Filter	On
! Intensity	Medium
Edge Enhancement	1
Smoothing	3
Unfiltered images	Off
Distortion Corr.	On
Mode	2D
Unfiltered images	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	ECG/Trigger
Average cycle	No Signal ms
Captured cycle	-not set-
Acquisition window	30000 ms
Trigger pulse	1
Trigger delay	0 ms
Segments	84
Phases	16
Tagging	None

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Dark blood	Off
Cine	On
Dummy heartbeats	1
Inline ventricular function	Off
<hr/>	
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	

Sequence

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Allowed
Bandwidth	149 Hz/Px
Flow comp.	No
Optimization	Min. TR
Allowed delay	0 s
Echo spacing	21.6 ms
Sequence type	Gre
<hr/>	
Define	Segments
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\VASO_118_broken

TA: 13:16 PAT: 2 Voxel size: 0.7x0.7x1.8 mm Rel. SNR: 1.00 USER: VASO_118

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	120.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	32.8 mm
FoV phase	300.0 %
Slice thickness	1.80 mm
TR	1648.90 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	

Contrast

Perfusion mode	SS-SI VASO
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	30 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	483
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	210.00 deg
R >> L	99 mm
A >> P	33 mm
F >> H	18 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1042 Hz/Px
Free echo spacing	Off
Echo spacing	1.08 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Excitation	Slab-sel.
RF spoiling	On
<hr/>	
Ampl	110
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	5.2
Renzo: Delta TI	72 ms
EFFECTIVE TR	16488 ms
PatPartitions	10
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	44
FlashRef BW	1000 Hz/px
FlashRef TE	4800 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\pgrs3d_ey

TA: 6.0 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: pgrs3d_ey

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	16
FoV read	184 mm
FoV phase	19.1 %
Slice thickness	0.8 mm
TR	2000 ms
TE	54.1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Off

Resolution

Base resolution	230
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	35 mm
Position	Isocenter
Orientation	Coronal

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	184 mm
A >> P	36 mm
F >> H	13 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1358 Hz/Px
Echo spacing	1 ms
Turbo factor	10
EPI factor	44
RF pulse type	Normal
Gradient mode	Fast*
Adjust flipangles	Off
Crusher Momentum	40000
Crusher Time	2000
FLIP ANGLES[1]	180 degrees
FLIP ANGLES[2]	180 degrees
FLIP ANGLES[3]	180 degrees
FLIP ANGLES[4]	180 degrees
FLIP ANGLES[5]	180 degrees
FLIP ANGLES[6]	180 degrees
FLIP ANGLES[7]	180 degrees
FLIP ANGLES[8]	180 degrees
FLIP ANGLES[9]	180 degrees
FLIP ANGLES[10]	180 degrees

FLIP ANGLES[11]	180 degrees
FLIP ANGLES[12]	180 degrees
FLIP ANGLES[13]	180 degrees
FLIP ANGLES[14]	180 degrees
FLIP ANGLES[15]	180 degrees
FLIP ANGLES[16]	180 degrees
FLIP ANGLES[17]	180 degrees
FLIP ANGLES[18]	180 degrees
FLIP ANGLES[19]	180 degrees
FLIP ANGLES[20]	180 degrees
FLIP ANGLES[21]	180 degrees
FLIP ANGLES[22]	180 degrees
FLIP ANGLES[23]	180 degrees
FLIP ANGLES[24]	180 degrees
FLIP ANGLES[25]	180 degrees
FLIP ANGLES[26]	180 degrees
FLIP ANGLES[27]	180 degrees
FLIP ANGLES[28]	180 degrees
FLIP ANGLES[29]	180 degrees
FLIP ANGLES[30]	180 degrees
FLIP ANGLES[31]	180 degrees
FLIP ANGLES[32]	180 degrees
Phase Encoding	On
Measurement Number	0
Inversion Time	1500000 us
Inversion Flag	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\VASO_118

TA: 8.3 s PAT: 2 Voxel size: 0.7x0.7x1.8 mm Rel. SNR: 1.00 USER: VASO_118

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	30.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	32.8 mm
FoV phase	300.0 %
Slice thickness	1.80 mm
TR	1659.70 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Perfusion mode	SS-SI VASO
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	30 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	120.00 deg
A >> P	99 mm
R >> L	33 mm
F >> H	18 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1042 Hz/Px
Free echo spacing	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Echo spacing	1.08 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	100
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	5.2
Renzo: Delta TI	73 ms
EFFECTIVE TR	16597 ms
PatPartitions	10
EPI phase correction	local
PAT refscan mode	segm LIN->PAR
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO\VASO_118_150V

TA: 2:44

PAT: 2

Voxel size: 0.7x0.7x1.8 mm

Rel. SNR: 1.00

USER: VASO_118

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R65.4 P0.0 H4.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	30.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	10
FoV read	32.8 mm
FoV phase	300.0 %
Slice thickness	1.80 mm
TR	2000.00 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Perfusion mode	SS-SI VASO
T12	1200 ms
T11	50 ms
T11s	50 ms
Flip angle	30 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	82
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	1200.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	150.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R65.4 P0.0 H4.0
Orientation	Transversal
Rotation	120.00 deg
A >> P	99 mm
R >> L	33 mm
F >> H	18 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1042 Hz/Px
Free echo spacing	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Echo spacing	1.08 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	100
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	5.2
Renzo: Delta TI	73 ms
EFFECTIVE TR	20000 ms
PatPartitions	10
EPI phase correction	local
PAT refscan mode	segm LIN->PAR
use CAIPI	Off

Table of contents

\\USER

Feinberglab

Tanja

GRASE_VASO

BP_grase_clean_VASO_V04_pF58_V170

BP_grase_clean_VASO_V04_pFoff_iPAT2x1_V170

BP_grase_clean_VASO

BP_grase_clean_VASO_V05_pF58_V170

pgrs3d_ey

fl2d10_retro_1slice_30s

fl3d_retro_pt5mm_FC

fl3d_retro_pt5mm_smallFlip

VASO_118_broken

pgrs3d_ey

VASO_118

VASO_118_150V