

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\localizer_200V_nova

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
----------	-----

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_24SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat.

Coronal

None

Table position

H

Table position

0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	24
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
actual ETL	10
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	98
Variable Flip Angle 02	58

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	53
Variable Flip Angle 04	51
Variable Flip Angle 05	55
Variable Flip Angle 06	60
Variable Flip Angle 07	67
Variable Flip Angle 08	76
Variable Flip Angle 09	90
Variable Flip Angle 10	133
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_24SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat.

Coronal

None

Table position

H

Table position

0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	24
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
actual ETL	10
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	82
Variable Flip Angle 02	47

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	43
Variable Flip Angle 04	42
Variable Flip Angle 05	45
Variable Flip Angle 06	49
Variable Flip Angle 07	57
Variable Flip Angle 08	66
Variable Flip Angle 09	82
Variable Flip Angle 10	130
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_24SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat.

Coronal

None

Table position

H

Table position

0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	24
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
actual ETL	10
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	68
Variable Flip Angle 02	38

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	36
Variable Flip Angle 04	35
Variable Flip Angle 05	38
Variable Flip Angle 06	42
Variable Flip Angle 07	48
Variable Flip Angle 08	58
Variable Flip Angle 09	75
Variable Flip Angle 10	126
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_CFA_24SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	24
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
actual ETL	10
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	180
Variable Flip Angle 02	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_36SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat.

Coronal

None

Table position

H

Table position

0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
actual ETL	14
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	83
Variable Flip Angle 02	46

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	42
Variable Flip Angle 04	40
Variable Flip Angle 05	41
Variable Flip Angle 06	43
Variable Flip Angle 07	46
Variable Flip Angle 08	50
Variable Flip Angle 09	54
Variable Flip Angle 10	59
Variable Flip Angle 11	65
Variable Flip Angle 12	74
Variable Flip Angle 13	89
Variable Flip Angle 14	135
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_36SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
actual ETL	14
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	69
Variable Flip Angle 02	38

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	35
Variable Flip Angle 04	33
Variable Flip Angle 05	34
Variable Flip Angle 06	35
Variable Flip Angle 07	38
Variable Flip Angle 08	41
Variable Flip Angle 09	46
Variable Flip Angle 10	51
Variable Flip Angle 11	58
Variable Flip Angle 12	67
Variable Flip Angle 13	82
Variable Flip Angle 14	130
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_VFA_36SL

TA: 0:00

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	20.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0
actual ETL	14
excite duration	2560
refoc duration	2560
excite BWTP	12
refoc BWTP	8
Variable Flip Angle 01	57
Variable Flip Angle 02	31

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	29
Variable Flip Angle 04	27
Variable Flip Angle 05	29
Variable Flip Angle 06	30
Variable Flip Angle 07	32
Variable Flip Angle 08	35
Variable Flip Angle 09	40
Variable Flip Angle 10	45
Variable Flip Angle 11	53
Variable Flip Angle 12	64
Variable Flip Angle 13	82
Variable Flip Angle 14	136
Crusher Gr	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_CS_CFA_36SL

TA: 4:36

PAT: Off

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: BP_grase_IV_CS_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	92
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat.

Coronal

None

Table position

H

Table position

0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	Mosaic
prep scans	0
actual ETL	6
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 02	180
Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Crusher Gr	40000

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\For_SetupFOV_36Slices

TA: 6.0 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	32.1 %
Slice thickness	0.80 mm
TR	3000 ms
TE	45.34 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	110 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Pause after meas. 1	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	25 mm
Position	Isocenter

Orientation
Special sat. Coronal
None

Table position
Table position
Inline Composing
H
0 mm
Off

System

T1 On
M2 On
B4 On
M3 On
V32 Off
Positioning mode REF
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Standard
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 220.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 90 mm
A >> P 29 mm
F >> H 29 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction Off
Dimension 3D
Reordering Centric
Contrasts 1
Bandwidth 1144 Hz/Px
Echo spacing 1 ms
Turbo factor 36
EPI factor 36
RF pulse type Normal
Gradient mode Fast
refocussing type sinc 2560
flip angle excit 90
phase encoding ON
Maxwell compensation Off
ICE program single
prepscans 0
excite duration 2560
refoc duration 2560
excite BWTP 12.0
refoc BWTP 8.0
T2 Validation Off
pre-crusher 40000

post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	180
Variable Flip Angle 02	180
Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_Regular

TA: 10:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	37.74 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	200
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation
Special sat. Coronal
None

Table position
Table position
Inline Composing
H
0 mm
Off

System

T1 On
M2 On
B4 On
M3 On
V32 Off
Positioning mode REF
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Standard
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 220.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 90 mm
A >> P 23 mm
F >> H 7 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction Off
Dimension 3D
Reordering Centric
Contrasts 1
Bandwidth 1144 Hz/Px
Echo spacing 1 ms
Turbo factor 6
EPI factor 28
RF pulse type Normal
Gradient mode Fast
refocussing type sinc 2560
flip angle excit 90
phase encoding ON
Maxwell compensation Off
ICE program single
prepscans 0
excite duration 2560
refoc duration 2560
excite BWTP 12.0
refoc BWTP 8.0
T2 Validation Off
pre-crusher 40000

post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	180
Variable Flip Angle 02	180
Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_grase_IV_VFA

TA: 10:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	37.74 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	200
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	15 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1 ms

Turbo factor	14
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast

refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
T2 Validation	Off
pre-crusher	40000

post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	90
Variable Flip Angle 02	49
Variable Flip Angle 03	45
Variable Flip Angle 04	43
Variable Flip Angle 05	43
Variable Flip Angle 06	43
Variable Flip Angle 07	42
Variable Flip Angle 08	43
Variable Flip Angle 09	46
Variable Flip Angle 10	51
Variable Flip Angle 11	58
Variable Flip Angle 12	67
Variable Flip Angle 13	77
Variable Flip Angle 14	90
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_VFA_24SL

TA: 10:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	200
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	22 mm
Position	Isocenter

Orientation
Special sat. Coronal
None

Table position H
Table position 0 mm
Inline Composing Off

System

T1 On
M2 On
B4 On
M3 On
V32 Off
Positioning mode REF
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Standard
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 220.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 90 mm
A >> P 23 mm
F >> H 20 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction Off
Dimension 3D
Reordering Centric
Contrasts 1
Bandwidth 1144 Hz/Px
Echo spacing 1.2 ms

Turbo factor 24
EPI factor 12
RF pulse type Normal
Gradient mode Fast

refocussing type variable sinc
flip angle excit 90
phase encoding ON
Maxwell compensation Off
ICE program single
prepscans 0
excite duration 2560
refoc duration 2560
excite BWTP 12.0
refoc BWTP 8.0
Variable Flip Angle 01 82
Variable Flip Angle 02 47

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	43
Variable Flip Angle 04	42
Variable Flip Angle 05	45
Variable Flip Angle 06	49
Variable Flip Angle 07	57
Variable Flip Angle 08	66
Variable Flip Angle 09	82
Variable Flip Angle 10	130
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	10
Which areas?	Motor Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_VFA_36SL

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	2000 ms
TE	24.1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	60
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Pause after meas. 40	0.0 s
Pause after meas. 41	0.0 s
Pause after meas. 42	0.0 s
Pause after meas. 43	0.0 s
Pause after meas. 44	0.0 s
Pause after meas. 45	0.0 s
Pause after meas. 46	0.0 s
Pause after meas. 47	0.0 s
Pause after meas. 48	0.0 s
Pause after meas. 49	0.0 s
Pause after meas. 50	0.0 s
Pause after meas. 51	0.0 s
Pause after meas. 52	0.0 s
Pause after meas. 53	0.0 s
Pause after meas. 54	0.0 s
Pause after meas. 55	0.0 s
Pause after meas. 56	0.0 s
Pause after meas. 57	0.0 s
Pause after meas. 58	0.0 s
Pause after meas. 59	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	22 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

System		Variable Flip Angle 12	74
T1		Variable Flip Angle 13	89
M2		Variable Flip Angle 14	135
B4		Variable Flip Angle 15	180
M3		Variable Flip Angle 16	180
V32		Variable Flip Angle 17	180
		Variable Flip Angle 18	180
		Variable Flip Angle 19	180
		Variable Flip Angle 20	180
Positioning mode		Regular or CS	CS
MSMA		actual ETL	14
Sagittal		Which areas?	Motor Cortex
Coronal			
Transversal			
Save uncombined			
Coil Combine Mode			
AutoAlign			
Auto Coil Select			
Shim mode			
Adjust with body coil			
Confirm freq. adjustment			
Assume Silicone			
! Ref. amplitude 1H			
Adjustment Tolerance			
Adjust volume			
Position			
Orientation			
Rotation			
R >> L			
A >> P			
F >> H			
Physio			
1st Signal/Mode			
Composing			
Sequence			
Introduction			
Dimension			
Reordering			
Contrasts			
Bandwidth			
Echo spacing			
Turbo factor			
EPI factor			
RF pulse type			
Gradient mode			
refocussing type			
flip angle excit			
phase encoding			
Maxwell compensation			
ICE program			
prepscans			
excite duration			
refoc duration			
excite BWTP			
refoc BWTP			
Variable Flip Angle 01			
Variable Flip Angle 02			
Variable Flip Angle 03			
Variable Flip Angle 04			
Variable Flip Angle 05			
Variable Flip Angle 06			
Variable Flip Angle 07			
Variable Flip Angle 08			
Variable Flip Angle 09			
Variable Flip Angle 10			
Variable Flip Angle 11			

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_CFA_18SL

TA: 8:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	24.1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	160
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation

Special sat. Coronal

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	15 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
Turbo factor	18
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	6
Which areas?	Motor Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_Regular

TA: 6.0 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	32.1 %
Slice thickness	0.80 mm
TR	3000 ms
TE	49.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Pause after meas. 1	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	28 mm
Position	Isocenter

Orientation
Special sat. Coronal
None

Table position H
Table position 0 mm
Inline Composing Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	29 mm
F >> H	29 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.1 ms

Turbo factor	27
EPI factor	36
RF pulse type	Normal
Gradient mode	Fast

refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	Regular
actual ETL	6
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_Regular

TA: 3.0 s PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	22
FoV read	112.0 mm
FoV phase	25.0 %
Slice thickness	1.00 mm
TR	1500 ms
TE	36.94 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Pause after meas. 1	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm
Position	Isocenter

Orientation

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	112 mm
R >> L	28 mm
F >> H	22 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1 ms
Turbo factor	17
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
T2 Validation	Off
pre-crusher	40000

post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	180
Variable Flip Angle 02	180
Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_CFA_22SL_0.8mm

TA: 9:48 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	22
FoV read	82.0 mm
FoV phase	37.3 %
Slice thickness	0.80 mm
TR	1500 ms
TE	32.72 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	160 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	392
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	102
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	30 mm
Position	Isocenter

Orientation

Special sat. Sagittal

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	82 mm
R >> L	31 mm
F >> H	18 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1140 Hz/Px
Echo spacing	1.1 ms
Turbo factor	22
EPI factor	20
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	8
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_CFA_22SL_1.0mm

TA: 9:48 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	22
FoV read	102.0 mm
FoV phase	37.3 %
Slice thickness	1.00 mm
TR	1500 ms
TE	32.08 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	160 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	392
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	102
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	37 mm
Position	Isocenter

Orientation

Special sat. Sagittal

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	102 mm
R >> L	38 mm
F >> H	22 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1140 Hz/Px
Echo spacing	1.1 ms
Turbo factor	22
EPI factor	20
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepsans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	8
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_CFA_22SL_1.2mm

TA: 9:48 PAT: Off Voxel size: 1.2x1.2x1.2 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	22
FoV read	122.4 mm
FoV phase	37.3 %
Slice thickness	1.20 mm
TR	1500 ms
TE	31.52 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	160 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	392
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	102
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	45 mm
Position	Isocenter

Orientation

Special sat. Sagittal

Table position	None
Table position	H
Inline Composing	0 mm
	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	123 mm
R >> L	46 mm
F >> H	27 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1140 Hz/Px
Echo spacing	1.1 ms
Turbo factor	22
EPI factor	20
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	8
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_VFA_36SL_1.2mm

TA: 9:48 PAT: Off Voxel size: 1.2x1.2x1.2 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	122.4 mm
FoV phase	37.3 %
Slice thickness	1.20 mm
TR	1500 ms
TE	31.52 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	160 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	392
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	102
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	45 mm
Position	Isocenter

Orientation

Special sat. Sagittal

Table position	None
Table position	H
Inline Composing	0 mm
	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	123 mm
R >> L	46 mm
F >> H	44 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1140 Hz/Px
Echo spacing	1.1 ms
Turbo factor	36
EPI factor	20
RF pulse type	Normal
Gradient mode	Fast
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180

Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	14
Which areas?	Motor Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_Regular_PA

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	37.74 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	180.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	7 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1 ms
<hr/>	
Turbo factor	6
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepScans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
T2 Validation	Off
pre-crusher	40000
post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	90
Variable Flip Angle 02	49
Variable Flip Angle 03	45
Variable Flip Angle 04	43
Variable Flip Angle 05	43
Variable Flip Angle 06	43
Variable Flip Angle 07	42
Variable Flip Angle 08	43
Variable Flip Angle 09	46
Variable Flip Angle 10	51
Variable Flip Angle 11	58
Variable Flip Angle 12	67
Variable Flip Angle 13	77
Variable Flip Angle 14	90
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_grase_IV_VFA_PA

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_Regular_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	37.74 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	180.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	15 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1 ms
<hr/>	
Turbo factor	14
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepScans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
T2 Validation	Off
pre-crusher	40000
post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000
Variable Flip Angle 01	90
Variable Flip Angle 02	49
Variable Flip Angle 03	45
Variable Flip Angle 04	43
Variable Flip Angle 05	43
Variable Flip Angle 06	43
Variable Flip Angle 07	42
Variable Flip Angle 08	43
Variable Flip Angle 09	46
Variable Flip Angle 10	51
Variable Flip Angle 11	58
Variable Flip Angle 12	67
Variable Flip Angle 13	77
Variable Flip Angle 14	90
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_CFA_18SL_PA

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	18
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	180.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	15 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	18
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	180
Variable Flip Angle 02	180
Variable Flip Angle 03	180
Variable Flip Angle 04	180
Variable Flip Angle 05	180
Variable Flip Angle 06	180
Variable Flip Angle 07	180
Variable Flip Angle 08	180
Variable Flip Angle 09	180
Variable Flip Angle 10	180
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	6
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_VFA_24SL_PA

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	24
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	180.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	24
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	82
Variable Flip Angle 02	47
Variable Flip Angle 03	43
Variable Flip Angle 04	42
Variable Flip Angle 05	45
Variable Flip Angle 06	49
Variable Flip Angle 07	57
Variable Flip Angle 08	66
Variable Flip Angle 09	82
Variable Flip Angle 10	130
Variable Flip Angle 11	180
Variable Flip Angle 12	180
Variable Flip Angle 13	180
Variable Flip Angle 14	180
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	10
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\GRASE w/ CS\BP_GRASE_CS_VFA_36SL_PA

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	180.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	variable sinc
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	83
Variable Flip Angle 02	46
Variable Flip Angle 03	42
Variable Flip Angle 04	40
Variable Flip Angle 05	41
Variable Flip Angle 06	43
Variable Flip Angle 07	46
Variable Flip Angle 08	50
Variable Flip Angle 09	54
Variable Flip Angle 10	59
Variable Flip Angle 11	65
Variable Flip Angle 12	74
Variable Flip Angle 13	89
Variable Flip Angle 14	135
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	14
Which areas?	Visual Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_VFA_36SL

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	22.82 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	22 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	83
Variable Flip Angle 02	46
Variable Flip Angle 03	42
Variable Flip Angle 04	40
Variable Flip Angle 05	41
Variable Flip Angle 06	43
Variable Flip Angle 07	46
Variable Flip Angle 08	50
Variable Flip Angle 09	54
Variable Flip Angle 10	59
Variable Flip Angle 11	65
Variable Flip Angle 12	74
Variable Flip Angle 13	89
Variable Flip Angle 14	135
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	12
Which areas?	Motor Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_VFA_36SL

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	24.1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	22 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	3840
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	83
Variable Flip Angle 02	46
Variable Flip Angle 03	42
Variable Flip Angle 04	40
Variable Flip Angle 05	41
Variable Flip Angle 06	43
Variable Flip Angle 07	46
Variable Flip Angle 08	50
Variable Flip Angle 09	54
Variable Flip Angle 10	59
Variable Flip Angle 11	65
Variable Flip Angle 12	74
Variable Flip Angle 13	89
Variable Flip Angle 14	135
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	12
Which areas?	Motor Cortex

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\GRASE w/ CS\BP_GRASE_CS_VFA_36SL

TA: 2:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_GRASE_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	36
FoV read	89.6 mm
FoV phase	25.0 %
Slice thickness	0.80 mm
TR	3000 ms
TE	25.38 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	40
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s

Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Pause after meas. 30	0.0 s
Pause after meas. 31	0.0 s
Pause after meas. 32	0.0 s
Pause after meas. 33	0.0 s
Pause after meas. 34	0.0 s
Pause after meas. 35	0.0 s
Pause after meas. 36	0.0 s
Pause after meas. 37	0.0 s
Pause after meas. 38	0.0 s
Pause after meas. 39	0.0 s
Multiple series	Off

Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	22 mm
Position	Isocenter
Orientation	Coronal
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	90 mm
A >> P	23 mm
F >> H	29 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1144 Hz/Px
Echo spacing	1.2 ms
<hr/>	
Turbo factor	36
EPI factor	12
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	5120
excite BWTP	12.0
refoc BWTP	8.0
Variable Flip Angle 01	83
Variable Flip Angle 02	46
Variable Flip Angle 03	42
Variable Flip Angle 04	40
Variable Flip Angle 05	41
Variable Flip Angle 06	43
Variable Flip Angle 07	46
Variable Flip Angle 08	50
Variable Flip Angle 09	54
Variable Flip Angle 10	59
Variable Flip Angle 11	65
Variable Flip Angle 12	74
Variable Flip Angle 13	89
Variable Flip Angle 14	135
Variable Flip Angle 15	180
Variable Flip Angle 16	180
Variable Flip Angle 17	180
Variable Flip Angle 18	180
Variable Flip Angle 19	180
Variable Flip Angle 20	180
Regular or CS	CS
actual ETL	12
Which areas?	Motor Cortex

Table of contents

\\USER

Feinberglab

Suhjung

GRASE w/ CS

localizer_200V_nova

For Effect of VFA Schemes

BP_grase_IV_CS_VFA_24SL

BP_grase_IV_CS_VFA_24SL

BP_grase_IV_CS_VFA_24SL

BP_grase_IV_CS_CFA_24SL

BP_grase_IV_CS_VFA_36SL

BP_grase_IV_CS_VFA_36SL

BP_grase_IV_CS_VFA_36SL

BP_grase_IV_CS_CFA_36SL

For Comparisons with Current Techinques in M1

For_SetupFOV_36Slices

BP_grase_IV_Regular

BP_grase_IV_VFA

BP_GRASE_CS_VFA_24SL

BP_GRASE_CS_VFA_36SL

BP_GRASE_CS_CFA_18SL

BP_GRASE_Regular

For STG Functional Data

BP_grase_IV_Regular

BP_GRASE_CS_CFA_22SL_0.8mm

BP_GRASE_CS_CFA_22SL_1.0mm

BP_GRASE_CS_CFA_22SL_1.2mm

BP_GRASE_CS_VFA_36SL_1.2mm

--flipped PE dir--

BP_grase_IV_Regular_PA

BP_grase_IV_VFA_PA

BP_GRASE_CS_CFA_18SL_PA

BP_GRASE_CS_VFA_24SL_PA

BP_GRASE_CS_VFA_36SL_PA

BP_GRASE_CS_VFA_36SL

BP_GRASE_CS_VFA_36SL

BP_GRASE_CS_VFA_36SL