

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\SMVMS\gre_PC_SH_goldstandard

TA: 8:40

PAT: Off

Voxel size: 1.1x1.1x5.0 mm

Rel. SNR: 1.00

USER: gre_PC_SH

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-----------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 800 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 2 | |
| Slices | 1 |
| Dist. factor | 800 % |
| Position | L0.0 P0.0 H45.0 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 3 | |
| Slices | 1 |
| Dist. factor | 800 % |
| Position | L0.0 P0.0 H90.0 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 10.0 ms |
| TE | 5.20 ms |
| Averages | 1 |
| Concatenations | 3 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| TD | 0 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 192 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|------------|
| Multi-slice mode | Sequential |
| Series | Ascending |
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ! Ref. amplitude 1H | 150.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|--------------------|---------------|
| 1st Signal/Mode | Pulse/Trigger |
| Average cycle | No Signal ms |
| Captured cycle | -not set- |
| Acquisition window | 900 ms |
| Trigger pulse | 1 |
| Trigger delay | 0 ms |
| Segments | 1 |
| Phases | 16 |
| Tagging | None |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|---------------|-----|
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|------|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 2000 us |
| Slice Distance | 45 mm |
| VENC | 80 cm/s |
| FOV Shift: FOV/ | 1 |
| Random Undersampling | Off |
| Calibration Mode | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\SMVMS\localizer_100V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

| | |
|-----------------------|------|
| Phase resolution | 90 % |
| Phase partial Fourier | 6/8 |
| Interpolation | On |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Sequential |
| Series | Interleaved |
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |
| Tim CT mode | Off |

Routine

| | |
|--------------------|----------------|
| Slice group 1 | |
| Slices | 5 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Sagittal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 2 | |
| Slices | 5 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Coronal |
| Phase enc. dir. | R >> L |
| Rotation | 0.00 deg |
| Slice group 3 | |
| Slices | 5 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 280 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR | 10.0 ms |
| TE | 3.00 ms |
| Averages | 1 |
| Concatenations | 15 |
| Filter | None |
| Coil elements | L1-8;R1-8;T1-8 |

System

| | |
|----|----|
| T1 | On |
| T2 | On |
| T3 | On |
| T4 | On |
| T5 | On |
| T6 | On |
| T7 | On |
| T8 | On |
| L1 | On |
| L2 | On |
| L3 | On |
| L4 | On |
| L5 | On |
| L6 | On |
| L7 | On |
| L8 | On |
| R1 | On |
| R2 | On |
| R3 | On |
| R4 | On |
| R5 | On |
| R6 | On |
| R7 | On |
| R8 | On |

Contrast

| | |
|-------------------|------------------|
| TD | 0 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 10 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

| | |
|-------------------|----------------|
| Positioning mode | FIX |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | On |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |
| Auto Coil Select | Off |

Resolution

| | |
|-----------------|-----|
| Base resolution | 256 |
|-----------------|-----|

| | |
|--------------------------|----------|
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ! Ref. amplitude 1H | 50.000 V |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|----------------------|-------------|
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|------|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|---------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Allowed |
| Bandwidth | 320 Hz/Px |
| Flow comp. | No |
| RF pulse type | Normal |
| Gradient mode | Whisper |
| Excitation | Slice-sel. |
| RF spoiling | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\SMVMS\localizer_100V_newcoil_Calibration

TA: 9.2 s PAT: Off Voxel size: 1.1x1.1x5.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|------------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 700 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 2 | |
| Slices | 1 |
| Dist. factor | 700 % |
| Position | L0.0 P0.0 H40.0 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 3 | |
| Slices | 1 |
| Dist. factor | 700 % |
| Position | L0.0 P0.0 H80.0 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Slice group 4 | |
| Slices | 1 |
| Dist. factor | 700 % |
| Position | L0.0 P0.0 H120.0 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR | 10.0 ms |
| TE | 3.03 ms |
| Averages | 1 |
| Concatenations | 4 |
| Filter | None |
| Coil elements | L1-8;R1-8;T1-8 |

Contrast

| | |
|-------------------|--------|
| TD | 0 ms |
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 10 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |

Averaging mode

| | |
|-----------------|------------------|
| Reconstruction | Short term |
| Measurements | Magnitude |
| Multiple series | 1 |
| | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 192 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | On |

| | |
|----------|------|
| PAT mode | None |
|----------|------|

| | |
|-------------------|-----|
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Sequential |
| Series | Interleaved |

| | |
|-----------------|----------|
| Saturation mode | Standard |
| Special sat. | None |

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

| | |
|-------------|-----|
| Tim CT mode | Off |
|-------------|-----|

System

| | |
|----|----|
| T1 | On |
| T2 | On |
| T3 | On |
| T4 | On |
| T5 | On |
| T6 | On |
| T7 | On |
| T8 | On |
| L1 | On |
| L2 | On |
| L3 | On |
| L4 | On |
| L5 | On |
| L6 | On |
| L7 | On |
| L8 | On |
| R1 | On |
| R2 | On |
| R3 | On |
| R4 | On |
| R5 | On |
| R6 | On |
| R7 | On |
| R8 | On |

| | |
|-------------------|----------------|
| Positioning mode | FIX |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | On |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|--------------------------|-------------|
| Auto Coil Select | Off |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ! Ref. amplitude 1H | 50.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|------|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|---------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Allowed |
| Bandwidth | 320 Hz/Px |
| Flow comp. | No |
| RF pulse type | Normal |
| Gradient mode | Whisper |
| Excitation | Slice-sel. |
| RF spoiling | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS

TA: 2:54 PAT: Off Voxel size: 1.1x1.1x5.0 mm Rel. SNR: 1.00 USER: gre_SMVMS

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR 1 | 13.5 ms |
| TR 2 | 4.5 ms |
| TE | 4.00 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 192 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|-----------------|----------|
| Saturation mode | Standard |
| Special sat. | None |

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|--------------------|---------------|
| 1st Signal/Mode | Pulse/Trigger |
| Average cycle | No Signal ms |
| Captured cycle | -not set- |
| Acquisition window | 900 ms |
| Trigger pulse | 1 |
| Trigger delay | 0 ms |
| Segments | 1 |
| Phases | 16 |

| | |
|------------|------|
| Tagging | None |
| Dark blood | Off |

| | |
|---------------|-----|
| Resp. control | Off |
|---------------|-----|

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|------------|-----|
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| SMS Pulse | On |
| SMS Factor | 4 Bands |
| RF Duration | 2048 us |
| Slice Distance | 40 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Slis. | 90 deg |
| Blips | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhying\SMVMS\gre_SMVMS_IPAT2

TA: 1:28 PAT: 2 Voxel size: 1.1x1.1x5.0 mm Rel. SNR: 1.00 USER: gre_SMVMS

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR 1 | 13.5 ms |
| TR 2 | 4.5 ms |
| TE | 4.00 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|----------|
| Base resolution | 192 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | GRAPPA |
| Accel. factor PE | 2 |
| Ref. lines PE | 24 |
| Reference scan mode | Separate |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |

Elliptical filter

Off

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
| Series | Interleaved |
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|--------------------|---------------|
| 1st Signal/Mode | Pulse/Trigger |
| Average cycle | No Signal ms |
| Captured cycle | -not set- |
| Acquisition window | 900 ms |
| Trigger pulse | 1 |
| Trigger delay | 0 ms |
| Segments | 1 |
| Phases | 16 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|--------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|----------------------|------|
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| <hr/> | |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |
| <hr/> | |
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 4 Bands |
| RF Duration | 2048 us |
| Slice Distance | 40 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\SMVMS\gre_PC_SH_goldstandard

TA: 2:54 PAT: Off Voxel size: 1.1x1.1x5.0 mm Rel. SNR: 1.00 USER: gre_PC_SH

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 700 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 10.0 ms |
| TE | 5.20 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 192 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|------------|
| Multi-slice mode | Sequential |
|------------------|------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ! Ref. amplitude 1H | 150.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|--------------------|---------------|
| 1st Signal/Mode | Pulse/Trigger |
| Average cycle | No Signal ms |
| Captured cycle | -not set- |
| Acquisition window | 900 ms |
| Trigger pulse | 1 |
| Trigger delay | 0 ms |
| Segments | 1 |
| Phases | 16 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|------------|-----|
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 2000 us |
| Slice Distance | 45 mm |
| VENC | 80 cm/s |
| FOV Shift: FOV/ | 1 |
| Random Undersampling | Off |
| Calibration Mode | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\SMVMS\ep2d_bold_SliceAcc_770B

TA: 0:29 PAT: 3 Voxel size: 3.9x3.9x5.0 mm Rel. SNR: 1.00 USER: ep2d_bold_SliceAcc_770B

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 5 |
| Dist. factor | 800 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 500 mm |
| FoV phase | 100.0 % |
| Slice thickness | 5.0 mm |
| TR | 1200 ms |
| TE | 98 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-----------------|-----------|
| MTC | Off |
| Flip angle | 90 deg |
| Fat suppr. | Fat sat. |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 20 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution

| | |
|-----------------------|----------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | GRAPPA |
| Accel. factor PE | 3 |
| Ref. lines PE | 36 |
| Reference scan mode | Separate |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Raw filter | On |
| Elliptical filter | Off |
| Hamming | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
| Series | Interleaved |

Special sat. None

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|-----|-----|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |

| | |
|-------------------|----------------|
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |
| Auto Coil Select | Default |

| | |
|--------------------------|-------------|
| Shim mode | Standard |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 500 mm |
| A >> P | 500 mm |
| F >> H | 185 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
|-----------------|------|

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | On |
| Dynamic t-maps | Off |
| Starting ignore meas | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 20 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-------------------|------------|
| Meas[20] | Active |
| Motion correction | On |
| Interpolation | 3D-K-space |
| Spatial filter | Off |

Sequence

| | |
|---------------------|-----------|
| Introduction | Off |
| Bandwidth | 752 Hz/Px |
| Free echo spacing | Off |
| Echo spacing | 1.4 ms |
| <hr/> | |
| EPI factor | 128 |
| RF pulse type | Normal |
| Gradient mode | Fast |
| <hr/> | |
| Slice acc. factor | 5 |
| RF clip | 0 |
| VERSE factor | 1.00 |
| PE shift factor | 1 |
| Fat sat. flip angle | 110 deg |
| FFT factor | 1.00 |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res35_000

TA: 5.1 s PAT: Off Voxel size: 3.4x3.4x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 64 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | Off |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res35_100

TA: 5.1 s PAT: Off Voxel size: 3.4x3.4x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 64 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res35_110

TA: 5.1 s PAT: Off Voxel size: 3.4x3.4x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 64 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | On |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res23_000

TA: 6.9 s PAT: Off Voxel size: 2.3x2.3x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 96 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|-----------------|----------|
| Saturation mode | Standard |
| Special sat. | None |

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | Off |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res23_100

TA: 6.9 s PAT: Off Voxel size: 2.3x2.3x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 96 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res23_110

TA: 6.9 s PAT: Off Voxel size: 2.3x2.3x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 96 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | On |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res17_000

TA: 8.7 s PAT: Off Voxel size: 1.7x1.7x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | Off |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res17_100

TA: 8.7 s PAT: Off Voxel size: 1.7x1.7x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sl. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | Off |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res17_110

TA: 8.7 s PAT: Off Voxel size: 1.7x1.7x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | On |
| Phase Enc. off | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhjung\SMVMS\gre_SMVMS_NonPhCy_res17_001

TA: 8.7 s PAT: Off Voxel size: 1.7x1.7x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sls. | 90 deg |
| Blips | On |
| RF Spoil | Off |
| VarFA | Off |
| Phase Enc. off | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Suhyung\SMVMS\gre_SMVMS_NonPhCy_res17_101

TA: 8.7 s PAT: Off Voxel size: 1.7x1.7x3.0 mm Rel. SNR: 1.00 USER: gre_SMVMS_NonPhCy

Properties

| | |
|---|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-------------|
| Slice group 1 | |
| Slices | 1 |
| Dist. factor | 20 % |
| Position | Isocenter |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 220 mm |
| FoV phase | 100.0 % |
| Slice thickness | 3.0 mm |
| TR 1 | 14.0 ms |
| TR 2 | 3.5 ms |
| TE | 3.80 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-------------------|------------------|
| MTC | Off |
| Magn. preparation | None |
| Flip angle | 14 deg |
| Fat suppr. | None |
| Water suppr. | None |
| SWI | Off |
| Averaging mode | Short term |
| Reconstruction | Magnitude |
| Measurements | 1 |
| Multiple series | Each measurement |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Image Filter | Off |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Normalize | Off |
| B1 filter | Off |
| Raw filter | Off |
| Elliptical filter | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
|------------------|-------------|

Series

| | |
|------------------|----------|
| Saturation mode | Standard |
| Special sat. | None |
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Tune up |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | Isocenter |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 350 mm |
| A >> P | 263 mm |
| F >> H | 350 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
| Segments | 1 |
| Tagging | None |
| Dark blood | Off |
| Resp. control | Off |

Inline

| | |
|----------------------|-----|
| Subtract | Off |
| Liver registration | Off |
| Std-Dev-Sag | Off |
| Std-Dev-Cor | Off |
| Std-Dev-Tra | Off |
| Std-Dev-Time | Off |
| MIP-Sag | Off |
| MIP-Cor | Off |
| MIP-Tra | Off |
| MIP-Time | Off |
| Save original images | On |
| Wash - In | Off |
| Wash - Out | Off |
| TTP | Off |
| PEI | Off |
| MIP - time | Off |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------|------|
| MapIt | None |
| Contrasts | 1 |

Sequence

| | |
|----------------------|------------|
| Introduction | On |
| Dimension | 2D |
| Phase stabilisation | Off |
| Asymmetric echo | Off |
| Bandwidth | 260 Hz/Px |
| Flow comp. | No |
| <hr/> | |
| RF pulse type | Normal |
| Gradient mode | Fast |
| Excitation | Slice-sel. |
| RF spoiling | On |
| <hr/> | |
| SMS Pulse | On |
| SMS Factor | 1 Bands |
| RF Duration | 1600 us |
| Slice Distance | 30 mm |
| High VENC | 80 cm/s |
| Low VENC | 40 cm/s |
| Random UnderSampling | Off |
| Phs Diff. Bet. Sl. | 90 deg |
| Blips | On |
| RF Spoil | On |
| VarFA | Off |
| Phase Enc. off | On |

Table of contents

\\USER

Feinberglab

Suhjung

SMVMS

gre_PC_SH_goldstandard
 localizer_100V_newcoil
 localizer_100V_newcoil_Calibration
 gre_SMVMS
 gre_SMVMS_IPAT2
 gre_PC_SH_goldstandard
 ep2d_bold_SliceAcc_770B
 gre_SMVMS_NonPhCy_res35_000
 gre_SMVMS_NonPhCy_res35_100
 gre_SMVMS_NonPhCy_res35_110
 gre_SMVMS_NonPhCy_res23_000
 gre_SMVMS_NonPhCy_res23_100
 gre_SMVMS_NonPhCy_res23_110
 gre_SMVMS_NonPhCy_res17_000
 gre_SMVMS_NonPhCy_res17_100
 gre_SMVMS_NonPhCy_res17_110
 gre_SMVMS_NonPhCy_res17_001
 gre_SMVMS_NonPhCy_res17_101
 gre_SMVMS_NonPhCy_res17_111

b-SMS

gre_PC_SH_calibration
 gre_PC_SH_mb5_fov3
 gre_PC_SH_mb5_fov3_pi3_rand
 gre_PC_SH_mb5_fov3_pi4_rand
 ep2d_bold_SliceAcc_770B

SMS_EPI

localizer_50V_newcoil
 AV_ep2d_bold_sd1ipat3mb4_pt75mm_tSNR_shimWholeVol
 AV_ep2d_bold_sd1ipat3mb5_pt75mm_tSNR_shimWholeVol
 AV_ep2d_bold_sd1ipat3mb6_pt75mm_tSNR_shimWholeVol

GRASE

localizer_50V_newcoil
 grase3d_reference
 grase3d_cs_SH
 grase3d_cspf_SH
 grase3d_cs_ms_SH
 grase3d_cspf_ms_SH
 grase3d_reference_res15
 grase3d_cspf_ms_res15_SH

GRASE_IV

localizer_200V_nova
 b1map_200V_TR1000_nova
 Visual Cortex
 measure point spread functions
 BP_grase_clean_IV_Regular_PSF_Par8_SH
 BP_grase_clean_IV_Regular_PSF_Par20_SH
 BP_grase_clean_IV_VFL_PSF_SH
 BP_grase_clean_IV_CS1D_Par20_NonSFL_PSF_SH
 BP_grase_clean_IV_CS2D_Par20_NonSFL_PSF_SH
 resolution: iso-0.8mm3
 BP_grase_clean_IV_TE37_Regular_SH
 BP_grase_clean_IV_TE37_VFL_SH
 BP_grase_clean_IV_TE37_Par20_STE_SH
 BP_grase_clean_IV_TE27_Par20_STE_SH
 2D MultiContrast GRASE for T2 Map
 BP_grase_clean_IV_Regular_MC_CGR_FA180
 BP_grase_clean_IV_Regular_MC_VGR_FA180
 BP_grase_clean_IV_Regular_MC_CGR_FA140
 BP_grase_clean_IV_Regular_MC_VGR_FA140
 Comparisons: STE vs. nSTE

\\USER

Feinberglab

Suhung

GRASE_IV

BP_grase_clean_IV_TE41_FA180_Regular_CGR
 BP_grase_clean_IV_TE41_FA180_Regular_VGR
 BP_grase_clean_IV_FA_TE32_FA180_CS_CGR
 BP_grase_clean_IV_FA_TE32_FA180_CS_VGR
 Pause

BP_grase_clean_IV_TE41_FA140_Regular_CGR
 BP_grase_clean_IV_TE41_FA180_Regular_VGR
 BP_grase_clean_IV_TE41_FA140_Regular_VGR
 for tSNR Test

BP_grase_clean_IV_TE32_FA180_CS_CGR
 BP_grase_clean_IV_TE32_FA180_CS_CGR
 BP_grase_clean_IV_FA_TE32_FA180_CS_CGR
 etc.....

AV_ep2d_bold_sd1ipat2mb1_1mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb1_p8mm_tSNR_32
 cmrr_mbep2d_se

--M1 STE--

BP_grase_clean_IV_TE37_STE_Regular_SH
 BP_grase_clean_IV_TE37_nSTE_Regular_SH
 BP_grase_clean_IV_TE27_Par20_STE_SH
 BP_grase_clean_IV_TE27_Par20_nSTE_SH

--- for nSTE Test

BP_grase_clean_IV_FA180_Regular_SH
 BP_grase_clean_IV_FA180_nSTE_CurrentPT_SH
 BP_grase_clean_IV_FA140_Regular_SH
 BP_grase_clean_IV_FA140_nSTE_CurrentPT_SH
 BP_grase_clean_IV_FA160_Regular_SH
 BP_grase_clean_IV_FA160_nSTE_CurrentPT_SH

Pause

Motor Cortex

BP_grase_clean_IV_Reg_ForComparison
 BP_grase_clean_IV_Reg_nSTE_T2Map
 BP_grase_clean_IV_Reg
 BP_grase_clean_IV_Reg_nSTE
 BP_grase_clean_IV_Reg_VFA
 BP_grase_clean_IV_Reg_nSTE_VFA
 BP_grase_clean_IV_CSv1
 BP_grase_clean_IV_CSv1_nSTE

Motor Cortex for PSF

BP_grase_clean_IV_Reg_PSF_20slc
 BP_grase_clean_IV_Reg_PSF
 BP_grase_clean_IV_Reg_VFA_PSF
 BP_grase_clean_IV_CSv1_PSF
 BP_grase_clean_IV_Rand_SH

GRASE_IV_TSNR

resolution: iso-0.8mm3

BP_grase_clean_IV_TE41_PF6/8_Res08_SH
 BP_grase_clean_IV_TE41_PF6/8_Res08_VFA_SH
 BP_grase_clean_IV_TE41_PF6/8_Res08_VFA_SH
 BP_grase_clean_IV_CS_TE41_PF6/8_Res08_SH
 BP_grase_clean_IV_CS_TE35_PF6/8_Res08_SH
 BP_grase_clean_IV_CS_TE38_PF6/8_Res08_fov48_SH
 pgrs3d_ey_20180816_p8mm

CS-M1

BP_grase_clean_IV_Regular_SH_ob
 BP_grase_clean_IV_CS2D_Par20_NonSFL_SH_Ob
 --time permitting--
 BP_grase_clean_VASO_V10t_noClip
 BP_grase_clean_IV_Regular_TE43_SH_TR3000
 --02/18--

\\USER

Feinberglab

Suhjung

CS-M1

BP_grase_clean_IV_Regular_SH_dbIOb
BP_grase_clean_IV_CS2D_Par20_NonSFL_SH_dbIOb

GRASE_IV_CS_GRASE

localizer_200V_nova
t1_mpr_sag_p9mm_iso
b1map_200V_TR1000_nova
mp2rage_0.8mmiso_TR4500
----- For STE vs nSTE
BP_grase_IV_Regular
BP_grase_IV_Regular_nSTE
BP_grase_IV_Regular
BP_grase_IV_Regular_nSTE
BP_grase_IV_Regular_nSTE_T2map_ICEOff
se_mc
----- For Point Spread Functions -----
BP_grase_IV_Regular
BP_grase_IV_VFA
BP_grase_IV_CS_CFA_SL18_ETL06
BP_grase_IV_CS_VFA_SL24_ETL10
BP_grase_IV_CS_VFA_SL36_ETL14
----- For Visual Cortex -----
BP_grase_IV_Regular
BP_grase_IV_VFA
BP_grase_IV_CS_CFA_SL18_ETL06
BP_grase_IV_CS_VFA_SL24_ETL10
BP_grase_IV_CS_VFA_SL36_ETL14
BP_grase_IV_CS_VFA_SL48_ETL14
-----For Motor Cortex -----
BP_grase_IV_Regular
BP_grase_IV_VFA
BP_grase_IV_CS_CFA_SL18_ETL06
BP_grase_IV_CS_VFA_SL24_ETL10
BP_grase_IV_CS_VFA_SL36_ETL14
----- For STG -----
BP_grase_IV_CS_VFA_SL36_TE22
BP_grase_IV_CS_VFA_SL36_TE27
For 3D EPI
BP_ep3D_bold_CAIPi_ICEOff_ChannelM2
BP_ep3D_bold_multiecho_new_SH_CAIPi_noICE
- For Motor Cortex: turn off phase encode
-----Point Spread Function for Alex- turn off phase encode
BP_grase_clean_VASO_V10t_noClip
VASO_116_phantom
Enhance Tissue Contrast For Sam
BP_grase_IV_CS_VFA_SL36_ETL14_TE22_TR12
BP_grase_IV_CS_VFA_SL36_ETL14_TE22_TR20
-----3DEPI CAIPi: turn on phase encode, ice off for last
GE_p8mm_MB1IPAT3_pf6_te23_tr1505_sat_PIFrSn_391i
BP_ep3D_bold_multiecho_new_SH_noCAIPi

GRASE VASO w/ CS

BP_grase_IV_Regular_VASO_SH
BP_GRASE_VASO_Short
BP_GRASE_VASO_Long
Pause
BP_GRASE_VASO_BloodNull_IR100ms
BP_GRASE_VASO_BloodNull_IR200ms
BP_grase_IV_VFA_VASO_SH
BP_grase_IV_CS_VASO_SH_24SL
BP_grase_IV_CS_VFA_VASO_GMS_SH
BP_grase_IV_CS_VASO_SH

\\USER

Feinberglab

Suhyung

GRASE w/ CS

localizer_200V_nova
 For Effect of VFA Schemes
 BP_grase_IV_CS_VFA_24SL
 BP_grase_IV_CS_VFA_24SL
 BP_grase_IV_CS_VFA_24SL
 BP_grase_IV_CS_CFA_24SL
 BP_grase_IV_CS_VFA_36SL
 BP_grase_IV_CS_VFA_36SL
 BP_grase_IV_CS_VFA_36SL
 BP_grase_IV_CS_CFA_36SL
 For Comparisons with Current Techinques in M1
 For_SetupFOV_36Slices
 BP_grase_IV_Regular
 BP_grase_IV_VFA
 BP_GRASE_CS_VFA_24SL
 BP_GRASE_CS_VFA_36SL
 BP_GRASE_CS_CFA_18SL
 BP_GRASE_Regular
 For STG Functional Data
 BP_grase_IV_Regular
 BP_GRASE_CS_CFA_22SL_0.8mm
 BP_GRASE_CS_CFA_22SL_1.0mm
 BP_GRASE_CS_CFA_22SL_1.2mm
 BP_GRASE_CS_VFA_36SL_1.2mm
 --flipped PE dir--
 BP_grase_IV_Regular_PA
 BP_grase_IV_VFA_PA
 BP_GRASE_CS_CFA_18SL_PA
 BP_GRASE_CS_VFA_24SL_PA
 BP_GRASE_CS_VFA_36SL_PA
 BP_GRASE_CS_VFA_36SL
 BP_GRASE_CS_VFA_36SL
 BP_GRASE_CS_VFA_36SL

GRASE w/ CS for STG

BP_RGRASE_24SL_1.0mm
 BP_AccVGRASE_CFA_24SL_1.0mm_EPI20
 BP_AccVGRASE_CFA_24SL_1.0mm_EPI16
 BP_AccVGRASE_24SL_1.0mm_EPI20
 BP_AccVGRASE_24SL_1.0mm_EPI16
 BP_AccVGRASE_24SL_1.0mm_EPI20_RefDuration2560
 BP_AccVGRASE_36SL_1.0mm_EPI20
 BP_AccVGRASE_36SL_0.8mm
 Functional Acquisitions
 BP_AccCGRASE_24SL_1.0mm_TE23ms
 BP_AccVGRASE_24SL_1.0mm_TE23ms
 BP_AccVGRASE_24SL_1.0mm_TE21ms
 BP_AccVGRASE_36SL_1.0mm_TE23ms
 Outer VolumeSuppressions
 ep2d_M1P2f1_iso150
 ep2d_venc_ms_sbmb_SAT_flashref
 ep2d_venc_ms_sbmb_SAT_flashref_SH