

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\EPI_VASO_BOLD\ep2d_fid_VASO

TA: 1:12 PAT: 2 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: ep2d_fid_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.0 mm
TR	4000.0 ms
TE 1	42 ms
TE 2	42 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	Non-sel. IR
T1	1450 ms
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	15
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode Interleaved
Series Interleaved

Special sat. None
Table position H
Table position 0 mm
Inline Composing Off

System

T1 On
M2 On
B4 On
M3 On
V32 Off
Positioning mode REF
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Sum of Squares
AutoAlign ---
Auto Coil Select Default
Shim mode Standard
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
? Ref. amplitude 1H 0.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 192 mm
A >> P 192 mm
F >> H 1 mm

Physio

1st Signal/Mode None

Perf

GBP On
PBP On
TTP On
Original images On
Starting ignore meas 2

Sequence

Introduction Off
Contrasts 2
Bandwidth 766 Hz/Px
Free echo spacing Off
Echo spacing 1.39 ms
EPI factor 192
RF pulse type Fast
Gradient mode Fast

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\EPI_VASO_BOLD\BP_grase_clean_VASO_V07_func_10132017

TA: 5:24 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V07_101320

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L0.8 A12.2 H37.8
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	162 mm
FoV phase	37.0 %
Slice thickness	1.0 mm
TR	4000 ms
TE	69.4 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1450 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	81
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm

Position	L0.8 A12.2 H37.8
Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.8 A12.2 H37.8
Orientation	Transversal
Rotation	0.00 deg
R >> L	162 mm
A >> P	60 mm
F >> H	8 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1342 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	60
RF pulse type	Normal
Gradient mode	Fast

flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0

Table of contents

\\USER

Feinberglab

Tanja

EPI_VASO_BOLD

ep2d_fid_VASO

BP_grase_clean_VASO_V07_func_10132017