

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO_FOCI_20170831\BP_grase_clean_VASO_V08

TA: 0:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V08

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L3.8 P36.0 F12.8
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.9 %
Slice thickness	0.8 mm
TR	3000 ms
TE	59 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
TI 1	1030 ms
TI 2	2530 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	2
Pause after meas. 1	0 s
Multiple series	Off

Resolution

Base resolution	162
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	

Thickness	40 mm
Position	L5.4 P36.0 H20.9
Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	225.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L3.8 P36.0 F12.8
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	41 mm
F >> H	7 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1144 Hz/Px

Turbo factor	5
EPI factor	50
RF pulse type	Normal
Gradient mode	Fast

Phase skip	0
Ampl	1
BWDTH	1 3.1kHz
thickness	0
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	Mosaic
Phase skip	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_VASO_FOCI_20170831\BP_grase_clean_VASO_V08_170912

TA: 3.0 s PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V08_170912

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	130 mm
FoV phase	30.8 %
Slice thickness	1.0 mm
TR	3000 ms
TE	52.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
TI 1	1400.0 ms
TI 2	2600 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	130
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	40 mm

Position	Isocenter
Orientation	Coronal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	130 mm
A >> P	40 mm
F >> H	8 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1040 Hz/Px
Echo spacing	1.1 ms
Turbo factor	5
EPI factor	40
RF pulse type	Normal
Gradient mode	Fast
Phase skip	30
Ampl	100
BWDTH	300 3.1kHz
thickness	100
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
Phase skip	30

Table of contents

\\USER

Feinberglab

Tanja

GRASE_VASO_FOCI_20170831

BP_grase_clean_VASO_V08

BP_grase_clean_VASO_V08_170912