

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_EPI_VASO_BIR4\BP_grase_clean_VASO_V10

TA: 3.0 s PAT: Off Voxel size: 7.8x7.8x3.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V10

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	3000 ms
TE	42.1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

Contrast

Magn. preparation	Non-sel. IR
T1	1200 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	66 mm
Position	Isocenter

Orientation
Special sat. Transversal
None

Table position H
Table position 0 mm
Inline Composing Off

System

T1 On
M2 Off
B4 Off
M3 Off
V32 Off
Positioning mode REF
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Standard
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
? Ref. amplitude 1H 0.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 500 mm
A >> P 500 mm
F >> H 24 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction Off
Dimension 3D
Reordering Centric
Contrasts 2
Bandwidth 2004 Hz/Px
Echo spacing 0.5 ms

Turbo factor 5
EPI factor 64
RF pulse type Normal
Gradient mode Fast

IR-RF: 2nd segm phase 0
flip angle excit 90
phase encoding ON
Maxwell compensation Off
ICE program single
prepscans 0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Tanja\GRASE_EPI_VASO_BIR4\ep_seg_fid_VASO_01

TA: 0:40

Voxel size: 4.0x3.9x1.0 mm

Rel. SNR: 1.00

USER: ep_seg_fid_VASO_01

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	4000 ms
TE 1	91 ms
TE 2	1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

Contrast

MTC	Off
Magn. preparation	Non-sel. IR
T1	1600 ms
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	98 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
------------------	-------------

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	500 mm
A >> P	500 mm
F >> H	8 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Contrasts	2
Bandwidth	752 Hz/Px
Free echo spacing	Off
Echo spacing	1.38 ms
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
IR-RF: 2nd segm phase	0

Table of contents

\\USER

Feinberglab

Tanja

GRASE_EPI_VASO_BIR4

BP_grase_clean_VASO_V10

ep_seg_fid_VASO_01