

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Suhyung\CS-M1\BP\_grase\_clean\_IV\_Regular\_SH\_ob

TA: 0:00 PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_IV\_SH

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	25.0 %
Slice thickness	0.8 mm
TR	3000 ms
TE	37.7 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	121
Pause after meas.	0 s
Multiple series	Off

## Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation Coronal  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode REF  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
? Ref. amplitude 1H 0.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position Isocenter  
Orientation Transversal  
Rotation 0.00 deg  
R >> L 94 mm  
A >> P 24 mm  
F >> H 7 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 1  
Bandwidth 1144 Hz/Px  
Turbo factor 5  
EPI factor 28  
RF pulse type Normal  
Gradient mode Fast  
refocussing type sinc 2560  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0  
excite duration 2560  
refoc duration 2560  
excite BWTP 16  
refoc BWTP 8  
Opposite Polarity Crusher Off  
pre-crusher 0  
post-crusher1 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

post-crusher2	0
post-crusher3	0
post-crusher4	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Suhyung\CS-M1\BP\_grase\_clean\_VASO\_V10t\_noClip

TA: 12:03 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_VASO\_V10\_noClip

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R27.3 P0.8 H16.6
Orientation	T > S24.2 > C-18.7
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	99 mm
FoV phase	25.8 %
Slice thickness	1.5 mm
TR	3000 ms
TE	45.9 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	241
Pause after meas.	0.0 s
Multiple series	Off

## Resolution

Base resolution	132
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm

Position	L10.8 P8.8 F15.4
Orientation	C > S26.5 > T17.3
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	215.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R27.3 P0.8 H16.6
Orientation	T > S24.2 > C-18.7
Rotation	-20.00 deg
R >> L	99 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode	None
-----------------	------

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1052 Hz/Px
Echo spacing	1.1 ms

Turbo factor	5
EPI factor	34
RF pulse type	Normal
Gradient mode	Fast

BIR4: 2nd segm phase	338
BIR4: duration	5120
excite duration	2560
refoc duration	2560
excite BWTP	10.4
refoc BWTP	5.2
phase encoding	ON
Maxwell compensation	Off
ICE program	single

## SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Suhyung\CS-M1\BP\_grase\_clean\_IV\_Regular\_TE43\_SH\_TR3000

TA: 0:00 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_IV\_SH

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R27.3 P0.8 H16.6
Orientation	T > S24.2 > C-18.7
Phase enc. dir.	A >> P
Rotation	-20 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	99 mm
FoV phase	25.8 %
Slice thickness	1.5 mm
TR	3000 ms
TE	45.92 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	241
Pause after meas.	0 s
Multiple series	Off

## Resolution

Base resolution	132
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm
Position	L10.8 P8.8 F15.4

Orientation C > S26.5 > T17.3  
Special sat. NoneTable position H  
Table position 0 mm  
Inline Composing Off

## System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	215.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R27.3 P0.8 H16.6
Orientation	T > S24.2 > C-18.7
Rotation	-20.00 deg
R >> L	99 mm
A >> P	26 mm
F >> H	12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1052 Hz/Px

Turbo factor	5
EPI factor	34
RF pulse type	Normal
Gradient mode	Fast

refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	10.4
refoc BWTP	5.2
Opposite Polarity Crusher	Off
pre-crusher	0
post-crusher1	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

post-crusher2	0
post-crusher3	0
post-crusher4	0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

\\USER\Feinberglab\Suhyung\CS-M1\BP\_grase\_clean\_IV\_Regular\_SH\_dbIOb

TA: 0:00 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP\_grase\_clean\_IV\_SH

## Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

## Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	94 mm
FoV phase	25.0 %
Slice thickness	1.5 mm
TR	3000 ms
TE	37.54 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

## Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	241
Pause after meas.	0 s
Multiple series	Off

## Resolution

Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

## Geometry

Series	Interleaved
Sat. region 1	
Thickness	20 mm
Position	Isocenter

Orientation Coronal  
Special sat. None

Table position H  
Table position 0 mm  
Inline Composing Off

## System

T1 On  
M2 On  
B4 On  
M3 On  
V32 Off  
Positioning mode REF  
MSMA S - C - T  
Sagittal R >> L  
Coronal A >> P  
Transversal F >> H  
Save uncombined Off  
Coil Combine Mode Adaptive Combine  
AutoAlign ---  
Auto Coil Select Default

Shim mode Standard  
Adjust with body coil Off  
Confirm freq. adjustment Off  
Assume Silicone Off  
? Ref. amplitude 1H 0.000 V  
Adjustment Tolerance Auto  
Adjust volume  
Position Isocenter  
Orientation Transversal  
Rotation 0.00 deg  
R >> L 94 mm  
A >> P 24 mm  
F >> H 12 mm

## Physio

1st Signal/Mode None

## Composing

## Sequence

Introduction Off  
Dimension 3D  
Reordering Centric  
Contrasts 1  
Bandwidth 1144 Hz/Px  
Turbo factor 5  
EPI factor 28  
RF pulse type Normal  
Gradient mode Fast  
refocussing type sinc 2560  
flip angle excit 90  
phase encoding ON  
Maxwell compensation Off  
ICE program single  
prepscans 0  
excite duration 2560  
refoc duration 2560  
excite BWTP 16  
refoc BWTP 8  
Opposite Polarity Crusher Off  
pre-crusher 0  
post-crusher1 0

# SIEMENS MAGNETOM Investigational\_Device\_7T syngo MR B17

post-crusher2	0
post-crusher3	0
post-crusher4	0

## Table of contents

\\USER

Feinberglab

Suhung

CS-M1

BP\_grase\_clean\_IV\_Regular\_SH\_ob  
BP\_grase\_clean\_IV\_CS2D\_Par20\_NonSFL\_SH\_Ob  
--time permitting--  
BP\_grase\_clean\_VASO\_V10t\_noClip  
BP\_grase\_clean\_IV\_Regular\_TE43\_SH\_TR3000  
--02/18--  
BP\_grase\_clean\_IV\_Regular\_SH\_dbIOb  
BP\_grase\_clean\_IV\_CS2D\_Par20\_NonSFL\_SH\_dbIOb