

**Table of contents**

\\USER

AMRIT

Alex

UC\_Berkeley\_Comparison

ep2d\_bold\_pat3\_sms3\_1.3x1.3x1.3\_dpg  
ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV3\_dpg  
ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV2\_dpg  
ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_dpg  
ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_68PF\_dpg  
ep2d\_bold\_pat5\_sms3\_0.6x0.6x0.6\_dpg  
ep2d\_bold\_pat4\_sms3\_0.6x0.6x0.6\_68pf\_dpg  
ep2d\_bold\_pat4\_sms2\_0.6x0.6x0.6\_68pf\_dpg  
ep2d\_bold\_pat4\_sms3\_0.5x0.5x0.6\_fov3\_58PF\_dpg

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_1.3x1.3x1.3\_dpg

TA: 1:57 PM: FIX Voxel size: 1.3x1.3x1.3 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	111
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
TR	2160 ms
TE	23.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	2160 ms
TE	23.0 ms
MTC	Off
Flip angle exc	78 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
Base resolution	168
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	111
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.3 mm
TR	2160 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	210 mm
R >> L	210 mm
F >> H	145 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2160 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.73 ms
Bandwidth	1654 Hz/Px

**Sequence - Part 2**

EPI factor	168
RF pulse type	Fast
Gradient mode	Normal

**Sequence - Special**

DPG	On
FLEET	On

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV3\_dpg

TA: 3:10 PM: FIX Voxel size: 0.8x0.8x0.8 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	135
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	3960 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	3960 ms
TE	30.0 ms
MTC	Off
Flip angle exc	87 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
Base resolution	260
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	135
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	3960 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	205 mm
R >> L	208 mm
F >> H	108 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3960 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1202 Hz/Px

**Sequence - Part 2**

EPI factor	256
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat3\_sms3\_0.8x0.8x0.8\_FOV2\_dpg

TA: 3:10 PM: FIX Voxel size: 0.8x0.8x0.8 mmPAT: 9 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	135
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	3960 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	3960 ms
TE	30.0 ms
MTC	Off
Flip angle exc	87 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
Base resolution	260
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	3

**Resolution - iPAT**

Ref. lines PE	66
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	135
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	208 mm
FoV phase	98.5 %
Slice thickness	0.8 mm
TR	3960 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	205 mm
R >> L	208 mm
F >> H	108 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3960 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1 ms
Bandwidth	1202 Hz/Px

**Sequence - Part 2**

EPI factor	256
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_dpg

TA: 4:57 PM: FIX Voxel size: 0.5x0.5x0.6 mmPAT: 16 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	4720 ms
TE	27.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	4720 ms
TE	27.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	88
Accel. factor slice	4
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	4720 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L



**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	104 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	4720 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms4\_0.5x0.5x0.6\_fov3\_68PF\_dpg

TA: 5:40 PM: FIX Voxel size: 0.5x0.5x0.6 mmPAT: 16 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	5390 ms
TE	42.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	5390 ms
TE	42.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	88
Accel. factor slice	4
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	5390 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	104 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5390 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat5\_sms3\_0.6x0.6x0.6\_dpg

TA: 5:11 PM: FIX Voxel size: 0.6x0.6x0.6 mmPAT: 15 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	5180 ms
TE	30.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	5180 ms
TE	30.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	5

**Resolution - iPAT**

Ref. lines PE	70
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	5180 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	103 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	5180 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms3\_0.6x0.6x0.6\_68pf\_dpg

TA: 5:29 PM: FIX Voxel size: 0.6x0.6x0.6 mmPAT: 12 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	6090 ms
TE	35.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	6090 ms
TE	35.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	56
Accel. factor slice	3
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	6090 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	103 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	6090 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms2\_0.6x0.6x0.6\_68pf\_dpg

TA: 6:53 PM: FIX Voxel size: 0.6x0.6x0.6 mmPAT: 8 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	9180 ms
TE	35.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	9180 ms
TE	35.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	334
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	56
Accel. factor slice	2
FOV shift factor	2
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	172
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	9180 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L



**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	104 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	9180 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.22 ms
Bandwidth	936 Hz/Px

**Sequence - Part 2**

EPI factor	334
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off

\\USER\AMRIT\Alex\UC\_Berkeley\_Comparison\ep2d\_bold\_pat4\_sms3\_0.5x0.5x0.6\_fov3\_58PF\_dpg

TA: 5:35 PM: FIX Voxel size: 0.5x0.5x0.6 mmPAT: 12 Rel. SNR: 1.00 : WIPep

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	6200 ms
TE	26.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	A32

**Contrast - Common**

TR	6200 ms
TE	26.0 ms
MTC	Off
Flip angle exc	80 deg
Flip angle fat sat	110 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
Base resolution	382
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off

**Resolution - iPAT**

Accel. mode	Slice accel.
Accel. factor PE	4

**Resolution - iPAT**

Ref. lines PE	88
Accel. factor slice	3
FOV shift factor	3
Reference scan mode	EPI/separate

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	171
Dist. factor	0 %
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.6 mm
TR	6200 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

**Geometry - AutoAlign**

Slice group	1
Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R4.4 A19.1 H0.6
R	4.4 mm
A	19.1 mm
H	0.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-7.5
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L

**System - Miscellaneous**

Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R4.4 A19.1 H0.6 mm
Orientation	T > C-7.5
Rotation	0.00 deg
A >> P	200 mm
R >> L	200 mm
F >> H	103 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	297.210676 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	6200 ms
Concatenations	1

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

**BOLD**

Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	25
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	1.35 ms
Bandwidth	818 Hz/Px

**Sequence - Part 2**

EPI factor	382
RF pulse type	Fast
Gradient mode	Fast*

**Sequence - Special**

DPG	On
FLEET	Off