

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\24ch fMRI\localizer_50V_newcoil

TA: 0:27

PAT: Off

Voxel size: 1.2x1.1x3.0 mm

Rel. SNR: 1.00

SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	L10-24;LV1-9

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Resolution

Base resolution	256
-----------------	-----

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	50.000 V

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\24ch fMRI\b1map_100V_TR1000_RL

TA: 1:09

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 P37.8 F5.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1000 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On

L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	100.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Composing

Sequence

Contrasts	2
Bandwidth	260.416667 Hz/Px
T1 Compensation	Mean T1
Mean T1	1000.0 ms
Angles	1
Amplitude Weighting	Linear
Scale Bar	Enabled
Raw Data	Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\24ch fMRI\lgFactorMap_100V

TA: 3:13

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 P37.8 F5.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	12
Filter	None
Coil elements	L10-24;LV1-9

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	496
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position

H

Table position

0 mm

Inline Composing

Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode

FIX

MSMA

S - C - T

Sagittal

R >> L

Coronal

A >> P

Transversal

F >> H

Save uncombined

Off

Coil Combine Mode

Adaptive Combine

AutoAlign

Auto Coil Select

Default

Shim mode

Standard

Adjust with body coil

Off

Confirm freq. adjustment

Off

Assume Silicone

Off

! Ref. amplitude 1H

100.000 V

Adjustment Tolerance

Auto

Adjust volume

! Position

Isocenter

! Orientation

Transversal

! Rotation

0.00 deg

! R >> L

350 mm

! A >> P

263 mm

! F >> H

350 mm

Physio

1st Signal/Mode

None

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off

MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\24ch fMRI\AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_shimWholeVol

TA: 1:47 PAT: 2 Voxel size: 0.7x0.7x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	R6.0 A6.3 F32.8
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	172 mm
FoV phase	100.0 %
Slice thickness	0.75 mm
TR	3000 ms
TE	24.4 ms
Multi-band accel. factor	2
Filter	None
Coil elements	L1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	230
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	56
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

L1	On
L2	On
L3	On
L4	On
L5	On
L6	On
L7	On
L8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	120.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R6.0 A6.3 F32.8
Orientation	Transversal
Rotation	0.00 deg
R >> L	172 mm
A >> P	172 mm
F >> H	45 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1144 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.05 ms
SIR accel. factor	1
EPI factor	230
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.20
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\flow\fl_fq_mb2f2_gre_3D_seg

TA: 9:52 PAT: 2 Voxel size: 1.5x1.5x2.0 mm Rel. SNR: 1.00 USER: fl_fq_mb_gre_3D_seg

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	2
Dist. factor	100 %
Position	L0.0 A29.7 F58.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	37.65 ms
TE	5.91 ms
Averages	1
Concatenations	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	15 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	Pulse/Trigger
Average cycle	No Signal ms
Captured cycle	-not set-
Acquisition window	770 ms
Trigger pulse	1
Trigger delay	0 ms
Segments	1
Phases	20

Angio

Flow mode	Single vel.
Encodings	3
Velocity enc.	90 cm/s
Direction 1	Through plane
Direction 2	A >> P
Direction 3	R >> L
Rephased images	On
Magnitude images	On
Magnitude sum	Off
Phase images	On
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
MB Number	2
FOV Shift	2

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\flow\fl_fq_mb2f2_gre_3D_seg2

TA: 5:27 PAT: 2 Voxel size: 1.5x1.5x2.0 mm Rel. SNR: 1.00 USER: fl_fq_mb_gre_3D_seg

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	2
Dist. factor	100 %
Position	L0.0 A29.7 F58.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	2.00 mm
TR	75.15 ms
TE	5.91 ms
Averages	1
Concatenations	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	15 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	Separate
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	Pulse/Trigger
Average cycle	No Signal ms
Captured cycle	-not set-
Acquisition window	849 ms
Trigger pulse	1
Trigger delay	0 ms
Segments	2
Phases	11

Angio

Flow mode	Single vel.
Encodings	3
Velocity enc.	90 cm/s
Direction 1	Through plane
Direction 2	A >> P
Direction 3	R >> L
Rephased images	On
Magnitude images	On
Magnitude sum	Off
Phase images	On
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
MB Number	2
FOV Shift	2

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\VASO\BOLD_local_1p5

TA: 7:24 PAT: Off Voxel size: 1.5x1.5x1.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.1 A20.2 H6.1
Orientation	T > C-30.6
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	4000 ms
TE	19.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	110
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	245.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.1 A20.2 H6.1
Orientation	T > C-30.6
Rotation	0.00 deg
R >> L	210 mm
A >> P	210 mm
F >> H	45 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1786 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.66 ms
<hr/>	
SIR accel. factor	1
EPI factor	140
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	1.00
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\VASO\ep2d_fid_VASO-1.5x1.5x1.5_1450_noIPAT

TA: 7:24 PAT: Off Voxel size: 1.5x1.5x1.5 mm Rel. SNR: 1.00 USER: ep2d_fid_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	L0.1 A16.8 H6.7
Orientation	T > C-30.6
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.5 mm
TR	4000.0 ms
TE 1	19 ms
TE 2	19 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	Non-sel. IR
T1	1450 ms
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	110
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	245.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.1 A20.2 H6.1
! Orientation	T > C-30.6
! Rotation	0.00 deg
! R >> L	210 mm
! A >> P	210 mm
! F >> H	45 mm

Physio

1st Signal/Mode	None
-----------------	------

Perf

GBP	Off
PBP	Off
TTP	Off
Original images	On

Sequence

Introduction	Off
Contrasts	2
Bandwidth	1786 Hz/Px
Free echo spacing	Off
Echo spacing	0.66 ms

EPI factor	140
RF pulse type	Fast
Gradient mode	Normal

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\VASO\ep2d_fid_VASO-1.5x1.5x1.5_1350_noIPAT

TA: 7:24 PAT: Off Voxel size: 1.5x1.5x1.5 mm Rel. SNR: 1.00 USER: ep2d_fid_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	L0.1 A16.8 H6.7
Orientation	T > C-30.6
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.5 mm
TR	4000.0 ms
TE 1	19 ms
TE 2	19 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	Non-sel. IR
T1	1350 ms
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	110
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	140
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	245.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.1 A20.2 H6.1
! Orientation	T > C-30.6
! Rotation	0.00 deg
! R >> L	210 mm
! A >> P	210 mm
! F >> H	45 mm

Physio

1st Signal/Mode	None
-----------------	------

Perf

GBP	Off
PBP	Off
TTP	Off
Original images	On

Sequence

Introduction	Off
Contrasts	2
Bandwidth	1786 Hz/Px
Free echo spacing	Off
Echo spacing	0.66 ms
EPI factor	140
RF pulse type	Fast
Gradient mode	Normal

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\VASO\BP_grase_clean_VASO_V07_func_1450

TA: 7:20 PAT: Off Voxel size: 1.5x1.5x1.5 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V07_101320

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R0.4 P38.2 F25.8
Orientation	T > C-30.6
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	192 mm
FoV phase	25.0 %
Slice thickness	1.5 mm
TR	4000 ms
TE	40.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1450 ms
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	110
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	48 mm

Position	L0.0 P39.2 F23.4
Orientation	C > T30.8
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R0.4 P38.2 F25.8
Orientation	T > C-30.6
Rotation	0.00 deg
R >> L	192 mm
A >> P	48 mm
F >> H	12 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1148 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	32
RF pulse type	Normal
Gradient mode	Fast

flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\localizer_50V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B1-8

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	50.000 V
Adjustment Tolerance	Auto
Adjust volume	

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\b1map_100V_TR1000_RL

TA: 1:09

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A12.8 F16.2
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1000 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B1-8

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P75.6 F0.7
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On

B4	On
B5	On
B6	On
B7	On
B8	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	100.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Composing

Sequence

Contrasts	2
Bandwidth	260.416667 Hz/Px
T1 Compensation	Mean T1
Mean T1	1000.0 ms
Angles	1
Amplitude Weighting	Linear
Scale Bar	Enabled
Raw Data	Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\localizer_200V_nova

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
----------	-----

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\b1map_200V_32

TA: 2:10

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1938 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P35.8 F18.2
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On

M3 On
V32 Off

Positioning mode FIX
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Tune up
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 100.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 350 mm
A >> P 263 mm
F >> H 350 mm

Composing

Sequence

Contrasts 2
Bandwidth 260.416667 Hz/Px
T1 Compensation Mean T1
Mean T1 1000.0 ms
Angles 1
Amplitude Weighting Linear
Scale Bar Enabled
Raw Data Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\gFactorMap_32

TA: 3:13

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	12
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	496
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_32ch

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	12
Meas[1]	Active
Meas[2]	Active
Meas[3]	Active
Meas[4]	Active
Meas[5]	Active
Meas[6]	Active
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No

Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_32

TA: 3:21 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L3.4 P45.8 F11.0
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.2 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L3.4 P45.8 F11.0
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_32

TA: 2:15 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_32

TA: 1:55 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_32

TA: 1:24 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.41 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_32

TA: 1:12 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_1mm_tSNR_32

TA: 1:12 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_1mm_tSNR_32

TA: 1:06 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.33 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\localizer_50V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	L0.0 P77.6 H18.2
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	L10-24;LV1-9

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	50.000 V

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\b1map_100V_24

TA: 1:09

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A10.8 F4.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1000 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P54.0 F18.2
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On

LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	100.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Composing

Sequence

Contrasts	2
Bandwidth	260.416667 Hz/Px
T1 Compensation	Mean T1
Mean T1	1000.0 ms
Angles	1
Amplitude Weighting	Linear
Scale Bar	Enabled
Raw Data	Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\gFactorMap_24

TA: 3:13

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A24.3 F8.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	12
Filter	None
Coil elements	L10-24;LV1-9

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	496
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off

MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px

Gradient mode	Fast
RF spoiling	On

ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_24

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L0.1 P67.4 F10.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.1 P67.4 F10.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_24

TA: 3:21 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L3.4 P67.4 F12.4
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.2 ms
Multi-band accel. factor	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L3.4 P67.4 F12.4
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_24

TA: 2:15 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	R1.3 P83.7 H0.0
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R1.3 P83.7 H0.0
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_24

TA: 1:55 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	R1.3 P83.7 H0.0
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R1.3 P83.7 H0.0
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_24

TA: 1:24 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P80.9 H0.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P80.9 H0.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.41 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_24

TA: 1:12 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P80.9 H0.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P80.9 H0.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20

Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat1mb2_1mm_tSNR_24

TA: 1:12 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties		Inline Composing	Off
Prio Recon		System	
Before measurement		LV1	On
After measurement		LV2	On
Load to viewer		LV3	On
Inline movie		LV4	On
Auto store images		LV5	On
Load to stamp segments		LV6	On
Load images to graphic segments		LV7	On
Auto open inline display		LV8	On
Start measurement without further preparation		LV9	On
Wait for user to start		L10	On
Start measurements		L11	On
		L12	On
		L13	On
		L14	On
		L15	On
		L16	On
		L17	On
		L18	On
		L19	On
		L20	On
		L21	On
		L22	On
		L23	On
		L24	On
Routine		Positioning mode	FIX
Slice group 1		MSMA	S - C - T
Slices		Sagittal	R >> L
Dist. factor		Coronal	A >> P
Position		Transversal	F >> H
Orientation		Coil Combine Mode	Sum of Squares
Phase enc. dir.		AutoAlign	---
Rotation		Auto Coil Select	Default
Phase oversampling		Shim mode	Standard
FoV read		Adjust with body coil	Off
FoV phase		Confirm freq. adjustment	On
Slice thickness		Assume Silicone	Off
TR		! Ref. amplitude 1H	0.000 V
TE		Adjustment Tolerance	Auto
Multi-band accel. factor		Adjust volume	
Filter		Position	L0.0 P80.2 H3.6
Coil elements		Orientation	Coronal
		Rotation	90.00 deg
		R >> L	90 mm
		F >> H	72 mm
		A >> P	30 mm
Contrast		Physio	
MTC		1st Signal/Mode	None
Magn. preparation		BOLD	
Flip angle		GLM Statistics	Off
Fat suppr.		Dynamic t-maps	Off
		Starting ignore meas	0
Averaging mode		Ignore after transition	0
Reconstruction		Model transition states	On
Measurements		Temp. highpass filter	On
Delay in TR		Threshold	4.00
Multiple series		Paradigm size	20
		Meas[1]	Baseline
		Meas[2]	Baseline
		Meas[3]	Baseline
Resolution			
Base resolution			
Phase resolution			
Phase partial Fourier			
Interpolation			
PAT mode			
Distortion Corr.			
Prescan Normalize			
Raw filter			
Elliptical filter			
Hamming			
Geometry			
Multi-slice mode			
Series			
Special sat.			
Table position			
Table position			

Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat1mb1_1mm_tSNR_24

TA: 1:06 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P80.2 H3.6
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	1
Filter	None
Coil elements	L10-24;LV1-9

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

LV1	On
LV2	On
LV3	On
LV4	On
LV5	On
LV6	On
LV7	On
LV8	On
LV9	On
L10	On
L11	On
L12	On
L13	On
L14	On
L15	On
L16	On
L17	On
L18	On
L19	On
L20	On
L21	On
L22	On
L23	On
L24	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P80.2 H3.6
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline

Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.33 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\localizer_50V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	L0.0 A8.8 F17.5
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	L0.0 P68.8 H0.7
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	L0.0 P68.8 H0.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B1-8

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	50.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\b1map_100V_8

TA: 1:09

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A12.1 F2.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1000 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B1-8

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L2.7 P68.1 H0.7
Orientation	Transversal
Rotation	0.00 deg
Base size phase	23 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On

B4	On
B5	On
B6	On
B7	On
B8	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	100.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Composing

Sequence

Contrasts	2
Bandwidth	260.416667 Hz/Px
T1 Compensation	Mean T1
Mean T1	1000.0 ms
Angles	1
Amplitude Weighting	Linear
Scale Bar	Enabled
Raw Data	Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\gFactorMap_8

TA: 3:13

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A12.1 F2.0
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	12
Filter	None
Coil elements	B1-8

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	496
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Rx coil diode switching	On
coil channel reordering	Off
<hr/>	
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_8

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L8.8 P72.2 F11.7
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L8.8 P72.2 F11.7
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_8

TA: 3:21 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L3.4 P72.2 F9.0
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.2 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L3.4 P72.2 F9.0
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_8

TA: 2:15 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	R1.3 P88.1 F1.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R1.3 P88.1 F1.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_8

TA: 1:55 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	R1.3 P88.1 F1.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R1.3 P88.1 F1.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_8

TA: 1:24 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P85.6 F2.1
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P85.6 F2.1
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.41 ms
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_8

TA: 1:12 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P85.6 F2.1
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P85.6 F2.1
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active

Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat1mb2_1mm_tSNR_8

TA: 1:12 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P84.8 H1.1
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P84.8 H1.1
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilTest2018\AV_ep2d_bold_sd1ipat1mb1_1mm_tSNR_8

TA: 1:06 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P84.8 H1.1
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P84.8 H1.1
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.33 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\localizer_200V_nova

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
----------	-----

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\b1map_200V_32

TA: 2:10

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1938 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P35.8 F18.2
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On

M3 On
V32 Off

Positioning mode FIX
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Tune up
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 100.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 350 mm
A >> P 263 mm
F >> H 350 mm

Composing

Sequence

Contrasts 2
Bandwidth 260.416667 Hz/Px
T1 Compensation Mean T1
Mean T1 1000.0 ms
Angles 1
Amplitude Weighting Linear
Scale Bar Enabled
Raw Data Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\gFactorMap_32

TA: 3:13

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	12
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	496
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	Isocenter
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	350 mm
! A >> P	263 mm
! F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\tSNR_32_sd1ipat2mb2_pt5mm_AV_ep2d_bold

TA: 2:15 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\tSNR_32_sd1ipat2mb1_pt5mm_AV_ep2d_bold

TA: 1:55 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	0 %
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	5000 ms
TE	26.0 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P66.1 F12.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.43 ms
<hr/>	
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\tsNR_32_sd1ipat2mb2_pt75mm_AV_ep2d_bold

TA: 1:24 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.41 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\tsNR_32_sd1ipat2mb1_pt75mm_AV_ep2d_bold

TA: 1:12 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.8 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P64.1 F11.5
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\tsnr_32_sd1ipat2mb2_1mm_AV_ep2d_bold

TA: 1:12 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\SNR_32_sd1ipat2mb1_1mm_AV_ep2d_bold

TA: 1:06 PAT: Off Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	80.0 %
Slice thickness	1.00 mm
TR	3000 ms
TE	22.6 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	90
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P65.4 F10.8
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	72 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.33 ms
<hr/>	
SIR accel. factor	1
EPI factor	72
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	12
Meas[1]	Active
Meas[2]	Active
Meas[3]	Active
Meas[4]	Active
Meas[5]	Active
Meas[6]	Active
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No

Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	12
Meas[1]	Active
Meas[2]	Active
Meas[3]	Active
Meas[4]	Active
Meas[5]	Active
Meas[6]	Active
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No

Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold

TA: 3:33 PAT: 2 Voxel size: 0.5x0.5x0.5 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	60
Dist. factor	0 %
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.9 %
Slice thickness	0.50 mm
TR	3000 ms
TE	26.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	180
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 P45.2 F6.3
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	12
Meas[1]	Active
Meas[2]	Active
Meas[3]	Active
Meas[4]	Active
Meas[5]	Active
Meas[6]	Active
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No

Free echo spacing	Off
Echo spacing	1.43 ms
SIR accel. factor	1
EPI factor	160
Gradient mode	Normal
RF spoiling	Off
Excite pulse duration	5820 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\CoilComparison_JT\visLoc_32_sd1ipat2mb1_pt75mm_AV_ep2d_bold

TA: 3:21 PAT: 2 Voxel size: 0.7x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L3.4 P45.8 F11.0
Orientation	Coronal
Phase enc. dir.	F >> H
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	90 mm
FoV phase	88.3 %
Slice thickness	0.75 mm
TR	3000 ms
TE	23.2 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	63
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	120
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L3.4 P45.8 F11.0
Orientation	Coronal
Rotation	90.00 deg
R >> L	90 mm
F >> H	80 mm
A >> P	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	816 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.37 ms
<hr/>	
SIR accel. factor	1
EPI factor	106
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	5820 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\localizer_200V_nova

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	500 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
----------	-----

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\b1map_200V_32

TA: 2:10

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1938 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P35.8 F18.2
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On

M3 On
V32 Off

Positioning mode FIX
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Tune up
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 200.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 350 mm
A >> P 263 mm
F >> H 350 mm

Composing

Sequence

Contrasts 2
Bandwidth 260.416667 Hz/Px
T1 Compensation Mean T1
Mean T1 1000.0 ms
Angles 1
Amplitude Weighting Linear
Scale Bar Enabled
Raw Data Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\mp2rage_1mm_TR4000

TA: 5:26 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: mp2rage_wip602B

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L1.9 A29.6 F31.5
Orientation	Sagittal
Phase enc. dir.	H >> F
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	11.1 %
Slices per slab	144
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	1.00 mm
TR	4000 ms
TE	3.23 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	240.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 A24.9 F9.3
! Orientation	Sagittal
! Rotation	0.00 deg
! F >> H	108 mm
! A >> P	160 mm
! R >> L	127 mm

Contrast

Magn. preparation	Non-sel. IR
T1 1	1000 ms
T1 2	3200 ms
Flip angle 1	4 deg
Flip angle 2	4 deg
Fat suppr.	Water excit. fast
Water suppr.	None
2nd Inversion-Contrast	On
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Resolution

Base resolution	192
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	200 Hz/Px
Flow comp.	Slice
Echo spacing	7.8 ms
<hr/>	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
<hr/>	
FFT Scale Factor	200 %
Line/Partition Swap	Off
Homodyne Phase Filter	Off
Flat Image	On
T1 Map	On
Division Image	Off
ExtInvPulseOn	On
OffResFreqInv	0
Invflipangle	970

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_pt8mm_MB2IPAT3_wGap_pf5

TA: 1:16

PAT: 3

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.80 mm
TR	2400 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	1200 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1220 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.03 ms
<hr/>	
SIR accel. factor	1
EPI factor	232
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_pt8mm_MB2IPAT3_wGap_pf6

TA: 1:15 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	28
Dist. factor	0 %
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.80 mm
TR	2400 ms
TE	28.2 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	1200 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	23 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1220 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.03 ms
<hr/>	
SIR accel. factor	1
EPI factor	232
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_pt8mm_MB2IPAT4_wGap_pf7

TA: 1:16 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	24
Dist. factor	0 %
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.80 mm
TR	2400 ms
TE	31.6 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	1200 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	64
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1220 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.04 ms
<hr/>	
SIR accel. factor	1
EPI factor	232
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_pt8mm_MB2IPAT4_wGap_pf7_altFOV

TA: 1:16 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	24
Dist. factor	0 %
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	153 mm
FoV phase	100.0 %
Slice thickness	0.80 mm
TR	2400 ms
TE	28.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	1200 ms
Multiple series	Off

Resolution

Base resolution	192
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	64
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L2.1 A33.1 F11.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	153 mm
A >> P	153 mm
F >> H	20 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1184 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.07 ms
<hr/>	
SIR accel. factor	1
EPI factor	192
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GRASE_pt8mm_wGap

TA: 6.0 s PAT: Off Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: pgrs3d_ey_20180816

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	16
FoV read	200 mm
FoV phase	11.2 %
Slice thickness	0.80 mm
TR	2000 ms
TE	35.66 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	180 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Multiple series	Off

Resolution

Base resolution	250
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Raw filter	Off

Geometry

Series	Ascending
Sat. region 1	
Thickness	24 mm
Position	Isocenter
Orientation	Coronal

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	100.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	23 mm
F >> H	13 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1112 Hz/Px
Echo spacing	1.06 ms

Slice turbo factor	10
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast

Adjust flipangles	Off
FLIP ANGLES[1]	180 deg
FLIP ANGLES[2]	180 deg
FLIP ANGLES[3]	180 deg
FLIP ANGLES[4]	180 deg
FLIP ANGLES[5]	180 deg
FLIP ANGLES[6]	180 deg
FLIP ANGLES[7]	180 deg
FLIP ANGLES[8]	180 deg
FLIP ANGLES[9]	180 deg
FLIP ANGLES[10]	180 deg

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Crusher Factor	7.00
Spoiler Factor	1.00
RF02 BWT Factor	1.00
RF02 time	2560 ms
RF Scaling[1]	1.00
Phase Encoding PE	On
Phase Encoding 3D	On
Measurement Number	1023
Inversion Flag	Off
FFT Scale Factor	0.10
ACROSSSEGMENTS	Off
PRIMARYMODE	On
AUTOCORR	Off
CROSSCORR	Off
FILTERED	On
FatSat FlipAngle	110 deg

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1pt5mm_MB5IPAT3_500ms

TA: 0:31

PAT: 3

Voxel size: 1.5x1.5x1.5 mm

Rel. SNR: 1.00

USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	500 ms
TE	18.8 ms
Multi-band accel. factor	5
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	40 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	60 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1446 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	128
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	5 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1pt5mm_MB1IPAT3_2000ms

TA: 1:03

PAT: 3

Voxel size: 1.5x1.5x1.5 mm

Rel. SNR: 1.00

USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	2000 ms
TE	18.8 ms
Multi-band accel. factor	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	60 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1446 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	0.8 ms
<hr/>	
SIR accel. factor	1
EPI factor	128
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Fake MB factor for SB	1
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1pt5mm_MB5IPAT3_2000ms

TA: 1:13

PAT: 3

Voxel size: 1.5x1.5x1.5 mm

Rel. SNR: 1.00

USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	40
Dist. factor	0 %
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	1.50 mm
TR	2000 ms
TE	18.8 ms
Multi-band accel. factor	5
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.3 A29.0 H0.7
Orientation	Transversal
Rotation	0.00 deg
R >> L	192 mm
A >> P	192 mm
F >> H	60 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1446 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	128
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	5 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT4_pf7

TA: 0:17 PAT: 4 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	24.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	3
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	64
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf5_420i_p+f

TA: 7:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	420
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf5_612i_p+f+s

TA: 10:29 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	612
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf5_360i_TIMITb5

TA: 6:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf5_360i_TIMITb6

TA: 6:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_p8mm_MB2IPAT3_pf5_360i_TIMITb5

TA: 6:18 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160 mm
FoV phase	91.0 %
Slice thickness	0.80 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	160 mm
A >> P	146 mm
F >> H	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.01 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\mp2rage_0.7mm_TR4500

TA: 7:36 PAT: 3 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: mp2rage_wip602B

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved
<hr/>	
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L1.2 A28.3 F30.2
Orientation	Sagittal
Phase enc. dir.	H >> F
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	8.3 %
Slices per slab	192
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.80 mm
TR	4500 ms
TE	3.33 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

System

T1	On
M2	On
B4	On
M3	On
V32	Off
<hr/>	
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
<hr/>	
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 A24.9 F9.3
! Orientation	Sagittal
! Rotation	0.00 deg
! F >> H	108 mm
! A >> P	160 mm
! R >> L	127 mm

Contrast

Magn. preparation	Non-sel. IR
T1 1	1000 ms
T1 2	3200 ms
Flip angle 1	4 deg
Flip angle 2	4 deg
Fat suppr.	Water excit. fast
Water suppr.	None
2nd Inversion-Contrast	On
<hr/>	
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Physio

1st Signal/Mode	None
<hr/>	
Dark blood	Off
<hr/>	
Resp. control	Off

Resolution

Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	6/8
Interpolation	Off
<hr/>	
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	200 Hz/Px
Flow comp.	Slice
Echo spacing	8.1 ms
<hr/>	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
<hr/>	
FFT Scale Factor	200 %
Line/Partition Swap	Off
Homodyne Phase Filter	Off
Flat Image	On
T1 Map	On
Division Image	Off
ExtInvPulseOn	On
OffResFreqInv	0
Invflipangle	970

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf5_420i_p+f

TA: 0:40 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT3_pf6_420i_p+f

TA: 0:40 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	23.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_PotentialScans\GE_1mm_MB2IPAT4_pf7_420i_p+f

TA: 0:37 PAT: 4 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	24.0 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	23
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	64
Reference scan mode	Segmented
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\BWTP\ep3d_epi_backup_CAnew

TA: 2.2 s PAT: Off Voxel size: 7.8x7.8x5.0 mm Rel. SNR: 1.00 USER: ep3d_epi_backup_CAnew

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	10 ms
TE	1.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	500 mm
A >> P	500 mm
F >> H	5 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	On
Interpolation	3D-K-space

Spatial filter	Off
----------------	-----

Sequence

Introduction	Off
Dimension	2D
Bandwidth	752 Hz/Px
Free echo spacing	Off
Echo spacing	1.4 ms
<hr/>	
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast
<hr/>	
RF90 duration	5120
MB Number	1
DummyScan Number	1
TR Extent(us)	0
BwTimeProd	52
PhaseOffset	0
FOV Shift	0
Polarity	0
Interleaved	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\localizer_200V_nova

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	500 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
----------	-----

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\b1map_200V_32

TA: 2:10

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	100 %
Position	R0.7 A30.3 F0.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	1938 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	L0.0 P35.8 F18.2
Orientation	Transversal
Rotation	0.00 deg
Base size phase	50 mm
Base size read	50 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On

M3 On
V32 Off

Positioning mode FIX
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Tune up
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 200.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 350 mm
A >> P 263 mm
F >> H 350 mm

Composing

Sequence

Contrasts 2
Bandwidth 260.416667 Hz/Px
T1 Compensation Mean T1
Mean T1 1000.0 ms
Angles 1
Amplitude Weighting Linear
Scale Bar Enabled
Raw Data Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\mp2rage_0.7mm_TR4500

TA: 7:36

PAT: 3

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: mp2rage_wip602B

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved
<hr/>	
Table position	H
Table position	0 mm
Inline Composing	Off

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L1.2 A28.3 F30.2
Orientation	Sagittal
Phase enc. dir.	H >> F
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	8.3 %
Slices per slab	192
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.80 mm
TR	4500 ms
TE	3.33 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

System

T1	On
M2	On
B4	On
M3	On
V32	Off
<hr/>	
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
<hr/>	
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 A24.9 F9.3
! Orientation	Sagittal
! Rotation	0.00 deg
! F >> H	108 mm
! A >> P	160 mm
! R >> L	127 mm

Contrast

Magn. preparation	Non-sel. IR
T1 1	1000 ms
T1 2	3200 ms
Flip angle 1	4 deg
Flip angle 2	4 deg
Fat suppr.	Water excit. fast
Water suppr.	None
2nd Inversion-Contrast	On
<hr/>	
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Physio

1st Signal/Mode	None
<hr/>	
Dark blood	Off
<hr/>	
Resp. control	Off

Resolution

Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	6/8
Interpolation	Off
<hr/>	
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	200 Hz/Px
Flow comp.	Slice
Echo spacing	8.1 ms
<hr/>	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
<hr/>	
FFT Scale Factor	200 %
Line/Partition Swap	Off
Homodyne Phase Filter	Off
Flat Image	On
T1 Map	On
Division Image	Off
ExtInvPulseOn	On
OffResFreqInv	0
Invflipangle	970

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\GE_1mm_MB2IPAT3_pf5_420i_p+f

TA: 7:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	420
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\GE_1mm_MB2IPAT3_pf5_612i_p+f+s

TA: 10:29 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	612
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\GE_1mm_MB2IPAT3_pf5_360i_TIMITb5

TA: 6:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\GE_1mm_MB2IPAT3_pf5_360i_TIMITb6

TA: 6:17 PAT: 3 Voxel size: 1.0x1.0x1.0 mm Rel. SNR: 1.00 USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	91.0 %
Slice thickness	1.00 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	182 mm
F >> H	30 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\GE_p8mm_MB2IPAT3_pf5_360i_TIMITb5

TA: 6:18

PAT: 3

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: AV_ep2d_bold_sd_20140727

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	30
Dist. factor	0 %
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160 mm
FoV phase	91.0 %
Slice thickness	0.80 mm
TR	1000 ms
TE	22.8 ms
Multi-band accel. factor	2
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Magn. preparation	None
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	360
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	200
Phase resolution	100 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	48
Reference scan mode	GRE
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved

Special sat.

None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	On
Assume Silicone	Off
! Ref. amplitude 1H	170.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.4 A31.7 F4.9
Orientation	Transversal
Rotation	0.00 deg
R >> L	160 mm
A >> P	146 mm
F >> H	25 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1250 Hz/Px
Flow comp.	No
Free echo spacing	Off
Echo spacing	1.01 ms
<hr/>	
SIR accel. factor	1
EPI factor	182
Gradient mode	Normal
RF spoiling	Off
<hr/>	
Excite pulse duration	3640 us
Slice multiplier	1
Multi-band PE shift	0 1/FoV
zBlip scheme	0
MB kernel size	0
MB knockout band	0
No. of interleaved TEs	0
RF pulse shape	1
EPI noise scans	0
EPI full reference scan	0
Single-band images	On
MB RF phase scramble	Off
SENSE1 coil combine	Off
Log physiology to file	Off
Invert RO/PE polarity	Off
Save reduced raw data	Off
Readout slice trace	Off
Disable ramp sampling	Off
PF omits higher k-space	Off
Online multi-band recon.	Online
FFT scale factor	0.02
GRE iPAT ref. FA	12.0 deg
Send B1 shim trigger	Never
Triggering scheme	Standard
Starting ignore meas	0
Paradigm size	2
Multiplier	1
Step [1]	1
Step [2]	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\VASO_116_phantom_420i_p+f
TA: 10:30 PAT: 2 Voxel size: 0.7x0.7x1.5 mm Rel. SNR: 1.00 USER: VASO_116

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R37.8 A6.9 H16.3
Orientation	T > S25.8 > C-11.1
Phase enc. dir.	R >> L
Rotation	70.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	32.8 mm
FoV phase	300.0 %
Slice thickness	1.50 mm
TR	1500.00 ms
TE	24 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Perfusion mode	Picore Q2TIPS
T12	900 ms
T11	50 ms
T11s	50 ms
Flip angle	26 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	420
Delay in TR	0 ms
Multiple series	Off
Perfusion mode	PICORE Q2T
Inversion time 1	50 ms
Saturation stop time	50 ms
Inversion time 2	900.0 ms
Flow limit	100.0 cm/s

Resolution

Base resolution	44
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	36
Accel. factor 3D	1
Ref. lines 3D	8
Reference scan mode	Separate

Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	Parallel F
Gap	25.0 mm
Thickness	100 mm

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R37.8 A6.9 H16.3
Orientation	T > S25.8 > C-11.1
Rotation	160.00 deg
R >> L	99 mm
A >> P	33 mm
F >> H	12 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

Motion correction	Off
Spatial filter	Off

Sequence

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	1
Bandwidth	1062 Hz/Px
Free echo spacing	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Echo spacing	1.07 ms
EPI factor	132
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Ampl	120
BWDTH	150 3.1kHz
ph.skip 4 Robert (the one)	30
use Ernst angle	Off
Maxwell Correction	Off
log physio files	Off
FFT scale	1.00
dummy prepscan time	3 s
z shim	0.00 mT/m*ms
RF duration	2560 us
RF BWTP	25.0
Renzo: Delta TI	71 ms
EFFECTIVE TR	12000 ms
PatPartitions	8
EPI phase correction	local
PAT refscan mode	Flash
FlashRef BaseRes	44
FlashRef BW	136 Hz/px
FlashRef TE	5000 us
FlashRef FA	5 deg
use CAIPI	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\Auditory_Pilot1\mp2rage_1mm_TR4500

TA: 4:16

PAT: 3

Voxel size: 0.8x0.8x1.0 mm

Rel. SNR: 1.00

USER: mp2rage_wip602B

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R3.0 A3.8 F29.8
Orientation	T > C34.8
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	11.1 %
Slices per slab	72
FoV read	154 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	3530 ms
TE	3.33 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
TI 1	1000 ms
TI 2	3200 ms
Flip angle 1	4 deg
Flip angle 2	4 deg
Fat suppr.	Water excit. fast
Water suppr.	None
2nd Inversion-Contrast	On
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	192
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	36
Accel. factor 3D	1
Reference scan mode	Integrated

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	230.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R3.0 A3.8 F29.8
Orientation	T > C34.8
Rotation	90.00 deg
A >> P	154 mm
R >> L	154 mm
F >> H	72 mm

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	200 Hz/Px
Flow comp.	Slice
Echo spacing	8.1 ms
<hr/>	
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
<hr/>	
FFT Scale Factor	200 %
Line/Partition Swap	Off
Homodyne Phase Filter	Off
Flat Image	On
T1 Map	On
Division Image	Off
ExtInvPulseOn	On
OffResFreqInv	0
Invflipangle	970

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\GV-0319\BP_grase_clean_VASO_V10t_noClip

TA: 12:03 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO_V10_noClip

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	99 mm
FoV phase	25.8 %
Slice thickness	1.5 mm
TR	3000 ms
TE	45.9 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1100 ms
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	241
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	132
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm

Position	L4.6 P8.1 H13.9
Orientation	C > S24.8 > T10.3
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Rotation	-20.00 deg
R >> L	99 mm
A >> P	26 mm
F >> H	12 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	2
Bandwidth	1052 Hz/Px
Echo spacing	1.1 ms
Turbo factor	5
EPI factor	34
RF pulse type	Normal
Gradient mode	Fast
BIR4: 2nd segm phase	338
BIR4: duration	5120
excite duration	2560
refoc duration	2560
excite BWTP	10.4
refoc BWTP	5.2
phase encoding	ON
Maxwell compensation	Off
ICE program	single

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\GV-0319\BP_grase_clean_IV_FA180_Regular_SH

TA: 12:03

PAT: Off

Voxel size: 0.8x0.8x1.5 mm

Rel. SNR: 1.00

USER: BP_grase_clean_IV_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	99 mm
FoV phase	25.8 %
Slice thickness	1.5 mm
TR	3000 ms
TE	39.32 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	180 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	241
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	132
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm
Position	L4.6 P8.1 H13.9

Orientation C > S24.8 > T10.3
Special sat. None

Table position H
Table position 0 mm
Inline Composing Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Rotation	-20.00 deg
R >> L	99 mm
A >> P	26 mm
F >> H	12 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1052 Hz/Px
Echo spacing	1.1 ms
Turbo factor	5
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0
excite duration	2560
refoc duration	2560
excite BWTP	16.0
refoc BWTP	8.0
Opposite Polarity Crusher	Off
pre-crusher	40000

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

post-crusher1	40000
post-crusher2	40000
post-crusher3	40000
post-crusher4	40000

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Alex\GV-0319\BP_grase_clean_IV_FA165_nSTE_CurrentPT_SH

TA: 12:03 PAT: Off Voxel size: 0.8x0.8x1.5 mm Rel. SNR: 1.00 USER: BP_grase_clean_IV_SH

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Phase enc. dir.	A >> P
Rotation	-20.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	8
FoV read	99 mm
FoV phase	25.8 %
Slice thickness	1.5 mm
TR	1500 ms
TE	36.52 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	None
Flip angle	165 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	482
Pause after meas.	0.0 s
Multiple series	Off

Resolution

Base resolution	132
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	5/8
Interpolation	Off
PAT mode	None
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	26 mm
Position	L4.6 P8.1 H13.9

Orientation C > S24.8 > T10.3
Special sat. NoneTable position H
Table position 0 mm
Inline Composing Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	220.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R40.0 A9.4 H29.7
Orientation	T > S40.5 > C-3.3
Rotation	-20.00 deg
R >> L	99 mm
A >> P	26 mm
F >> H	12 mm

Physio

1st Signal/Mode None

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	1184 Hz/Px
Echo spacing	1 ms

Turbo factor	5
EPI factor	28
RF pulse type	Normal
Gradient mode	Fast

refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0
excite duration	2560
refoc duration	2560
excite BWTP	16.0
refoc BWTP	8.0
Opposite Polarity Crusher	Off
pre-crusher	40000

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

post-crusher1	20000
post-crusher2	40000
post-crusher3	20000
post-crusher4	40000

Table of contents

\\USER

Feinberglab
Alex

24ch fMRI

localizer_50V_newcoil
 b1map_100V_TR1000_RL
 gFactorMap_100V
 AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_shimWholeVol

flow

fl_fq_mb2f2_gre_3D_seg
 fl_fq_mb2f2_gre_3D_seg2
 svb_mslaser_uhf_15uT

VASO

BOLD_local_1p5
 ep2d_fid_VASO-1.5x1.5x1.5_1450_noIPAT
 ep2d_fid_VASO-1.5x1.5x1.5_1350_noIPAT
 BP_grase_clean_VASO_V07_func_1450

CoilTest2018

localizer_50V_newcoil
 b1map_100V_TR1000_RL
 ---32---
 localizer_200V_nova
 b1map_200V_32
 gFactorMap_32
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_32ch
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_32
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb2_1mm_tSNR_32
 AV_ep2d_bold_sd1ipat2mb1_1mm_tSNR_32
 ---24---
 localizer_50V_newcoil
 b1map_100V_24
 gFactorMap_24
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_24
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_24
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_24
 AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_24
 AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_24
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_24
 AV_ep2d_bold_sd1ipat1mb2_1mm_tSNR_24
 AV_ep2d_bold_sd1ipat1mb1_1mm_tSNR_24
 ---8---
 localizer_50V_newcoil
 b1map_100V_8
 gFactorMap_8
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_visLoc_8
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_visLoc_8
 AV_ep2d_bold_sd1ipat2mb2_pt5mm_tSNR_8
 AV_ep2d_bold_sd1ipat2mb1_pt5mm_tSNR_8
 AV_ep2d_bold_sd1ipat2mb2_pt75mm_tSNR_8
 AV_ep2d_bold_sd1ipat2mb1_pt75mm_tSNR_8
 AV_ep2d_bold_sd1ipat1mb2_1mm_tSNR_8
 AV_ep2d_bold_sd1ipat1mb1_1mm_tSNR_8

CoilComparison_JT

---32---
 localizer_200V_nova
 b1map_200V_32
 gFactorMap_32
 tSNR_32_sd1ipat2mb2_pt5mm_AV_ep2d_bold
 tSNR_32_sd1ipat2mb1_pt5mm_AV_ep2d_bold
 tSNR_32_sd1ipat2mb2_pt75mm_AV_ep2d_bold

\\USER

Feinberglab
Alex

CoilComparison_JT

tSNR_32_sd1ipat2mb1_pt75mm_AV_ep2d_bold
 tSNR_32_sd1ipat2mb2_1mm_AV_ep2d_bold
 tSNR_32_sd1ipat2mb1_1mm_AV_ep2d_bold
 visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold
 visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold
 visLoc_32ch_sd1ipat2mb2_pt5mm_AV_ep2d_bold
 visLoc_32_sd1ipat2mb1_pt75mm_AV_ep2d_bold

Auditory_PotentialScans

localizer_200V_nova
 b1map_200V_32
 mp2rage_1mm_TR4000
 -- best match to Moerel 2018 --
 -- longer TE and/or PF necessary --
 -- due to slower gradients --
 GE_pt8mm_MB2IPAT3_wGap_pf5
 GE_pt8mm_MB2IPAT3_wGap_pf6
 GE_pt8mm_MB2IPAT4_wGap_pf7
 GE_pt8mm_MB2IPAT4_wGap_pf7_altFOV
 GRASE_pt8mm_wGap
 -- fast TR comparisons at 1p5 --
 -- cf Denison 2014 and Vu 2016 --
 -- tSNR comparisons --
 GE_1pt5mm_MB5IPAT3_500ms
 GE_1pt5mm_MB1IPAT3_2000ms
 GE_1pt5mm_MB5IPAT3_2000ms
 -- sub Second 1mm --
 GE_1mm_MB2IPAT4_pf7
 --for piloting--
 GE_1mm_MB2IPAT3_pf5_420i_p+f
 GE_1mm_MB2IPAT3_pf5_612i_p+f+s
 GE_1mm_MB2IPAT3_pf5_360i_TIMITb5
 GE_1mm_MB2IPAT3_pf5_360i_TIMITb6
 GE_p8mm_MB2IPAT3_pf5_360i_TIMITb5
 mp2rage_0.7mm_TR4500
 --tsnr--
 GE_1mm_MB2IPAT3_pf5_420i_p+f
 GE_1mm_MB2IPAT3_pf6_420i_p+f
 GE_1mm_MB2IPAT4_pf7_420i_p+f

BWTP

ep3d_epi_backup_CAnew

Auditory_Pilot1

localizer_200V_nova
 b1map_200V_32
 mp2rage_0.7mm_TR4500
 GE_1mm_MB2IPAT3_pf5_420i_p+f
 GE_1mm_MB2IPAT3_pf5_612i_p+f+s
 GE_1mm_MB2IPAT3_pf5_360i_TIMITb5
 GE_1mm_MB2IPAT3_pf5_360i_TIMITb6
 GE_p8mm_MB2IPAT3_pf5_360i_TIMITb5
 Pause
 VASO_116_phantom_420i_p+f
 mp2rage_1mm_TR4500

GV-0319

--check FOV is 99x25.8% and Matrix is 132x34 & b/w is 1052--
 BP_grase_clean_VASO_V10t_noClip
 BP_grase_clean_IV_FA180_Regular_SH
 BP_grase_clean_IV_FA165_nSTE_CurrentPT_SH
 BP_grase_clean_IV_TE27_Par20_STE_SH