

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\iceTesting\BP_grase_clean_sat_2mb_m1m1

TA: 0:20 PAT: Off Voxel size: 1.6x1.6x1.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_sat_2mb_new

Properties

| | |
|-----------------------------------------------|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|-----------------|
| Slab group 1 | |
| Slabs | 1 |
| Dist. factor | 200 % |
| Position | L0.0 P23.7 H6.8 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| Slice oversampling | 25.0 % |
| Slices per slab | 8 |
| FoV read | 160 mm |
| FoV phase | 32.0 % |
| Slice thickness | 1.0 mm |
| TR | 4000 ms |
| TE | 34.32 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|---------------------|-----------|
| Flip angle | 160 deg |
| Fat suppr. | Fat sat. |
| Fat sat. mode | Strong |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 5 |
| Pause after meas. 1 | 0.0 s |
| Pause after meas. 2 | 0.0 s |
| Pause after meas. 3 | 0.0 s |
| Pause after meas. 4 | 0.0 s |
| Multiple series | Off |

Resolution

| | |
|-----------------------|-------|
| Base resolution | 100 |
| Phase resolution | 100 % |
| Slice resolution | 100 % |
| Slice partial Fourier | Off |
| Interpolation | Off |
| PAT mode | None |
| Raw filter | Off |

Geometry

| | |
|---------------|-----------|
| Series | Ascending |
| Sat. region 1 | |
| Thickness | 50 mm |

| | |
|--------------|-----------------|
| Position | L0.0 P23.7 H0.0 |
| Orientation | Coronal |
| Special sat. | None |

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Save uncombined | Off |
| Coil Combine Mode | Adaptive Combine |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Standard |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | L0.0 P23.7 H6.8 |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 160 mm |
| A >> P | 52 mm |
| F >> H | 8 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
|-----------------|------|

Composing

Sequence

| | |
|----------------------|------------|
| Introduction | Off |
| Dimension | 3D |
| Reordering | Centric |
| Contrasts | 1 |
| Bandwidth | 1924 Hz/Px |
| Echo spacing | 0.7 ms |
| Turbo factor | 10 |
| EPI factor | 32 |
| RF pulse type | Normal |
| Gradient mode | Fast* |
| refocussing type | sinc 2560 |
| flip angle excit | 90 |
| Crusher Momentum | 40000 |
| Crusher Time | 2000 |
| phase encoding | ON |
| Maxwell compensation | Off |
| ICE program | single |
| prepsans | 0 |
| MB Number | 1 |
| Shift or Not | 0 |
| Dummy Number | 1 |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-----------------|------|
| MB Number2 | 1 |
| Dummy TR(ut:ms) | 2000 |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\iceTesting\ep2d_bold_m4p2

TA: 0:33 PAT: 2 Voxel size: 1.6x1.6x1.0 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asy

Properties

| | |
|-----------------------------------------------|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|------------------|
| Slice group 1 | |
| Slices | 96 |
| Dist. factor | 0 % |
| Position | L0.0 A27.8 H28.5 |
| Orientation | T > C-15.3 |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 200 mm |
| FoV phase | 100.0 % |
| Slice thickness | 1.0 mm |
| TR | 1500 ms |
| TE | 19 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-----------------|-----------|
| MTC | Off |
| Flip angle | 60 deg |
| Fat suppr. | Fat sat. |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 19 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution

| | |
|-----------------------|----------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Interpolation | Off |
| PAT mode | GRAPPA |
| Accel. factor PE | 2 |
| Ref. lines PE | 24 |
| Reference scan mode | Separate |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Raw filter | Off |
| Elliptical filter | Off |
| Hamming | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
| Series | Ascending |

Special sat. None

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|--------------------------|------------------|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |
| Auto Coil Select | Default |
| Shim mode | Standard |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | L0.0 A27.8 H28.5 |
| Orientation | T > C-15.3 |
| Rotation | 0.00 deg |
| R >> L | 200 mm |
| A >> P | 200 mm |
| F >> H | 96 mm |

Physio

1st Signal/Mode None

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | On |
| Dynamic t-maps | Off |
| Starting ignore meas | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 19 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Meas[17] | Active |
| Meas[18] | Active |
| Meas[19] | Active |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-------------------|------------|
| Motion correction | On |
| Interpolation | 3D-K-space |
| Spatial filter | Off |

Sequence

| | |
|-------------------|------------|
| Introduction | Off |
| Bandwidth | 1954 Hz/Px |
| Free echo spacing | Off |
| Echo spacing | 0.64 ms |
| <hr/> | |
| EPI factor | 128 |
| RF pulse type | Normal |
| Gradient mode | Fast |
| RF spoiling | On |
| <hr/> | |
| RF90 duration | 5120 |
| MB Number | 4 |
| DummyScan Number | 3 |
| FOV Shift Number | 3 |
| Shift K0 Center | 1 |
| Every Other Slice | 1 |
| SER Number | 1 |
| 2nd RFOff(1) | 0 |
| Polarity(1) | 0 |
| Dephase(0) | 0 |
| Echo Distance | 1.00 |
| MB Measurements | 4 |
| Ramp On | On |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\liceTesting\ep2d_bold_OVS2_m2

TA: 1:16 PAT: 2 Voxel size: 1.7x1.7x1.5 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS2

Properties

| | |
|-----------------------------------------------|--------|
| Prio Recon | Off |
| Before measurement | |
| After measurement | |
| Load to viewer | On |
| Inline movie | Off |
| Auto store images | On |
| Load to stamp segments | Off |
| Load images to graphic segments | Off |
| Auto open inline display | Off |
| Start measurement without further preparation | On |
| Wait for user to start | Off |
| Start measurements | single |

Routine

| | |
|--------------------|------------------|
| Slice group 1 | |
| Slices | 48 |
| Dist. factor | 0 % |
| Position | L0.0 A33.6 H27.7 |
| Orientation | Transversal |
| Phase enc. dir. | A >> P |
| Rotation | 0.00 deg |
| Phase oversampling | 0 % |
| FoV read | 212 mm |
| FoV phase | 62.5 % |
| Slice thickness | 1.50 mm |
| TR | 4760 ms |
| TE | 20 ms |
| Averages | 1 |
| Concatenations | 1 |
| Filter | None |
| Coil elements | B4;M2,3;T1 |

Contrast

| | |
|-----------------|-----------|
| MTC | Off |
| Flip angle | 90 deg |
| Fat suppr. | Fat sat. |
| Averaging mode | Long term |
| Reconstruction | Magnitude |
| Measurements | 13 |
| Delay in TR | 0 ms |
| Multiple series | Off |

Resolution

| | |
|-----------------------|----------|
| Base resolution | 128 |
| Phase resolution | 100 % |
| Phase partial Fourier | 6/8 |
| Interpolation | Off |
| PAT mode | GRAPPA |
| Accel. factor PE | 2 |
| Ref. lines PE | 24 |
| Reference scan mode | Separate |
| Distortion Corr. | Off |
| Prescan Normalize | Off |
| Raw filter | On |
| Elliptical filter | Off |
| Hamming | Off |

Geometry

| | |
|------------------|-------------|
| Multi-slice mode | Interleaved |
| Series | Ascending |

Sat. region 1

| | |
|-------------|------------------|
| Thickness | 50 mm |
| Position | L0.0 A116.8 H0.0 |
| Orientation | Coronal |

Sat. region 2

| | |
|--------------|-----------------|
| Thickness | 50 mm |
| Position | L0.0 P46.6 H0.0 |
| Orientation | Coronal |
| Special sat. | None |

| | |
|------------------|------|
| Table position | H |
| Table position | 0 mm |
| Inline Composing | Off |

System

| | |
|-----|-----|
| T1 | On |
| M2 | On |
| B4 | On |
| M3 | On |
| V32 | Off |

| | |
|-------------------|----------------|
| Positioning mode | REF |
| MSMA | S - C - T |
| Sagittal | R >> L |
| Coronal | A >> P |
| Transversal | F >> H |
| Coil Combine Mode | Sum of Squares |
| AutoAlign | --- |
| Auto Coil Select | Default |

| | |
|--------------------------|------------------|
| Shim mode | Standard |
| Adjust with body coil | Off |
| Confirm freq. adjustment | Off |
| Assume Silicone | Off |
| ? Ref. amplitude 1H | 0.000 V |
| Adjustment Tolerance | Auto |
| Adjust volume | |
| Position | L0.0 A33.6 H27.7 |
| Orientation | Transversal |
| Rotation | 0.00 deg |
| R >> L | 212 mm |
| A >> P | 133 mm |
| F >> H | 72 mm |

Physio

| | |
|-----------------|------|
| 1st Signal/Mode | None |
|-----------------|------|

BOLD

| | |
|-------------------------|----------|
| GLM Statistics | Off |
| Dynamic t-maps | Off |
| Starting ignore meas | 0 |
| Ignore after transition | 0 |
| Model transition states | On |
| Temp. highpass filter | On |
| Threshold | 4.00 |
| Paradigm size | 16 |
| Meas[1] | Baseline |
| Meas[2] | Baseline |
| Meas[3] | Baseline |
| Meas[4] | Baseline |
| Meas[5] | Baseline |
| Meas[6] | Baseline |
| Meas[7] | Baseline |
| Meas[8] | Baseline |
| Meas[9] | Baseline |
| Meas[10] | Baseline |
| Meas[11] | Active |

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

| | |
|-------------------|------------|
| Meas[12] | Active |
| Meas[13] | Active |
| Meas[14] | Active |
| Meas[15] | Active |
| Meas[16] | Active |
| Motion correction | On |
| Interpolation | 3D-K-space |
| Spatial filter | Off |

Sequence

| | |
|-------------------|------------|
| Introduction | Off |
| Bandwidth | 2056 Hz/Px |
| Free echo spacing | Off |
| Echo spacing | 0.59 ms |
| <hr/> | |
| EPI factor | 80 |
| RF pulse type | Normal |
| Gradient mode | Fast |
| RF spoiling | On |
| <hr/> | |
| RF90 duration | 5120 |
| MB Number | 2 |
| DummyScan Number | 3 |
| FOV Shift Number | 2 |
| Shift K0 Center | 1 |
| Every Other Slice | 1 |
| SER Number | 1 |
| Spoil factor | 5 |
| Skew Direction | 0 |
| Sat RF90 duration | 5120 |
| Dual On(1) | 1 |
| Echo Distance | 1.00 |
| MB Measurements | 4 |
| Ramp On | On |