

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2_OVS_flash_iso55_pkk

TA: 8:36

PAT: 2

Voxel size: 0.5x0.5x0.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	140 mm
FoV phase	50.0 %
Slice thickness	0.55 mm
TR	4000 ms
TE	28 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	126
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A43.2 H1.5
Orientation	C > T2.0

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L12.7 P37.1 H14.7
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	70 mm
! F >> H	48 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Allowed
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.08 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	7680
Dual On(1)	3
Echo Distance	1.00
MB Measurements	123
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2_OVS_flash_iso55

TA: 2:04

PAT: 2

Voxel size: 0.5x0.5x0.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	140 mm
FoV phase	50.0 %
Slice thickness	0.55 mm
TR	4000 ms
TE	34 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	28
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A43.2 H1.5
Orientation	C > T2.0

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L12.7 P37.1 H14.7
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	70 mm
! F >> H	48 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.4 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	25
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2_OVS_flash_iso75

TA: 2:04

PAT: 2

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	50.0 %
Slice thickness	0.75 mm
TR	4000 ms
TE	27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	28
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A43.2 H1.5
Orientation	C > T2.0

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L12.7 P37.1 H14.7
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	70 mm
! F >> H	48 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	25
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\fl_fq_mb

TA: 0:21

PAT: Off

Voxel size: 1.6x1.6x5.0 mm

Rel. SNR: 1.00

USER: fl_fq_mb

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	150.00 ms
TE	10.00 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	15 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Single dir.
Encodings	1
Velocity enc.	90 cm/s
Direction	Through plane
Rephased images	On
Magnitude images	On
Phase images	On
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Asymmetric echo	Off
Contrasts	1
Bandwidth	260 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
MB Number	2
FOV Shift	1

| FOV Shift

20

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\t2_fl2d_sag_hemo

TA: 0:48 PAT: 2 Voxel size: 1.1x0.9x4.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	35
Dist. factor	30 %
Position	L0.0 A13.6 H13.6
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	400 ms
TE	7.00 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	E01-20

Contrast

MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	80 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On

Mode

Inplane

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Bandwidth	200 Hz/Px
Flow comp.	Slice/Read
<hr/>	
RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\t2_fl2d_cor_hemo

TA: 0:48 PAT: 2 Voxel size: 1.1x0.9x2.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	35
Dist. factor	10 %
Position	R7.5 P35.3 H13.6
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	400 ms
TE	7.00 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	E01-20

Contrast

MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	80 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On

Mode

Inplane

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off

Shim mode

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Bandwidth	200 Hz/Px
Flow comp.	Slice/Read
<hr/>	
RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\t2_fl2d_cor_hemo

TA: 0:48 PAT: 2 Voxel size: 1.1x0.9x2.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	25
Dist. factor	10 %
Position	L0.0 A13.6 H13.6
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	2.0 mm
TR	400 ms
TE	7.00 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	E01-20

Contrast

MTC	Off
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	80 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On

Mode

Inplane

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	P >> A
Transversal	H >> F
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off

Shim mode

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Bandwidth	200 Hz/Px
Flow comp.	Slice/Read
<hr/>	
RF pulse type	Fast
Gradient mode	Fast*
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\NoiseMeasSensitivityMap

TA: 0:23

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	Off
GFactor	Off
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\localizer_200V

TA: 0:27

PAT: Off

Voxel size: 1.2x1.1x3.0 mm

Rel. SNR: 1.00

SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	B4;M2,3;T1

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Physio

1st Signal/Mode	None
Segments	1

Tagging	None
Dark blood	Off

Resp. control	Off
---------------	-----

Resolution

Base resolution	256
-----------------	-----

Inline

Subtract	Off
----------	-----

Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
<hr/>	
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
<hr/>	
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
<hr/>	
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\b1map_200V_TR100

TA: 0:32

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	150 %
Position	L0.0 A5.4 H40.7
Orientation	T > C-17.8
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	425 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	R1.4 A11.5 H37.3
Orientation	T > C-18.6
Rotation	0.00 deg
Base size phase	129 mm
Base size read	87 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On

M3 On
V32 Off

Positioning mode FIX
MSMA S - C - T
Sagittal R >> L
Coronal A >> P
Transversal F >> H
Save uncombined Off
Coil Combine Mode Adaptive Combine
AutoAlign ---
Auto Coil Select Default

Shim mode Tune up
Adjust with body coil Off
Confirm freq. adjustment Off
Assume Silicone Off
! Ref. amplitude 1H 200.000 V
Adjustment Tolerance Auto
Adjust volume
Position Isocenter
Orientation Transversal
Rotation 0.00 deg
R >> L 350 mm
A >> P 263 mm
F >> H 350 mm

Composing

Sequence

Contrasts 2
Bandwidth 260.416667 Hz/Px
T1 Compensation Mean T1
Mean T1 500.0 ms
Angles 1
Amplitude Weighting Linear
Scale Bar Enabled
Raw Data Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\NoiseCovMap

TA: 0:23

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	Off
GFactor	Off
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\gFactorMap

TA: 0:23

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\SnrMap

TA: 0:23

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None
Table position	H

Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
Gradient mode	Fast
RF spoiling	On
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	Off
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\localizer_200V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	E01-20

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On

PAT mode	None
----------	------

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Saturation mode	Standard
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

Tim CT mode	Off
-------------	-----

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal

Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\gFactorMap

TA: 0:46

Voxel size: 0.5x0.5x5.0 mm

Rel. SNR: 1.00

USER: NoiseMeasSensitivityMap

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	3
Dist. factor	400 %
Position	L0.0 A2.0 H2.7
Orientation	T > C-2.9
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	30 ms
TE	6.0 ms
Averages	1
Concatenations	3
Filter	None
Coil elements	E01-20

Contrast

TD	0 ms
MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	384
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Special sat.	None

Table position

H

Table position

0 mm

Inline Composing

Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

REF

MSMA

S - C - T

Sagittal

R >> L

Coronal

A >> P

Transversal

F >> H

Save uncombined

Off

Coil Combine Mode

Adaptive Combine

AutoAlign

Auto Coil Select

Default

Shim mode

Tune up

Adjust with body coil

Off

Confirm freq. adjustment

Off

Assume Silicone

Off

? Ref. amplitude 1H

0.000 V

Adjustment Tolerance

Auto

Adjust volume

! Position

Isocenter

! Orientation

Transversal

! Rotation

0.00 deg

! R >> L

350 mm

! A >> P

263 mm

! F >> H

350 mm

Physio

1st Signal/Mode	None
-----------------	------

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Introduction	Off
Dimension	2D
Contrasts	1
Bandwidth	200 Hz/Px
<hr/>	
Gradient mode	Fast
RF spoiling	On
<hr/>	
ICE program	CoilArrayUtil
number of noise lines	384 lines
Optimal SNR	On
GFactor	On
Condition number	Off
Rx coil diode switching	On
coil channel reordering	Off
<hr/>	
TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2f1_OVS_flash_iso75

TA: 0:28

PAT: 2

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P48.6 H1.1
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	50.0 %
Slice thickness	0.75 mm
TR	4000 ms
TE	27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	4
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A33.7 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 P34.6 H0.7
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	80 mm
! F >> H	56 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	1
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2_OVS_flash_iso55

TA: 0:28

PAT: 2

Voxel size: 0.5x0.5x0.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P48.6 H1.1
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	140 mm
FoV phase	50.0 %
Slice thickness	0.55 mm
TR	4000 ms
TE	34 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	4
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A33.7 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 P34.6 H0.7
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	80 mm
! F >> H	56 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.4 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	1
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2_OVS_flash_iso46

TA: 0:40

PAT: 2

Voxel size: 0.5x0.5x0.5 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P58.1 H0.5
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	147 mm
FoV phase	25.0 %
Slice thickness	0.45 mm
TR	4000 ms
TE	29 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	7
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	14
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A7.3 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L0.0 P34.6 H0.7
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	80 mm
! F >> H	56 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.56 ms
<hr/>	
EPI factor	80
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On
<hr/>	
RF90 duration	7680
MB Number	2
DummyScan Number	2
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	1
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2Px2_OVS_flash_iso46_fmri

TA: 4:36

PAT: 2

Voxel size: 0.5x0.5x0.5 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P75.8 F5.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	147 mm
FoV phase	25.0 %
Slice thickness	0.45 mm
TR	4000 ms
TE	25 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	66
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	320
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	14
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 P4.9 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 P64.2 F7.0
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	147 mm
! A >> P	54 mm
! F >> H	44 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline

Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Allowed
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.28 ms
<hr/>	
EPI factor	80
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On
<hr/>	
RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	63
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2P2f1_OVS_flash_iso75-fmri

TA: 4:36

PAT: 2

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P50.8 F5.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	50.0 %
Slice thickness	0.75 mm
TR	4000 ms
TE	27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	66
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A48.0 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 P64.2 F7.0
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	147 mm
! A >> P	54 mm
! F >> H	44 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline

Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms
<hr/>	
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	63
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\0317Test\ep2d_M2Px2_OVS_flash_iso55_fmri

TA: 4:36

PAT: 2

Voxel size: 0.5x0.5x0.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P65.8 F5.6
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	140 mm
FoV phase	50.0 %
Slice thickness	0.55 mm
TR	4000 ms
TE	28 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B1-8

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	66
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A20.9 H0.0
Orientation	Coronal

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

B1	On
B2	On
B3	On
B4	On
B5	On
B6	On
B7	On
B8	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.9 P64.2 F7.0
! Orientation	Transversal
! Rotation	0.00 deg
! R >> L	147 mm
! A >> P	54 mm
! F >> H	44 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Allowed
Bandwidth	850 Hz/Px
Free echo spacing	Off
Echo spacing	1.07 ms
<hr/>	
EPI factor	128
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On
<hr/>	
RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	63
Ramp On	On

Table of contents

\\USER

Feinberglab

Test

0317Test

```

ep2d_M2P2_OVS_flash_iso55_pbx
ep2d_M2P2_OVS_flash_iso55
ep2d_M2P2_OVS_flash_iso75
fl_fq_mb
t2_fl2d_sag_hemo
t2_fl2d_cor_hemo
t2_fl2d_cor_hemo
NoiseMeasSensitivityMap
-----run nova coil-----
localizer_200V
b1map_200V_TR100
NoiseCovMap
gFactorMap
SnrMap
-----run new coil-----
localizer_200V_newcoil
gFactorMap
ep2d_M2P2f1_OVS_flash_iso75
ep2d_M2P2_OVS_flash_iso55
-----run new coil-----
ep2d_M2P2_OVS_flash_iso46
ep2d_M2Px2_OVS_flash_iso46_fmri
ep2d_M2P2f1_OVS_flash_iso75_fmri
ep2d_M2Px2_OVS_flash_iso55_fmri

```