

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\M2F4P4_flashref_fa60

TA: 0:19 PAT: 4 Voxel size: 0.9x0.9x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flg

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	36
Dist. factor	50 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	232 mm
FoV phase	100.0 %
Slice thickness	0.8 mm
TR	1750 ms
TE	27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A21.2 H18.8
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1302 Hz/Px
Free echo spacing	Off
Echo spacing	0.91 ms
<hr/>	
EPI factor	256
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	4
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFOff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0

TA: 0:28 PAT: 4 Voxel size: 0.9x0.9x3.0 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	200 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	81.3 %
Slice thickness	3.00 mm
TR	2000 ms
TE	23 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	8
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	50 mm
Position	L0.0 A65.8 H0.0
Orientation	Coronal

Sat. region 2

Thickness	50 mm
Position	L0.0 P60.0 H0.0
Orientation	Coronal
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Rotation	0.00 deg
R >> L	220 mm
A >> P	179 mm
F >> H	102 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1220 Hz/Px
Free echo spacing	Off
Echo spacing	0.94 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	2
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	5
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	1
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_m3F3_pcasl_DE_BS_flash

TA: 0:42 PAT: 2 Voxel size: 3.1x3.1x5.0 mm Rel. SNR: 1.00 USER: ep2d_fid_mb_pcasl_DE_BS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	50 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000 ms
TE 1	14 ms
TE 2	31 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	11
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	200 mm
A >> P	200 mm
F >> H	88 mm

Physio

1st Signal/Mode	None
-----------------	------

Perf

GBP	Off
PBP	Off
TTP	Off
Original images	On

Sequence

Introduction	Off
Contrasts	2
Bandwidth	3004 Hz/Px
Free echo spacing	Off
Echo spacing	0.53 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
RF90 duration	5120
MB Number	3
DummyScan Number	2
FOV Shift Number	3
Shift K0 Center	1
Every Other Slice	1
BS Type	1
Background Suppr.	Off
BS parameter[1]	500 ms
BS parameter[2]	10 ms

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

BS parameter[3]	0 ms
FOCI parameter[1]	800
FOCI parameter[2]	24
FOCI parameter[3]	1.0
FOCI parameter[4]	2000
Label Offset	80 mm
Post Label Delay	1000000 us
MB Measurements	3
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\M2F2P4_flashref_fa60

TA: 0:20 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flk

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	48
Dist. factor	50 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	81.3 %
Slice thickness	0.8 mm
TR	1820 ms
TE	25 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A21.2 H18.8
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFoff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\M2F2P4_flashref_fa60

TA: 0:20 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flg

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	48
Dist. factor	50 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	81.3 %
Slice thickness	0.80 mm
TR	1820 ms
TE	25 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A21.2 H18.8
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	976 Hz/Px
Free echo spacing	Off
Echo spacing	1.15 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFOff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\M2F2P4_flashref_fa60-S36-TR1500

TA: 0:17 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flk

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	36
Dist. factor	100 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	81.3 %
Slice thickness	0.80 mm
TR	1500 ms
TE	25 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A21.2 H18.8
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	976 Hz/Px
Free echo spacing	Off
Echo spacing	1.15 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFOff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRMM2F2P4_flashref_fa60_S36_TR1500_AP

TA: 0:17 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flk

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	36
Dist. factor	100 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.80 mm
TR	1500 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	48
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat. None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Rotation	0.00 deg
R >> L	200 mm
A >> P	200 mm
F >> H	57 mm

Physio

1st Signal/Mode None

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.03 ms
<hr/>	
EPI factor	256
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFoff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0

TA: 0:22 PAT: 4 Voxel size: 0.9x0.9x3.0 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	200 %
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	220 mm
FoV phase	81.3 %
Slice thickness	3.00 mm
TR	2000 ms
TE	23 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	50 mm
Position	L0.0 A98.8 H25.9
Orientation	C > T14.7

Sat. region 2

Thickness	50 mm
Position	L0.0 P57.6 F15.1
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.2 A16.4 H26.2
Orientation	T > C-15.0
Rotation	0.00 deg
R >> L	220 mm
A >> P	179 mm
F >> H	102 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1220 Hz/Px
Free echo spacing	Off
Echo spacing	0.94 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	5
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	1
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\t1_mpr_tra_p2_iso

TA: 3:21

PAT: 2

Voxel size: 0.9x0.9x0.9 mm

Rel. SNR: 1.00

SIEMENS: tfl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L0.0 P12.9 F4.1
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	14 %
Slice oversampling	44.4 %
Slices per slab	144
FoV read	230 mm
FoV phase	87.5 %
Slice thickness	0.90 mm
TR	1900 ms
TE	3.32 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	900 ms
Flip angle	9 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	96 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off

Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Ascending
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Asymmetric echo	Allowed
Bandwidth	180 Hz/Px
Flow comp.	No
Echo spacing	8.1 ms
<hr/>	
RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\localizer_200V

TA: 0:27

PAT: Off

Voxel size: 1.2x1.1x3.0 mm

Rel. SNR: 1.00

SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	E01-20

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Off
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
-----------------	-----

Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\b1map_200V_TR100

TA: 0:32

Voxel size: 3.9x3.9x5.0 mm

Rel. SNR: 1.00

USER: b1map_658

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	150 %
Position	L0.0 A5.4 H40.7
Orientation	T > C-17.8
Phase enc. dir.	A >> P
Rotation	0.00 deg
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5 mm
TR	425 ms
TE 1	14 ms
TE 2	14 ms
Averages	1
Filter	None
Coil elements	

Contrast

Flip angle 1	90 deg
Flip angle 2	120 deg
Flip angle 3	60 deg
Flip angle 4	135 deg
Flip angle 5	45 deg
Measurements	1

Resolution

Base resolution	64
Phase resolution	100 %
Raw filter	Off

Geometry

Series	Interleaved
Navigator 1	
Position	R1.4 A11.5 H37.3
Orientation	T > C-18.6
Rotation	0.00 deg
Base size phase	129 mm
Base size read	87 mm
Thickness	50 mm
Table position	H
Table position	0 mm
Inline Composing	Off

System

Positioning mode	FIX
MSMA	S - C - T

Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	200.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Composing

Sequence

Contrasts	2
Bandwidth	260.416667 Hz/Px
T1 Compensation	Mean T1
Mean T1	500.0 ms
Angles	1
Amplitude Weighting	Linear
Scale Bar	Enabled
Raw Data	Disabled

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\M2F2P4_flashref_fa60

TA: 0:48 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_fl

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	48
Dist. factor	50 %
Position	L1.2 A17.8 H50.0
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	81.3 %
Slice thickness	0.8 mm
TR	4320 ms
TE	25 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline

Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms
<hr/>	
EPI factor	208
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFOff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRMM2F2P4_flashref_fa60_S36_TR1500_AP

TA: 0:36 PAT: 4 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_sbmb_cte_ipat_fov_asym_flk

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	36
Dist. factor	100 %
Position	L1.2 A17.8 H50.0
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	0.80 mm
TR	3240 ms
TE	26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	36
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Special sat.

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline

Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.03 ms
<hr/>	
EPI factor	256
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
2nd RFOff(1)	0
Polarity(1)	0
Dephase(0)	0
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0

TA: 0:11

PAT: 4

Voxel size: 0.8x0.8x0.8 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	200 %
Position	L1.2 A17.8 H50.0
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	40.6 %
Slice thickness	0.75 mm
TR	1000 ms
TE	28 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	50 mm
Position	L0.0 A98.8 H25.9
Orientation	C > T14.7

Sat. region 2

Thickness	50 mm
Position	L0.0 P57.6 F15.1
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1220 Hz/Px
Free echo spacing	Off
Echo spacing	1.02 ms

EPI factor	104
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	9120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	5
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	1
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0_iso6

TA: 0:22 PAT: 4 Voxel size: 0.6x0.6x0.6 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	24
Dist. factor	115 %
Position	L0.5 P52.7 H18.1
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	160 mm
FoV phase	31.3 %
Slice thickness	0.60 mm
TR	2000 ms
TE	17 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	80 mm
Position	L0.0 A14.5 H3.8
Orientation	C > T14.7

Sat. region 2

Thickness	80 mm
Position	L0.0 P120.6 F31.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	976 Hz/Px
Free echo spacing	Off
Echo spacing	1.23 ms
<hr/>	
EPI factor	80
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	0
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0_iso7

TA: 0:22 PAT: 4 Voxel size: 0.7x0.7x0.7 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	24
Dist. factor	115 %
Position	L0.5 P45.3 H19.4
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	180 mm
FoV phase	31.3 %
Slice thickness	0.70 mm
TR	2000 ms
TE	16 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	50 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	5
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	80 mm
Position	L0.0 A26.1 H6.9
Orientation	C > T14.7

Sat. region 2

Thickness	80 mm
Position	L0.0 P120.6 F31.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.5 P45.3 H19.4
Orientation	T > C-15.0
Rotation	0.00 deg
R >> L	180 mm
A >> P	57 mm
F >> H	36 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.13 ms
<hr/>	
EPI factor	80
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	1
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	0
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	2
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_bold_OVS2_m2

TA: 1:11 PAT: 2 Voxel size: 0.7x0.7x0.6 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS2

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	400 %
Position	L1.2 P45.3 H18.1
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	180 mm
FoV phase	50.0 %
Slice thickness	0.6 mm
TR	4760 ms
TE	32 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	12
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A46.9 H12.3
Orientation	C > T14.7

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	16
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1028 Hz/Px
Free echo spacing	Off
Echo spacing	1.1 ms
<hr/>	
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	3
FOV Shift Number	2
Shift K0 Center	1
Every Other Slice	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	3
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P4_OVS_flash_D1K0_fsovs0

TA: 1:11

PAT: 2

Voxel size: 0.7x0.7x0.6 mm

Rel. SNR: 1.00

USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	12
Dist. factor	400 %
Position	L1.2 P45.3 H18.1
Orientation	T > C-15.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	180 mm
FoV phase	50.0 %
Slice thickness	0.60 mm
TR	4760 ms
TE	32 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	12
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A46.9 H12.3
Orientation	C > T14.7

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L1.2 A17.8 H47.9
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	192 mm
! A >> P	192 mm
! F >> H	72 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1028 Hz/Px
Free echo spacing	Off
Echo spacing	1.1 ms
<hr/>	
EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
<hr/>	
RF90 duration	5120
MB Number	2
DummyScan Number	3
FOV Shift Number	2
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	3
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P2_OVS_flash_iso55

TA: 1:37 PAT: 2 Voxel size: 0.5x0.5x0.6 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	140 mm
FoV phase	50.0 %
Slice thickness	0.55 mm
TR	3120 ms
TE	34 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	28
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A43.2 H1.5
Orientation	C > T2.0

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Rotation	0.00 deg
R >> L	140 mm
A >> P	70 mm
F >> H	41 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	782 Hz/Px
Free echo spacing	Off
Echo spacing	1.4 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	25
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_M2P2_OVS_flash_iso75

TA: 1:37 PAT: 2 Voxel size: 0.8x0.8x0.8 mm Rel. SNR: 1.00 USER: ep2d_bold_OVS_flash

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	50
Dist. factor	50 %
Position	L1.2 P44.6 H15.4
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	192 mm
FoV phase	50.0 %
Slice thickness	0.75 mm
TR	3120 ms
TE	27 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	E01-20

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	28
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Reference scan mode	Separate
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending

Sat. region 1

Thickness	110 mm
Position	L0.0 A43.2 H1.5
Orientation	C > T2.0

Sat. region 2

Thickness	110 mm
Position	L0.0 P136.2 F35.7
Orientation	C > T14.7
Special sat.	None

Table position

Table position	H
Table position	0 mm
Inline Composing	Off

System

E17	On
E18	On
E19	On
E20	On
E01	On
E02	On
E03	On
E04	On
E05	On
E06	On
E07	On
E08	On
E09	On
E10	On
E11	On
E12	On
E13	On
E14	On
E15	On
E16	On

Positioning mode

MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode

Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
! Position	L12.7 P37.1 H14.7
! Orientation	T > C-15.0
! Rotation	0.00 deg
! R >> L	100 mm
! A >> P	70 mm
! F >> H	48 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off

Sequence

Introduction	Off
Asymmetric echo	Off
Bandwidth	1086 Hz/Px
Free echo spacing	Off
Echo spacing	1.05 ms

EPI factor	128
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

RF90 duration	7680
MB Number	2
DummyScan Number	1
FOV Shift Number	1
SkewType(1ff)	0
OVS flash(1on)	1
SER Number	1
Spoil factor	1
Skew Direction	1
Sat RF90 duration	5120
Dual On(1)	3
Echo Distance	1.00
MB Measurements	25
Ramp On	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\ISMRM\ep2d_venc_ms_sbmb_SAT_flashref

TA: 0.2 s PAT: Off Voxel size: 7.8x3.9x5.0 mm Rel. SNR: 1.00 USER: ep2d_venc_ms_sbmb_SAT_flashref

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	200 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0 deg
Phase oversampling	0 %
FoV read	500 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	79 ms
TE	1 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	T1

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	20
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	50 %
Phase partial Fourier	5/8
Interpolation	Off
PAT mode	None
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm

Inline Composing

Off

System

T1	On
M2	Off
B4	Off
M3	Off
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	500 mm
A >> P	500 mm
F >> H	5 mm

Physio

1st Signal/Mode	None
-----------------	------

Angio

Flow mode	Single dir.
Encodings	1
Velocity enc.	90 cm/s
Direction	Through plane
Magnitude sum	Off

Sequence

Introduction	Off
Bandwidth	752 Hz/Px
Free echo spacing	Off
Echo spacing	0 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On
RF90 duration	5120
MB Number	1
DummyScan Number	1
FOV Shift Number	1
Shift K0 Center	1
Every Other Slice	1
SER Number	1
Venc Repetition	1
Spoil factor	5
Skew Direction	0
Dual On(1)	0
Venc Type(0off,1+,-,20+,3on)	0

Table of contents

\\USER

Feinberglab

Test

ISMIRM

M2F4P4_flashref_fa60
 ep2d_M2P4_OVS_flash_D1K0
 ep2d_m3F3_pcasl_DE_BS_flash
 =====nova coil=====
 M2F2P4_flashref_fa60
 M2F2P4_flashref_fa60
 M2F2P4_flashref_fa60-S36-TR1500
 M2F2P4_flashref_fa60_S36_TR1500_AP
 ep2d_M2P4_OVS_flash_D1K0
 t1_mpr_tra_p2_iso
 =====new coil=====
 localizer_200V
 b1map_200V_TR100
 M2F2P4_flashref_fa60
 M2F2P4_flashref_fa60_S36_TR1500_AP
 ep2d_M2P4_OVS_flash_D1K0
 =====
 ep2d_M2P4_OVS_flash_D1K0_iso6
 ep2d_M2P4_OVS_flash_D1K0_iso7
 =====compare flash/epi reference=====
 ep2d_bold_OVS2_m2
 ep2d_M2P4_OVS_flash_D1K0_fsovs0
 ep2d_M2P2_OVS_flash_iso55
 ep2d_M2P2_OVS_flash_iso75
 ep2d_venc_ms_sbmb_SAT_flashref