

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\Testing\localizer_100V_newcoil

TA: 0:27 PAT: Off Voxel size: 1.2x1.1x3.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Phase resolution	90 %
Phase partial Fourier	6/8
Interpolation	On
PAT mode	None
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off
Tim CT mode	Off

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Slice group 3	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	3.0 mm
TR	10.0 ms
TE	3.00 ms
Averages	1
Concatenations	15
Filter	None
Coil elements	L1-8;R1-8;T1-8

System

T1	On
T2	On
T3	On
T4	On
T5	On
T6	On
T7	On
T8	On
L1	On
L2	On
L3	On
L4	On
L5	On
L6	On
L7	On
L8	On
R1	On
R2	On
R3	On
R4	On
R5	On
R6	On
R7	On
R8	On

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Positioning mode	FIX
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	On
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Off

Resolution

Base resolution	256
-----------------	-----

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
! Ref. amplitude 1H	50.000 V

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Tagging	None
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
MapIt	None
Contrasts	1

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Allowed
Bandwidth	320 Hz/Px
Flow comp.	No
RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\Testing\BP_grase_clean_VASO

TA: 6.0 s PAT: 4 Voxel size: 2.7x2.0x3.0 mm Rel. SNR: 1.00 USER: BP_grase_clean_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	20
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	3.0 mm
TR	3000 ms
TE	26.36 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1000 ms
Flip angle	100 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	75 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	2
Ref. lines 3D	20
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off

Geometry

Series

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	90.00 deg
A >> P	256 mm
R >> L	192 mm
F >> H	60 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	3256 Hz/Px
Echo spacing	0.6 ms
Turbo factor	7
EPI factor	72
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	Mosaic
prepscans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\Testing\BP_grase_clean_VASO

TA: 0:15

PAT: 12

Voxel size: 2.0x2.0x3.0 mm

Rel. SNR: 1.00

USER: BP_grase_clean_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R1.3 P0.0 H31.7
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	12.5 %
Slices per slab	32
FoV read	256 mm
FoV phase	75.0 %
Slice thickness	3.0 mm
TR	3000 ms
TE	18.42 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	1000 ms
Flip angle	100 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	4
Ref. lines PE	24
Accel. factor 3D	3
Ref. lines 3D	33
Reference scan mode	Separate
Prescan Normalize	Off
Raw filter	Off

Geometry

Series

Interleaved

Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R1.3 P0.0 H31.7
Orientation	Transversal
Rotation	90.00 deg
A >> P	256 mm
R >> L	192 mm
F >> H	96 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	2170 Hz/Px
Echo spacing	0.6 ms
Turbo factor	12
EPI factor	96
RF pulse type	Normal
Gradient mode	Fast
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prep scans	0

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\Testing\pgrs3d_asl_collection_new_7T

TA: 3:20

Voxel size: 1.7x1.7x4.0 mm

Rel. SNR: 1.00

USER: pgrs3d_asl_collection_new_7T

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	R2.0 A29.7 H9.4
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	26
FoV read	212 mm
FoV phase	75.0 %
Slice thickness	4.0 mm
TR	5000 ms
TE	20.88 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Flip angle	100 deg
Fat suppr.	Fat sat.
Fat sat. mode	Strong
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	2
Pause after meas. 1	0.0 s
Multiple series	Each measurement

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	6/8
Interpolation	On
Prescan Normalize	Off
Raw filter	Off

Geometry

Series	Interleaved
Sat. region 1	
Thickness	17 mm
Position	L0.0 P0.0 F67.7
Orientation	Transversal
Special sat.	None

Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R2.0 A29.7 H9.4
Orientation	Transversal
Rotation	90.00 deg
A >> P	212 mm
R >> L	159 mm
F >> H	104 mm

Physio

1st Signal/Mode	None
Segments	10

Composing

Sequence

Introduction	Off
Dimension	3D
Reordering	Centric
Contrasts	1
Bandwidth	2056 Hz/Px
Echo spacing	0.7 ms
Turbo factor	10
EPI factor	20
RF pulse type	Normal
Gradient mode	Fast
ASL mode	TE-pCASL
Pre sat	Off
manual LABEL/SAT region	Off
Robust Hadamard	Off
Saturation mode	var. I
Background Suppr	2*T1-based BackgroundSupression
BS parameter_1	700 ms
BS parameter_2	100 ms
Thickness of BS slice	200 mm
FOCI parameter_0	600
FOCI parameter_1	12
FOCI parameter_2	1.0
FOCI parameter_3	300

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

FOCI parameter_4	0.500
Additional scaling factor	5.0
Distribution Mode	Equal
z-spoiling	equal
Fix label plane offset	Off
pCASL bolus length	1000 ms
Subbolus length	50 ms
Time encoding steps	2
RF gap	360 usec
RF FA	30 deg
post labeling delay	1000 ms
Number of echoes	1

SIEMENS MAGNETOM Investigational_Device_7T syngo MR B17

\\USER\Feinberglab\Test\Testing\BP_grase_clean_VASO2x2-ti1330

TA: 0:20

PAT: 4

Voxel size: 1.8x1.8x3.0 mm

Rel. SNR: 1.00

USER: BP_grase_clean_VASO

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	0 %
Position	L0.0 A30.8 H1.2
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	12
FoV read	234 mm
FoV phase	75.0 %
Slice thickness	3.0 mm
TR	4000 ms
TE	53.26 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	B4;M2,3;T1

Contrast

Magn. preparation	Non-sel. IR
T1	2000 ms
Flip angle	180 deg
Fat suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	4
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	100 %
Slice partial Fourier	Off
Interpolation	On
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	2
Ref. lines 3D	12
Reference scan mode	Separate
Prescan Normalize	Off

Raw filter

Off

Geometry

Series	Ascending
Special sat.	None
Table position	H
Table position	0 mm
Inline Composing	Off

System

T1	On
M2	On
B4	On
M3	On
V32	Off
Positioning mode	REF
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 A30.8 H1.2
Orientation	Transversal
Rotation	90.00 deg
A >> P	234 mm
R >> L	176 mm
F >> H	36 mm

Physio

1st Signal/Mode	None
-----------------	------

Composing

Sequence

Introduction	Off
Dimensioning	3D
Reordering	Centric
Contrasts	1
Bandwidth	3256 Hz/Px
Echo spacing	1 ms
Turbo factor	6
EPI factor	96
RF pulse type	Normal
Gradient mode	Fast*
refocussing type	sinc 2560
flip angle excit	90
phase encoding	ON
Maxwell compensation	Off
ICE program	single
prepscans	0

Table of contents

\\USER	Feinberglab	Test	Testing	localizer_100V_newcoil
				BP_grase_clean_VASO
				BP_grase_clean_VASO
				pgrs3d_asl_collection_new_7T
				BP_grase_clean_VASO2x2-ti1330